
NT-VAL-Jig

*Gebrauchsanweisung / Instructions for Use / Mode d'emploi /
Istruzioni per l'uso / Instrucciones de uso*

 **DEUTSCH** 02

 **ENGLISH** 09

 **FRANÇAIS** 16

 **ITALIANO** 23

 **ESPAÑOL** 30

DEUTSCH

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ENGLISH

Instructions for Use

SERIES COMPATIBILITY WITH THE IMPLANT SYSTEMS

SERIES (NT-TRADING)	COMPATIBLE IMPLANT SYSTEM
CAM	Altatec® / Camlog®
COL	Altatec® / Conelog®
E	Nobel Biocare™ / Replace Select®
F	Nobel Biocare™ / NobelActive™
H	Biomet 3i® / Certain®
K	Nobel Biocare™ / Brånemark®
L	Straumann / Bone Level®
N	Straumann / SynOcta®
R	Zimmer Dental / Tapered Screw-Vent®
S	DENTSPLY Implants® / Astra Tech Implant System OsseoSpeed® TX
T	DENTSPLY Implants® / Xive®

Table 1: Assignment of NT-VAL-jig to implant systems.

Item numbers for the different diameters can be found in the nt-trading product catalogue.

TORQUES

NCM	SERIES
10	A-X

 **Rework of any type, particularly in relation to the connection geometry to the implant, leads to fitting inaccuracies that can render the device unsuitable for further use.**

INDICATION

The **NT-VAL-Jig** is a position holder fixture for dental implants inside a patient's mouth with a support unit that is used during intraoral application for precise implant position transfer and subsequent precise position checking and correction of implant analogues in the working model and patient model.

The **NT-VAL-Jig** is made up of several parts, comprising transfer posts, fixing elements and a transfer key. All components of the **NT-VAL-Jig** are mechanically fixed below one another using screw connections.

MATERIAL

Ti base and abutment screw: Ti6Al4V, medical grade 5, ASTM F136

CONTRAINDICATION

The transfer posts of the **NT-VAL-Jig** in each series can only be used with the corresponding compatible implant system, e.g. the E series can only be used in combination with the Replace Select® implant system. They cannot be combined with implants or implant analogues of a different implant type or from a different implant manufacturer. Transfer posts with an incorrect diameter must not be used. This can lead to irritation of the peri-implant soft tissue and cause damage to the implant.

PROCESSING

The transfer post for the **NT-VAL-Jig** is compatible with the corresponding implant systems shown in Table 1. The product is suitable for disinfection and sterilisation. The expiration date of the product is printed on the label.

Always read these instructions for use before using the **NT-VAL-Jig**. The **NT-VAL-Jig** may be used only depending on the indication in accordance with the general principles of medical/surgical treatment and in consideration of the applicable occupational safety and accident prevention regulations. **NT-VAL-Jig** is only suitable for medical/dental use when used together with the relevant implant system. If any points relating to the indication or the type of application are unclear, do not use the product until all issues have been clarified.

As part of our terms of sale and delivery, we guarantee perfect quality of our products. The following descriptions are not sufficient to enable practitioners inexperienced in implant prosthetics to perform the work correctly. The practitioner must be experienced in dental surgery and prosthetics, including diagnostics and pre-operative planning and/or laboratory procedures. Because the way in which our product is used is beyond our control, all liability for resultant damage is excluded.


The practitioner bears sole responsibility for

the use of the product. The product must be inspected before use to ensure it is not damaged. In the case of damaged packaging, the product must not be used.

The prefabricated **NT-VAL-Jig** made from grade 5 titanium is used for precise implant position transfer and the subsequent precise position checking and correction of implant analogues in the individual work model and patient model. To perform transfer and checking, in phase 1 the individual components (transfer posts, fixing elements and a transfer key) must be mechanically assembled in a series of steps within the patient's mouth, using screw connections. Prosthetic tool, insert: W 11.VAL.G70, W 11.VAL.G75, W 11.VAL.G85

STEP 1**POSITIONING THE TRANSFER POST, STRAIGHT/ANGLED**


- Ensure good access to the implant site. Check the compatibility of the implant series and that the diameter is correct. Also check the implant angle, to enable selection of the correct transfer posts.
- Before positioning, check and clean the internal connection of the implant.
- Place the transfer post precisely in the implant and tighten the transfer screw so that is hand-tight. (10 Ncm)


In the same way as inserting an abutment, the correct position of the transfer post can be checked with an X-ray inspection.

STEP 2


- Connect the fixing post with a positive fit to the transfer post with the knurled screw (pink) aligned ventrally. Tighten the knurled screw so that it is hand-tight.
- Before connecting, check the cleanliness of the outer hexagon of the transfer post.

**STEP 3**

- Using forceps, push the left lateral end fixing element (blue) and the right lateral fixing element (green), starting from the end and moving forwards, onto the fixing posts with the clamping screw aligned ventrally.

**STEP 4**

- With the clamping screw aligned ventrally, slide the retral fixing element (yellow) with the long connecting piece onto the fixing posts over the lateral fixing elements (end blue/green elements).


**STEP 5**

- With the clamping screw aligned ventrally, slide the frontal fixing element (yellow) with the short connecting piece onto the fixing posts over the lateral fixing elements (end blue/green elements, frontal yellow).

The frontal fixing element can be fitted with a grip to enable better handling when disconnecting from the transfer posts.

Tighten all screws using light force in a clockwise direction. Gradually increase the torque to approx. 10 Ncm to achieve an even, tension-free clamping action on all fixing elements.

If all fixing elements are tightened to the maximum torque immediately after sliding onto the fixing posts, this can create tension, which can be avoided by gradually tightening each screw at intervals after the installation of all fixing elements.


In **phase 2**, the NT-VAL-jig can be removed by following the steps described below.

Prosthetic tool, insert: W 11.VAL.G70, W 11.VAL.G75, W 11.VAL.G85


STEP 1

Following longitudinal marking of the transfer posts and fixing posts across both elements, the fixing elements, which are now tightly connected to each other, can be removed from the transfer posts together after undoing the knurled screws (pink).

STEP 2

Undo the transfer screw, and remove the transfer posts individually from the implant. Insert them individually into the fixing posts (according to the marking), and secure in place with the knurled screw (pink).

To prevent unintentional changing or loosening of the fixing elements, a stable, cushioned container should be used for storage and transport of the assembled NT-VAL-jig.

The use of the NT-VAL-jig in **phase 3** for checking or correcting the position of implant analogues in the work and patient model is performed in the following steps:


STEP 1

Disinfection of the NT-VAL-jig and checking of the fixed elements. Checking and, if necessary, renewing the marking between the transfer posts and fixing posts.

STEP 2

Undo the knurled screws (pink) and remove all transfer posts. Screw a transfer post onto one of the implant analogues for which a position check on the implant model is to be performed. After briefly checking the fit between the transfer posts and the fixing post mounted on the jig, screw another loosened transfer post into an implant analogue to be inspected. This enables the tension-free, positive fit of the fixing posts mounted on the jig to be checked on the transfer post on the implant model.


If a deviation from the implant analogue position in the implant model is identified during referencing with the **NT-VAL-Jig**, a correction is not permitted; to do this, the VAL jig is detached from the implant model.

For a precisely fitting dental implant, it is essential that the dental implant analogues correspond exactly to the situation in the patient's mouth, particularly in terms of angular alignment and their alignment to other dental implant analogues.


STEP 3

For correction, the corresponding implant analogue must be removed from the implant model so that this can be attached to the previously loosened transfer post and connected to the fixing post according to the marking.

The opening in the implant model for mounting the implant analogue must be extended far enough to ensure that no pressure is required to insert the implant analogue using the **NT-VAL-Jig**.

The **NT-VAL-Jig** is screwed onto a correctly positioned implant analogue with installed transfer post using a transfer screw, so that the implant analogue to be corrected can be connected to the implant model in its new position.

Finally, check all implant analogues in the implant model, and correct if necessary.

GENERAL INFORMATION

CAUTION: The following precautionary measures should be taken before or during a treatment.

- Before every procedure, ensure that all required parts, resources and instruments are available in their entirety and in the correct quantity, and are fully functional.
- For your own safety, always wear suitable protective clothing.
- Position the patient to minimise the danger of aspiration of any parts. Secure all parts used inside the patient's mouth against the risk of aspiration and swallowing.

SIDE EFFECTS

Allergies or sensitivities related to the Ti6Al4V alloy, medical grade 5, cannot be excluded in very rare individual cases.

INTERACTIONS

Different types of alloys in the same oral cavity may lead to galvanic reactions in the event of occlusal or approximal contact.

WARRANTY

10 years on the mechanical stability of the **Val-Jig** components depending on correct processing according to our instructions for use. Independently of this, information given verbally, in writing or during practical seminars based upon tests and experience gathered can only be judged as a standard value.

Our products are subject to continual development. Due to this fact, we reserve the right to change or alter our products with regard to construction and constituents.

CLEANING, DISINFECTION AND STERILIZATION

The listed products in this IFU are supplied in non-sterile condition. The components should be cleaned, disinfected and in specific clinical procedures and cases be sterilized, prior to use after they are received from the dental laboratory (no liability on disregard). Effective cleaning and disinfection is an indispensable requirement for effective sterilization of the medical device. Prior to sterilization, please keep medical devices clean when handling in

the laboratory and operator.

Additionally, please pay attention to legal regulations valid for your local areas as well as to the hygienic instructions of your dental practice. This applies particularly to the different guidelines regarding the inactivation of prions. Please observe all manufacturers guidelines for disinfection and cleaning agents with special regard to the concentration, exposure time and temperature. Only pH-neutral disinfection solutions without chlorine, ammoniac and aldehydes and with a proven effectiveness against HBV, HCV, and HIV must be used. The products have to meet the respective national regulations for disinfectants. If disinfectants containing aldehydes are used, this might lead to a possible fixation of proteins. Please use only freshly prepared solutions.

1. PRE-DISINFECTION (AVOIDANCE OF CROSS CONTAMINATIONS)

Place the devices in a germicidal bath* immediately after use. Remove all residues and disassemble demountable products.

2. CLEANING

Please use distilled water and neutral cleaning agents only. The internal irrigation tube has to be cleaned with a Miller needle and must be

rinsed with distilled water at the beginning and end of the exposure time using a disposable syringe (min. 10 ml). The products must be cleaned with a plastic instrument cleaning brush and then rinsed with distilled water. Please control all products after cleaning in order to avoid either damaging or corrosion. Damaged products must be replaced.

3. RINSING AND DRYING

After removal of the products from the germicidal bath, all components must be rinsed 3 times with distilled water/ e.g. Aqua purificata). Please dry all components thoroughly with a lint-free disposable cloth. For the cleaning of the internal irrigation tube oil-free compressed air is mandatory. Please re-check all parts for damage or corrosion afterwards.

4. DISINFECTION

We recommend a high level disinfectant such as, Cidex OPA (Johnson & Johnson) for disinfection of the medical device.

- a) Soak the medical device in the disinfectant solution for the required amount of time. See IFU of Cidex OPA.
- b) Remove the medical device from the disinfectant solution.
- c) Rinse at least three times with highly purified water.
- d) Air dry and package the devices immediately.

5. STERILIZATION

If no sterilization device is available in the laboratory, this information should be forwarded to the dentist so proper sterilization can occur. Please use only validated sterilization procedures for the sterilization of the devices. Other sterilization procedures must not be used.

Reusability: You may only sterilize the medical device one time. In case of inadvertent contamination, you may re-sterilize one time after cleaning and disinfection.

6. STEAM STERILIZATION

- Fractionated vacuum procedure or gravity procedure (with sufficient product drying)
- Steam sterilizer according to ISO 17665: 2006 or EN 13060 and EN 285 respectively or equivalent national standards
- Validated according to EN ISO/ANSI AAMI 17665 (in past: EN 554/ANSI AAMI ISO 11134) (valid IQ/OQ) (commissioning and product specific performance qualification)
- Sterilization time not less than 15 minutes at 121°C (250°F) (listed exposure times are at sterilization temperature)

7. STORAGE

Store the sterilized parts dry and dust-free at room temperature (18-25°C / 64-77°F).

Application only by qualified personnel!

MANUFACTURER

nt-trading GmbH & Co. KG
G.-Braun-Str. 18
76187 Karlsruhe, Germany

Tel: +49-721-91547160
Fax: +49-721-91547161
E-mail: info@nt-trading.com

Products indicated with ® are registered brand names of the manufacturers

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