Scankörper 3D-Guide
Indikation

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Instructions for use [ENG]
Scan Body 3D-Guide und Implant Pick-Up
Scan Body 3D-Guide
Indication

The Scan Body 3D-Guide is used for demonstrating the position and alignment of dental tooth or implant analogs with intraoral and extraoral CAD/CAM scanning processes.

Composition

PEEK (Polyether-ether-ketone), stainless steel (1.4305), silicone.

Application

Scan Body 3D-Guides are fixed on the laboratory analog or compatible implant, which are shown on the label and which are applied in the dental working model or in the mandibles of the patient, by using the integrated fixing screw. In the individual scan processes, the scan body 3D guide shows the three-dimensional position of the individual implant or implant analog. This is necessary for the purpose of a correct 3D alignment for the following design process by using the CAD software (3Shape/Exocad/DentalWings).

Important information for Application

The Scan Body 3D-Guide is positioned in the pertaining dental implant or implant analog in the mandible of the patient or in the dental working model. Before inserting the scan body 3D guide, make sure that all components are clean and undamaged. In the case of intraoral application, the scan bodies must also be disinfected (according to professional standards) and protected against aspiration by using an appropriate tool (dental floss).

Check the correct position of the 3D guide in the tooth or implant analog and fix it hand-tight with the integrated fixing screw (max.10 Ncm).

Use the pertaining ISI shank screwdriver only for fixing the scan body 3D guide in the implant or in the implant analog.

Check the correct position of the scan body 3D guide in the tooth or implant analog again. If you feel any loosening in vertical or rotational direction, the scan body 3D guide is not correctly placed in the tooth or implant analog and must be positioned again and correctly.

Contraindication/Warnings

If you do not adhere to the course of actions described in this instruction for use, one or more of the following complications may occur:

- Injury to the patient, such as cross-contamination, aspiration of components, damage to or loss of the implant, etc.
- Erroneous scan result
- Incorrect final restoration or improper function of the crown, bridge or other prosthetic component
- Damage to the implant, secondary part, components or equipment

Only use scan body 3D guides with the compatible tooth or implant analog. If not, injury to patients may occur and/or implant, components or tools may be damaged.

The scan body 3D guide and all of its components are designed for single use. Multiple use of a scan body 3D-guide may cause inaccurate scanning results.

Implant Pick-Up
Open tray impression
Indication

High-precision impression components give a precise replica of the intraoral situation. Clear-cut tactile response from the prosthetic connection verifies proper seating of components.

Composition

TiAlAV, medical grade 5, ASTM F136, stainless steel (1.4305), silicone

Generals

The open tray impression post has the following characteristics for simplicity, reliability and application. Slender emergence profile accommodates space limitations. Guide screw can be tightened either by hand or with the screwdriver of the corresponding implant system.

Application

This open tray impression procedure requires a custom-made tray with perforations. Impression posts are intended for single use only to provide optimal fit and precise impression taking for each patient.

Open tray impression – Clinical procedure
Step 1 – Positioning the impression post

Ensure sufficient access to the implant site in order to avoid pinching in the gingival tissue. Be aware that the sulcus may collapse rapidly once the healing components have been removed.

Clean the internal configuration of the implant thoroughly from blood, tissue, etc. prior to the impression procedure.

Place the impression post accurately into the implant and hand-tighten the guide screw. In case of occlusal space limitation, the length of the impression post can be reduced by one retention ring after the knurled screw has been removed.

Step 2 – Impression taking

Adapt the custom-made impression tray according to the individual situation so that the positioning screw of the impression post sticks out.

Take the impression using an elastomerimpression material (polyvinyl siloxane or polyether-rubber). Uncover the screws before the material is cured.

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