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

Instructions for use [ENG] DIM Parallel Counterbore


Item name/Item description	Material	Item number
A DIM parallel counterbore	Stainless steel (1.4034, 1.4310)	W 16.DIM.000

COMPONENT DESCRIPTION

The DIM parallel counterbore is a one-piece tool with an integrated, spring-loaded guide pin.

Optional accessories:

- an ISO shaft insert made of stainless steel (1.4034) for insertion in a laboratory handpiece
- an insert (for an ISO shaft for the T, H, CAM and COL series, DIM analogue), 1.20 mm hex head, made of stainless steel (1.4034)
- and a screwdriver (including an ISO shaft for the T, H, CAM and COL series, DIM analogue), 1.20 mm hex head, made of stainless steel (1.4034)

Item No.	Name
1	DIM Parallel Counterbore
2	Spring-loaded guide pin
3	Internal springs
4	Safety cap
5	ISO shaft insert, laboratory handpiece

INDICATIONS

By using a DIM parallel counterbore, you can expand and deburr the opening to the knurled screw when using the DIM analogues; the positioning takes place via the guide pin of this tool.

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
GENERAL INFORMATION

The DIM parallel counterbore, or the optional components, can be used in the laboratory for the production of implant models. The use of interchangeable inserts allows loosening and tightening of the safety cap located at the end of the DIM parallel counterbore. In addition, there are DIM knurled screws and abutment screws for the respective implant systems listed (Item D). To tighten the screws, the torque recommended by the screw manufacturer must be strictly observed.

APPLICATIONS

DIM analogues are used for the production of implant models in the digital workflow and, using the DIM sleeve, also in the classic impression procedure. Both applications are described STEP BY STEP in the „DIM Analogue“ instructions and are supplemented by this detailed description of the DIM parallel counterbore, which is mainly used in classic model production. The knurled screw of the DIM analogue can be seen in the basal view of the manufactured implant model. If necessary, the opening to the knurled screw can be enlarged and deburred by means of the DIM parallel counterbore. To do this, the spring-loaded guide pin is inserted into the screwdriver holder of the knurled screw.

Guided by the pin, the parallel counterbore is pressed against the spring pressure in the opening to the knurled screw. By means of manual rotary movements, the opening is enlarged and deburred to such an extent that easy access to the knurled screw results without damaging it.

To use the optional ISO shaft insert for the laboratory handpiece, remove the safety cap on the end using a 1.20 mm hex head screwdriver. It should be noted that the opening, after unscrewing the cap, is permanently facing up to secure the internal springs. In the opening with internal threads, the threads for the laboratory handpiece insert can be screwed in. When using the torque ratchet for screwing in the insert, the insert can be connected to the ratchet via the hex head of the insert. To avoid damage and injury, the max. speed of 1500-2000/min for the laboratory handpiece should not be exceeded. In addition, the standard personal protective equipment used in a laboratory, e.g., safety glasses and gloves, is to be worn.

CONTRAINDICATIONS

According to current knowledge, there are no known contraindications.

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HAZARD AND SAFETY INSTRUCTIONS

The DIM parallel counterbore, or the individual components, may only be used by those with sufficient knowledge and qualifications in dental technology processing and dental applications (dentists, oral and maxillofacial surgeons, dental technicians, dental assistants, dental hygienists, dental assistance). All parts must be disinfected before use in accordance with hygiene guidelines. For disinfection, only instrument disinfectants may be used. The product may only be used as indicated. The torques recommended by the screw manufacturer must be strictly adhered to. Before use, check that the selected insert has the correct size and shape for matching the compatible implant system.

CLEANING, DISINFECTION AND STERILIZATION

The DIM parallel counterbore is supplied non-sterile in suitable packaging and must be cleaned and disinfected before use. (Failure to do so will void the warranty). Thorough cleaning and disinfection is an indispensable requirement for the instruments and tools that are used in the production of medical devices. In addition to the hygiene guidelines for dental practice, the statutory regulations that apply locally must be observed.

MANUFACTURER:

nt-trading GmbH & Co. KG	Tel: +49 - 721 - 91 54 71 60
G.-Braun-Straße 18	Fax: +49 - 721 - 91 54 71 61
76187 Karlsruhe, Germany	E-mail: info@nt-trading.com

Products indicated with * are registered brand names of the manufacturers.

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