INDIKATION

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Instructions for use: Dentokeep PEEK Disc [ENG]

INDICATION

nt-Dentokeep PEEK discs are blanks to be applied for the manufacturing of permanent dental prostheses by milling as well as for model casting technology in a CAM process. The prostheses made from these blanks can be fixed on remaining tooth stumps and/or on implant abutments and guarantee a fundamental esthetical and functional treatment.

CONTRAINDICATION

Based upon currently available information, there are no known contra-indications. If any incompatibility is suspected, this product may only be used after prior allergological clarification and verification that the product does not cause allergic reactions.

COMPOSITION

Polyetheretherketones (PEEK)

SAFETY INSTRUCTIONS

In order to guarantee safety and quality, the processor must adhere to all information and instructions included in these directions for use. During processing, the usual safety precautions (respirator, extraction) as well as personal protective equipment (protective goggles) must be used.

PROCESSING

- Minimum wall thickness: 0.6 mm
- Cross-sectional areas of connectors in the front teeth zone: 10 mm²
- Cross-sectional areas in the lateral teeth zone: 16 mm²

Do not position more than two links in the lateral teeth zone.

nt-Dentokeep discs must be processed by using the special nt-PEEK Burr carbide metal milling tools. The use of other milling tools may cause breaking of the milling tools and damage to the disc or CAM unit. Use a coarsely toothed carbide metal trepan to cut the milled prosthesis from the blank. Cut the holding frames without applying any pressure. Prepolishing of the nt-Dentokeep restorations can be carried out by using the appropriate silicone polishing tools and goat hair brushes.

FIXING

After fitting the prosthesis, a composite primer (Ambarino P60 Creamed, Visio.link Bredent, Z-Prime Plus Bisco) is applied to the inner surface of the crown, which has been sand-blasted with 220 mµ Al₂O₃. The primer establishes a firm connection between Nt-Dentokeep and the fixing composite.

VENEER

Suited materials for veneering like Gradia (from the company GC Europe) and Sinfony (from the company 3M Espe) and the primer.

AVAILABLE VERSIONS

Durchmesser (mm)	Dicke (mm)	Farbe	Artikel-Nr.
98.5	14	pearl white	12.000.P14
98.5	18	pearl white	12.000.P18
98.5	23	pearl white	12.000.P23
98.5	14	tooth-colored	12.00A.P14
98.5	18	tooth-colored	12.00A.P18
98.5	23	tooth-colored	12.00A.P23
98.5	14	gingiva-colored	12.GUM.P14
98.5	18	gingiva-colored	12.GUM.P18
98.5	23	gingiva-colored	12.GUM.P23

PLEASE NOTE

This information is a recommendation and shall be considered as guidelines. Changes in handling or composition reserved.

CLEANING, DISINFECTION AND STERILIZATION:

The listed products in this IFU are supplied in non-sterile condition. The components should be cleaned, disinfected and in specific clinical procedures and cases be sterilized, prior to use after they are received from the dental laboratory (no liability on disregard). Effective cleaning and disinfection is an indispensable requirement for effective sterilization of the medical device. Prior to sterilization, please keep medical devices clean when handling in the laboratory and operatory. Additionally, please pay attention to legal regulations valid for your local areas as well as to the hygienic instructions of your dental practice. This applies particularly to the different guidelines regarding the inactivation of prions. Please observe all manufacturers guidelines for disinfection and cleaning agents with special regard to the concentration, exposure time and temperature. Only pH-neutral disinfection solutions without chlorine, ammoniac and aldehydes and with a proven effectiveness against HBV, HCV, and HIV must be used. The products have to meet the respective national regulations for disinfectants. If disinfectants containing aldehydes are used, this might lead to a possible fixation of proteins. Please use only freshly prepared solutions.

1. Pre-Disinfection (Avoidance of cross Contaminations)

Place the devices in a germicidal bath immediately after use. Remove all residues and disassemble demountable products.

2. Cleaning

Please use distilled water and neutral cleaning agents only. The products must be cleaned with a plastic instrument cleaning brush and then rinsed with distilled water. Please control all products after cleaning in order to avoid either damaging or corrosion. Damaged products must be replaced.

3. Rinsing and Drying

After removal of the products from the germicidal bath, all components must be rinsed 3 times with distilled water (e.g. Aqua purificata). Please dry all components thoroughly with a lint-free disposable cloth. Please re-check all parts for damage or corrosion afterwards.

4. Disinfection

We recommend a high level disinfectant such as, Cidex OPA (Johnson & Johnson) for disinfection of the Dentokeep PEEK Disc
a) Soak the medical device in the disinfectant solution for the required amount of time. See instructions for use of Cidex OPA
b) Remove the medical device from the disinfectant solution
c) Rinse at least three times with highly purified water
d) Air dry and package the devices immediately

5. Sterilization:

If no sterilization device is available in the laboratory, this information should be forwarded to the dentist so proper sterilization can occur. Please use only validated sterilization procedures for the sterilization of the devices. Other sterilization procedures must not be used.

Reusability: You may only sterilize the medical device one time. In case of inadvertent contamination, you may re-sterilize one time after cleaning and disinfection.

6. Steam Sterilization:

- fractionated vacuum procedure or gravity procedure (with sufficient product drying)
- steam sterilizer according to ISO 17665: 2006 or EN 13060 and EN 285 respectively or equivalent national standards
- validated according to EN ISO / ANSI AAMI 17665 (in past: EN 554 / ANSI AAMI ISO 11134)
(valid IQ / OQ) (commissioning and product specific performance qualification)
- sterilization time not less than 15 minutes at 121 °C (250 °F) (listed exposure times are at sterilization temperature)

MILLING TOOLS

The carbide metal milling tool nt-PEEK burr is appropriate for the processing of PEEK materials.

nt-PEEK Burr

Ø	1.00 mm	2.50 mm
Item no.	11301	11302

STORAGE

Dentokeep may only be stored in original packaging dry and at normal room temperature (15°C - 25°C). Dentokeep must be protected against sunlight and UV-light.

⚠ Application only by qualified personnel!

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Products indicated with ® are registered brand names of the manufacturers.

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nt-trading GmbH & Co KG