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Processing Instructions of gold-plastic abutment

Cast-On Alloy with Residue-Free, Burn-Out Plastics Duct. During casting the non-oxidizable gold plastic abutment is wetted by the melt of a precious casting, universal or bonding alloy (Au / Pt min. 25%). The metallic compound is produced by diffusion of alloy components of both materials.

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Gold-plastic abutment of each series can only be combined with the matching implant. They cannot be combined with implants of a different implant type or manufacturer. The diameter of the Gold-plastic abutment must correspond in size to the used implant in order to prevent a peri-implant tissue irritation

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2. Scope of Application

For the production of customized build-ups made of precious metal alloys.

For the cast-on process, precious metal alloys type 4 are appropriate:
Fuse-on alloy standard ISO 22764 + 9693
Gold casting alloy standard ISO 22764
Precious metal red. alloy standard ISO 22764
Base metal alloys must not be used!

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4. Technical Data

Melting range 1400 – 1490 °C
Vickers hardness HV5
(w) soft-annealed 1000 °C/30 min, w HV5 145
(s) self-hardening by slow cooling of the casting, s HV5 225
(a) Hardening/tempering 700 °C/30 min, cooling in air, a HV5 240

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5. Processing

Modeling: The length of the residue-free burn-out plastics sleeve can be reduced to meet the anatomic requirements. Modeling is carried out according to the established rules of dental technology. The internal connection and the implant support must be absolutely free of: Plastics, wax and grease, to avoid a cast-on process in this area.

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6. Investing

Only phosphate-bonded (gypsum-free) investment compounds, which are recommended for precious metal casting, can be applied. In particular, much attention must be paid to bubble-free investing.

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7. Preheating

The specifications given by the manufacturer of the applied investment compound must be adhered to. The preheating temperatures specified by the alloy manufacturers must be increased by approx. 50°C, but a temperature of at least 750°C is required for accurate casting. The casting plates should be preheated for a longer period.

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Nicht steril

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8. Casting

The melting range of the alloy to be cast on must be considerably below the solidus temperature of the gold-plastics build-up to avoid partial melting of the cast-on alloy.

With centrifugal casting, the object should be placed in the casting cylinder contrary to the direction of rotation to guarantee accurate casting

9. Cooling

Cool the casting cylinder to room temperature to avoid stresses in the casting object.

10. Divesting/Cleaning

Divesting must be carried out with much care. The casting objects are blasted with high-brilliance blasting beads at low pressure (max. 2 bar). The internal connection must not be blasted to avoid the risk of endangering the high precision. For complete cleaning, the casting object can be placed into a hot mordanting solution or can be cleaned by steam blasting.

11. Elaboration

The alloy is elaborated in compliance with the recommendations given by the precious metal manufacturers. During the elaborating process, it must be made sure under all circumstances that the cast-on alloy is not uncovered in the veneered area. The cast-on alloy must have a minimum wall thickness of 0.3 mm in this area. For elaborating, fixing of the gold-plastics build-up in a laboratory implant is recommendable.

12. Oxidizing

Oxidizing is carried out according to the specifications given by the manufacturer of the alloy.

13. Veneering

The preparation of the veneering area and the veneering process (ceramics/plastics) are carried out according to the specifications given by the manufacturers. In case of ceramic veneering, the cast-on alloy must not be covered, since this may cause cracks or spalling due to the different thermal expansions coefficients of the alloys.

14. Polishing

For polishing, fixing of the gold-plastics build-up in a laboratory implant is recommendable.

15. Recommended torques for the abutment screw:

<u>Ncm</u>	<u>Abutment</u>
20	H-Serie
25	S-, T-Series, F-Series NP
30	R-Series
35	I-, K-, L-, E-, N-Series, F-Series RP

16. Cleaning, Disinfection and Sterilization

The abutments and abutment screws of the series E, I, K, N, R, S, T, H, L and F are supplied in non-sterile condition.

The abutments and abutment screws must be cleaned and disinfected according to professional standards before and after each use one a patient.

17. Explanations of symbols

Do not reuse

Lot number

REF

Article number

Manufacturer

Consult instructions for use

Non-sterile

Usable up to