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
Securing and assembling the nt-Preform® milling blank in the nt-trading preform basic holder, using the example of the VHF S1/S2 dental milling machine

1) Holder components

- Basic holder

The basic holder (W12.Vhf.000) is designed in such a way that it fits precisely into the circular blank holder of the corresponding milling machine and therefore ensures perfect and defined support during the milling operation in the machine. Using abutment grips, a maximum of six nt-Preform® blanks can be in clamped in one operation.

- Abutment grips

Three different abutment grips (yellow Ø3.5 mm/W12.100.000, red Ø4.3 mm/W12.200.000 and blue Ø5.0 mm/W12.300.000) are used to hold the milling analogs required for the chosen implant system. The nt-Preform® is fixed in a defined position where the implant connection is optimally protected during the milling operation.

You can find an overview of the classification for this purpose in our nt-trading product catalog.

- Securing frame

The securing frame (W12.F11.005) ensures that the nt-Preform® is positioned in the holder securely and without any vibration during the machining operation.

2) Milling components

- Milling blank (nt-Preform®)*

The industrially produced, high-precision implant connection and the prefabricated screw channel make it possible to manufacture individual titanium abutments without compromising on quality. The abutment that you have designed will only be milled in the cylindrical area of the nt-Preform®, meaning that critical areas of the implant connection or the screw channel do not have to be machined.

- Milling screw*

This item is used to screw the nt-Preform® to the corresponding milling analog. The milling screw may also be used as the working screw in the cast model. **This milling screw must not be used to integrate the abutment when treating a patient!**

- Milling analog*

The milling analog secures the nt-Preform® in the relevant abutment grip assigned to the implant system and implant diameter. It also protects the implant connection during the milling operation.

*In our nt-trading catalog, you can find the classification of the item numbers according to the implant systems and diameters

3) Tools

Every tool is clearly assigned its screw by the different head profiles and thread diameters in the holder. This ensures that there is no risk of confusion and the components can be assembled intuitively.

- Tool A (W12.WS0.001)

For securing the milling analog to the abutment grip

- Tool B (W12.SD0.000)

For securing the nt-Preform® to the milling analog

- Tool C (W12.WE1.000)

For securing the abutment grip to the holder

- Tool D (W12.WE2.000)

For securing the securing frame to the holder

4) Fitting the basic holder

- Visual inspection

Before assembling the holder, you should ensure that none of the required components are damaged or contaminated (e.g. by shavings from milling) so that the components can be precisely aligned with one another.


- Fitting the abutment grip

Depending on the implant system and diameter, the right color-coded abutment grip (see catalog for classification) should be fitted to the associated milling analog. The groove in the bottom section of the analog must correspond to the pin in the abutment grip.

By gently rotating, it is easy to determine the correct position and the milling analog should audibly click into the final position.


The associated set screw is screwed in at a max. of 60 Ncm (tool A/W12.WS0.001) and secures the milling analog in its distinct position.

The set screws which are in the positions on the abutment grip not occupied by the milling analogs should be screwed in slightly to avoid being lost during the milling operation. Under no circumstances should these screws damage the pins in the abutment grip for locking the milling analogs.


Warning:


The flat side on the upper section of the nt-Preform® must be aligned to the left of the grip!

The nt-Preform® to be milled, along with its associated milling screw, can now be secured finger-tight to the analog (tool B/W12.WS0.000).

Screwing in set screws A and B makes it easier to uncouple the abutment grip from the holder after the milling process (it may be held in place due to shavings from milling, cooling lubricant, etc.).

Before each milling operation, you should therefore ensure that these screws are loosened again after being screwed in and that they do not protrude from the rear side of the abutment grip.

- Securing the abutment grip

The abutment grip is secured to the holder finger-tight by the two screws provided for this purpose (tool C/W12.WE1.000). The upper section of the nt-Preform® that has already been fitted is clearly positioned in the central bar of the holder, with the flat area aligned to the left.

- Securing the securing frame

The securing frame is placed over the screws that are already fitted in the holder's central bar and is moved laterally. The screws are then tightened finger-tight (tool D/W12.WE2.000).

5) Replacing components

- Screws

A screw set (W12.SC0.000) can be reordered separately if necessary.

- Cover strip

The cover strip (W12.CS0.000) can also be reordered if it needs to be replaced. To replace the cover strip, the screws underneath the abutment grip must be undone (tool A/W12.WS0.001) and the strip must be removed and replaced with a new one. Ensure that the pins on the cover strip are aligned to the rear of the abutment grip.