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INSTRUCTIONS FOR USE [ENG] SCREWDRIWER-SET

ARTICLE NUMBERS

W11.000.SET • W11.LN2.COF (yellow) • W11.SEV.SMU (orange) • W11.IRS.000 (black) • W11.EFK.000 (green) • W11.2CO.NMO (violet)
W11.THCDIM (white) • W11.MEX.000 (grey) • W11.VAL.BAL (brown)
W11.Y00.000 (red) • W11.NEO.000 (magenta) • W11.THO.000 (light blue)

DESCRIPTION OF THE COMPONENTS

The screwdriver set consists of the following:

- A case made of Aluminium EN AW 6082 [AlSi1MgMn] with a lockpart consisting of a sleeve (1.4303), a ball (1.3541) and a spring (1.4568) made of stainless steel.
- eleven screwdrivers made out of stainless steel (1.4034) with a composite ceramic color marking (see article numbers / CAS No. 72869-86-4)

INDICATIONS

The screwdriver set is used for attaching prosthetic abutments to dental implants by means of abutment screws.

GENERAL

The screwdrivers in the screwdriver set can be used in the laboratory on the implant model as well as orally in the mouth of a patient. The use of replaceable screwdrivers allows the insertion or removal of a removable dental prosthesis by the loosening and tightening of abutment screws, also with different implant systems. To tighten the screws, it is essential that the torque recommended by the screw manufacturer be applied.

- ! Abutment screws should to be retightened after approximately 15 minutes with the same torque to achieve the maximum screw pretension.

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time and temperature. Only pH-neutral disinfection solutions without chlorine, ammoniac and aldehydes and with a proven effectiveness against HBV, HCV, and HIV must be used. The products have to meet the respective national regulations for disinfectants. If disinfectants containing aldehydes are used, this might lead to a possible fixation of proteins. Please use only freshly prepared solutions.

1. PRE-DESINFEXION

(avoidance of cross contaminations)

Place the devices in a germicidal bath immediately after use. Remove all residues and disassemble demountable products.

2. CLEANING

Please use distilled water and neutral cleaning agents only. The products must be cleaned with a plastic instrument cleaning brush and then rinsed with distilled water. Please control all products after cleaning in order to avoid either damaging or corrosion. Damaged products must be replaced.

3. RINSING AND DRYING

After removal of the products from the germicidal bath, all components must be rinsed 3 times with distilled water (e.g. Aqua purificata). Please dry all components thoroughly with a lint-free disposable cloth. Please re-check all parts for damage or corrosion afterwards.

4. DESINFEXION

We recommend a high level disinfectant such as, Cidex OPA (Johnson & Johnson) for disinfection of the abutments and screws.

- Soak the medical device in the disinfectant solution for the required amount of time. See instructions for use of Cidex OPA.
- Remove the medical device from the disinfectant solution.
- Rinse at least three times with highly purified water.
- Air dry and package the devices immediately.

5. STERILIZATION

If no sterilization device is available in the laboratory, this information should be forwarded to the dentist so proper sterilization can occur. Please use only validated sterilization procedures for the sterilization of the devices. Other sterilization procedures must not be used.

6. STEAM

- Fractionated vacuum procedure or gravity procedure (with sufficient product drying)
- Steam sterilizer according to ISO 17665: 2006 or EN 13060 and EN 285 respectively or equivalent national standards
- Validated according to EN ISO/ANSI AAMI 17665 (in past: EN 554/ANSI AAMI ISO 11134) (valid IQ/OQ)

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(commissioning and product specific performance qualification)
• Sterilization time not less than 15 minutes at 121 °C (250 °F)
(listed exposure times are at sterilization temperature).

STORAGE

Store the sterilized parts dry and dust-free at room temperature.

 nt-trading GmbH & Co KG	Fon: +49 - 721 - 91 54 71 60
G.-Braun-Straße 18	Fax: +49 - 721 - 91 54 71 61
76187 Karlsruhe, Germany	E-mail: info@nt-trading.com

Products indicated with * are registered brand names of the manufacturers.

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