
## GEBRAUCHSANWEISUNG [DEU] SURGICAL GUIDE SLEEVE

### INDIKATION

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
## INSTRUCTION FOR USE [ENG] SURGICAL GUIDE SLEEVE

### INDICATION

nt-trading Surgical Guide Sleeves support surgical drilling parallel to the axis and are, therefore, indexed for the application in CAD/CAM technology. Surgical drilling guides (drill-sleeve and guide-sleeve), planned and created with CAD/CAM, are used for guided insertion of prosthetically and anatomically ideal placed implants.

#### Implant-Planing-System:

nt-trading Surgical Guide Sleeves serve software-based presurgical planning and support surgical implementing of parallel to the axis planned implant insertions. Surgical Guide Sleeves are used during installation in surgical drilling guides made by CAD/CAM technology for use with surgical drills.

The outer diameter of these surgical drills must match the inner diameter of the corresponding Surgical Guide Sleeve or the outer diameter of a separate drill-key system. This ensures guidance and support of the drill. Surgical Guide Sleeves are used in the first steps for pilot drilling.

### GENERAL INDICATION AND ARTICEL NO. OF SURGICAL GUIDE SLEEVE

Surgical guide sleeves (SGS) are drill guide sleeves and made from titanium with different diameters in the same length.

Article-No. Guide Sleeve	Implant system	d	D	DD	L	SW
SGS 16.SD4.148	NobelBiocare NP	Ø 4,10 mm	4,80 mm	5,80 mm	4,00 mm	4,50 mm
SGS 16.SD5.060	NobelBiocare RP	Ø 5,00 mm	6,00 mm	7,00 mm	4,00 mm	5,70 mm
SGS 16.SD6.270	NobelBiocare WP	Ø 6,20 mm	7,00 mm	8,00 mm	4,00 mm	6,70 mm
SGS 16.SD3.548	Camlog	Ø 3,50 mm	4,80 mm	5,80 mm	3,00 mm	4,50 mm
SGS 16.SD4.563	Camlog	Ø 4,50 mm	6,00 mm	7,00 mm	3,00 mm	5,70 mm
SGS 16.SD4.560	Ankylos/Xive ND	Ø 4,48 mm	6,00 mm	7,00 mm	4,00 mm	5,70 mm
SGS 16.SD4.960	Ankylos WD	Ø 4,90 mm	6,00 mm	7,00 mm	4,00 mm	5,70 mm
SGS 16.SD5.260	SEV /Xive WD	Ø 5,20 mm	6,00 mm	7,00 mm	4,00 mm	5,70 mm
SGS 16.SD4.660	SEV ND	Ø 4,60 mm	6,00 mm	7,00 mm	4,00 mm	5,70 mm
SGS 16.SD5.060	Straumann	Ø 5,00 mm	6,00 mm	7,00 mm	4,00 mm	5,70 mm
SGS 16.SD4.760	Astra	Ø 4,70 mm	6,00 mm	7,00 mm	4,00 mm	5,70 mm
SGS 16.SD5.770	Astra	Ø 5,70 mm	7,00 mm	8,00 mm	4,00 mm	6,70 mm
Article-No. Individual Sleeve	Implant system	d	D	DD	L	SW
SGS 16.PD2.032	pilot drill-2,0 mm	2,00 mm	3,20 mm	4,20 mm	4,00 mm	2,90 mm
SGS 16.PD2.232	pilot drill-2,2 mm	2,20 mm	3,20 mm	4,20 mm	4,00 mm	2,90 mm
SGS 16.PD2.332	pilot drill-2,3 mm	2,35 mm	3,20 mm	4,20 mm	4,00 mm	2,90 mm
SGS 16.PD2.848	pilot drill-2,8 mm	2,80 mm	4,80 mm	5,80 mm	4,00 mm	4,50 mm

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### GARANTIE

nt-trading GmbH & Co. KG gives a 10-year warranty on the mechanical stability of the Surgical Guide Sleeve with professional treatment in compliance with our instruction for use. Irrespective of that, are all informations mediated verbally, written or in practical tutorials based on tests and our experience and can therefore only be seen as standard values. Our products are subject to continuous further development. In this context, we reserve for ourself, to make product modifications or changes relating to the construction and composition.

### CLEANING, DESINFECTION AND STERILIZATION

Surgical Guide Sleeves are supplied in non-sterile condition and should be cleaned, disinfected and in specific clinical procedures and cases be sterilized. No liability on disregard. Effective cleaning and disinfection is an indispensable requirement for effective sterilization of the medical device. Prior to sterilization, please keep medical devices clean when handling in the laboratory and operatory. Additionally, please pay attention to legal regulations valid for your local areas as well as to the hygienic instructions of your dental practice. This applies particularly to the different guidelines regarding the inactivation of prions. Please observe all manufacturers guidelines for disinfection and cleaning agents with special regard to the concentration, exposure time and temperature. Only pH-neutral disinfection solutions without chlorine, ammoniac and aldehydes and with a proven effectiveness against HBV, HCV, and HIV must be used. The products have to meet the respective national regulations for disinfectants. If disinfectants containing aldehydes are used, this might lead to a possible fixation of proteins. Please use only freshly prepared solutions.

#### 1. PRE-DESINFECTION (AVOIDANCE OF CROSS CONTAMINATIONS)

Place the devices in a germicidal bath immediately after use. Remove all residues.

#### 2. CLEANING

Please use distilled water and neutral cleaning agents only. The products must be cleaned with a plastic instrument cleaning brush and then rinsed with distilled water. Please control all products after cleaning in order to avoid either damaging or corrosion. Damaged products must be replaced.

#### 3. RISING AND DRYING

After removal of the products from the germicidal bath, all components must be rinsed 3 times with distilled water (e.g. Aqua purificata). Please dry all components thoroughly with a lint-free disposable cloth. Please re-check all parts for damage or corrosion afterwards.

#### 4. DISINFECTION

We recommend a high level disinfectant such as, Cidex OPA (Johnson & Johnson) for disinfection of the Surgical Guide Sleeves.

- Soak the medical device in the disinfectant solution for the required amount of time. See instructions for use of Cidex OPA
- Remove the medical device from the disinfectant solution
- Rinse at least three times with highly purified water
- Air dry and package the devices immediately

#### 5. STERILIZATION

Please use only validated sterilization procedures for the sterilization of the devices. Other sterilization procedures must not be used. The process described refers exclusively to the Surgical Guide Sleeves. The packaging is not suitable for sterilization. For Cleaning, Desinfection and Sterilisation instructions of the stencil material please contact the respective manufacturer.

You may only sterilize the medical device one time. In case of inadvertent contamination, you may re-sterilize one time after cleaning and disinfection.

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### 6. STEAM STERILIZATION

- Fractionated vacuum procedure or gravity procedure (with sufficient product drying)
- Steam sterilizer according to ISO 17665: 2006 or EN 13060 and EN 285 respectively or equivalent national standards
- Validated according to EN ISO/ANSI AAMI 17665 (in past: EN 554/ANSI AAMI ISO 11134) (valid IQ/OQ) (commissioning and product specific performance qualification)
- Sterilization time not less than 15 minutes at 121 °C (250 °F) (listed exposure times are at sterilization temperature).

### STORAGE

Store the sterilized party dry and dust-free at room temperature. Protected against direct sunlight or light irradiation.

nt-trading GmbH & Co KG	Fon: +49 - 721 - 91 54 71 60
G.-Braun-Straße 18	Fax: +49 - 721 - 91 54 71 61
76187 Karlsruhe, Germany	E-mail: info@nt-trading.com

Products indicated with \* are registered brand names of the manufacturers.

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