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Instructions for use of laboratory analogs [EN]

COMPONENT DESCRIPTION

In the case of NT-trading implant analogs, a distinction is made between Lab Analog (laboratory analog) and NT-Mill Lab Analog (milling analog). Both variants are used for simulating and transferring the implant position of the patient into a dental working model. The NT-Mill Lab Analog can also be used in combination with the Abutment Grip for mounting NT-Preforms during the milling process.

The laboratory analog, in combination with the impression post, serves to transfer the position of the implant from the patient to a master dental model. For this purpose, the laboratory analog is screwed onto the particular compatible impression post in the impression and then model material (plaster, plastic) is poured into the impression. After the hardening of the model material and the releasing of the impression post and the impression, a three-dimensional dental model of the dental patient situation is created, in which the laboratory analog is anchored in a stable position by its retentive outer geometry and can now be used to accommodate various compatible abutment applications or prosthetic constructions.

COMPATIBILITY OF THE LABORATORY ANALOGS

The laboratory analogs and milling analogs are compatible with the particular implant systems shown in Table 1.

Assignment of laboratory analogs and milling analogs to the implant systems:

Serie (Laboranalogue / Fräsanaloge)	Kompatibles Implantatsystem
BEG	BEGO Implant Systems/Semados [®]
BIO	BioHorizons [®] / internal
BPS	BEGO Implant Systems/Semados [®] SC/SCX-/ RS/R SX-/RI-Line mit PS Design
BRE	Bredent Medical [®] / Sky [®]
CAM	Altatec / Camlog [®]
COL	Altatec / Conelog [®]
E	Nobel Biocare [™] / Replace Select [®]
EVMU	DENTSPLY Implants [®] / Astra Tech Implant System [™] EV / Multi Unit
F	Nobel Biocare [™] / Nobel Active [™]
H	Biomet 3i [®] / Certain [®]
HIO	HIOSEN [®] / OSSTEM TS [®]
I	Biomet 3i [®] / Osseotite [®] External Hex [®]
ISY	Altatec / ISY [®]
K	Nobel Biocare [™] / Brånemark [®]
L	Straumann / Bone Level [®]
LOC	Locator [®] Halteelement
MEG	MEGAGEN [®] / AnyRidge [®]
MEX	Medentis Medical [®] / ICX [®]
MIC	MIS [®] / C1 [®]
MIS	MIS [®] / SEVEN [®]
N	Straumann / SynOcta [®]
NBMU	Nobel Biocare [™] / Multi Unit
NEO	NEOSS Implant System [®] / NEOSS [™]
R	Zimmer Dental / Tapered Screw-Vent [®]
S	DENTSPLY Implants [®] / Astra Tech Implant System OsseoSpeed [®] TX
SEV	DENTSPLY Implants [®] / Astra Tech Implant System [™] EV
SMU	DENTSPLY Implants [®] / Astra Tech Implant System OsseoSpeed [®] TX / Multi Unit
STRMU	Straumann [®] / Multi Unit
T	DENTSPLY Implants [®] / Xive [®]
TA	Tatum Implant System
THO	Thommen Medical [®] / SPI [®]
UPR	Sweden & Martina Implantology [™] / Premium [®]
UTG	Sweden & Martina Implantology [™] / Premium Kohno [®] / Premium TG [®]
Y	DENTSPLY Implants [®] ANKYLOS [®]
2-CON	Nur in Anwendung mit 2-CONnect Abutment BEGO Implant Systems/ Semados [®] ; BioHorizons [®] / internal; Bredent Medical [®] / Sky [®] ; Altatec [®] / Camlog [®] ; Altatec / Conelog [®] ; Nobel Biocare [™] / Replace Select [®] ; Nobel Biocare [™] / Nobel Active [™] ; Biomet 3i [®] / Certain [®] ; HIOSEN [®] / OSSTEM TS [®] ; Biomet 3i [®] / Osseotite [®] External Hex [®] ; Nobel Biocare [™] / Brånemark [®] ; Straumann / Bone Level [®] ; MEGAGEN [®] / Medentis Medical [®] / ICX [®] ; MIS [®] / C1 [®] ; MIS [®] / SEVEN [®] ; Straumann / SynOcta [®] ; NEOSS Implant System [®] / NEOSS [™] ; Zimmer Dental / Tapered Screw-Vent [®] ; DENTSPLY Implants [®] / Astra Tech Implant System OsseoSpeed [®] TX; DENTSPLY Implants [®] / Astra Tech Implant System [™] EV; Sweden & Martina Implantology [™] / Premium Kohno [®] / Premium TG [®] ; DENTSPLY Implants [®] / Xive [®] ; Thommen Medical [®] / SPI [®]

MATERIAL

The laboratory analogs and milling analogs are made of stainless steel (1.4305).

INDICATIONS

Laboratory analogs and milling analogs are used to simulate and transfer the implant position of the patient into a dental working model.

Milling analogs are also used in combination with the Abutment Grip for mounting NT-Preforms during the milling process.

ACCESSORIES

Milling analog Abutment Grip (1.4305) is a holder for the individual milling of NT-Preform abutments.

APPLICATIONS

- **Laboratory analogs**

To simulate and transfer the implant position of the patient, the laboratory analog is fixed on the correspondingly compatible impression posts in the impression and this is then filled with a suitable model material (plaster, plastic). After releasing the impression, the laboratory analog remains in the dental model.

- **Milling analogs**

The milling analog, in combination with the Abutment Grip, is used for the individual milling of an NT-Preform abutment. The milling analog is inserted into the analog holders of the Abutment Grip and fastened with the screws present in the Abutment Grip. The compatible NT-Preform is then fixed on the milling analog with a milling screw and machined in the milling holder.

The use as an aid for simulating and transferring the patient situation into a dental model corresponds to the sequence of the laboratory analog described above.

APPLICATION NOTES

- Laboratory analogs and milling analogs are designed for a single use. Repeated use increases the risk of mechanical failure due to material fatigue.
- Before using a laboratory analog or milling analog, check whether the selected analog is compatible with the implant system.
- Before use, check the components for contamination and damage.
- Milling analogs: The connection component of the NT-Preform abutment must be clamped with the flattened side always on the left on the Abutment Grip, so that an optimal hold is ensured during the milling process.
- The milling analog must always be inserted with the groove down in the Abutment Grip in order to occupy its final position in the Abutment Grip.

STORAGE

The parts must be stored dry and dust-free at room temperature.

nt-trading GmbH & Co KG
G.-Braun-Str. 18
76187 Karlsruhe, Germany
Tel: +49-721-91547160
Fax: +49-721-91547161
E-mail: info@nt-trading.com

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