


DENTOKEEP PEEK DISC

Dentokeep PEEK discs for milling fabrication of dentures in CAM process



BENEFITS:

- biocompatible
- MMA and benzoyl peroxide-free
- sterilizable
- high fracture strength
- low thermo conductivity
- low plaque affinity
- for composite veneering
- high abrasion resistance



APPLICATION:

- monolithic crowns (alternative to metal anatomical crowns)
- implant fixed structures (bars, abutments, bridges combined with titanium bases) - dental telescopes, bar riders
- individual gingiva forming
- tabletops
- temporary prosthetics
- crowns / bridges (for composite veneering)
- occlusal splints



TECHNICAL DETAILS:

- tensile modulus of elasticity: ≥ 3800 MPa
- density: 1,30 – 1,50 g/cm³
- color: pearl white, A, GUM
- melting temperature: 343 °C
- 14 mm, 18 mm, 23 mm thickness
- diameter 98,5 mm