



NT

nt-Preform[®]

Milling blank for individual one-piece titanium abutments in CAM process



BENEFITS:

- precise prefabricated implant interface
- optimal gingiva forming by two different gingiva heights
- individual abutment design by prosthetical terms (diameter, etc.)
- optimal interface protection by bilateral fixing in preform holder during milling process



APPLICATION:

- primary telescopes on implants
- mesio-structures for crowns and bridges
- occlusal screwable structures



TECHNICAL DETAILS:

- Titanium Grade 5 ELI

