


UTG-SERIES

compatible to:

*Sweden & Martina
Implantology™
Premium TG®*



PLATFORMS:



3,3 - 4,25 mm



5,0 mm



TITANIUM BASE

incl. screw

UTG 80334200



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

UTG 9.S3D3.342



NT-MILL LAB ANALOG

UTG 5.La3.342




IMPLANT PICK-UP

for open-tray with integrated screw

UTG TR-00013.3



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

UTG 603342



DIM ANALOG®

UTG 5.DIM.334













SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A20





<p>NT-PREFORM® incl. coloured milling screw</p>	<p>UTG 10.Mi3.342.05</p>	
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>	<p>UTG 9.S3D3.342</p>	
<p>NT-MILL LAB ANALOG for NT-Preform®</p>	<p>UTG 5.La3.342</p>	
<p>COLOURED MILLING SCREW for NT-Preform®</p>	<p>UTG 6.Sc3.342</p>	
<p>ABUTMENT SCREW</p>  <p>Overview of all screwdrivers on pages 344 and 345</p>	<p>UTG 603342</p>	
<p>DIM ANALOG®</p>	<p>UTG 5.DIM.334</p>	
<p>SLEEVE (belongs to DIM ANALOG®) plaster model</p>	<p>DIM 5.XSL.A20</p>	
<p>ABUTMENT GRIP for NT-Holder fixture</p>	<p>W 12.200.000</p>	
<p>I-BUTMENT individual milled by NT-iDent</p> <p style="color: orange;">Not available in the USA</p>	<p>TZ15.TiO.Si1</p>	



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm UTG-2CON 8033421S

GH 2,5 mm UTG-2CON 8033422S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm UTG-2CON 8033421M

GH 2,5 mm UTG-2CON 8033422M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 20 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800






2-CONNECT®-ABUTMENT SET (ANGLED)

incl. 2-CONnect®-Coping, mounting aid,
Coping-screw and Male Part Screw, angled 20°,
Gingiva Height 2,5 mm


UTG-2CON 80334A3S



2-CONNECT®-ABUTMENT

angled 20°, Gingiva Height 2,5 mm
Male Part inclusive Screw,
mounting aid

UTG-2CON 80334A3M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning
KS 62 = 1,6 mm

2CON 9.S3D1.600



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 15 Ncm

KS 62



MALE PART SCREW

for angled 20° 2-CONnect® Abutment
Recommended tightening torques: 20 Ncm

UTG 6.Ms3.342



2-CONNECT®-DIM ANALOG®

KS 62 = 1,6 mm

2CON 5.DIM.160



SLEEVE (belongs to DIM ANALOG®)
plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw,
KS 62 = 1,6 mm

2CON TR-00031.6



2-CONNECT®-LAB ANALOG

KS 62 = 1,6 mm

2CON 5.La1.600





TITANIUM BASE

incl. screw

UTG 80500000



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

UTG 9.S3D5.000



NT-MILL LAB ANALOG

UTG 5.La5.000




IMPLANT PICK-UP

for open-tray with integrated screw

UTG TR-00025.0



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

UTG 605000



DIM ANALOG®

UTG 5.DIM.500













SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30





<p>NT-PREFORM® incl. coloured milling screw</p>	<p>UTG 10.Mi5.000.05</p>	
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>	<p>UTG 9.S3D5.000</p>	
<p>NT-MILL LAB ANALOG for NT-Preform®</p>	<p>UTG 5.La5.000</p>	 <p>groove</p>
<p>COLOURED MILLING SCREW for NT-Preform®</p>	<p>UTG 6.Sc5.000</p>	
<p>ABUTMENT SCREW</p>  <p>Overview of all screwdrivers on pages 344 and 345</p>	<p>UTG 605000</p>	
<p>DIM ANALOG®</p>	<p>UTG 5.DIM.500</p>	
<p>SLEEVE (belongs to DIM ANALOG®) plaster model</p>	<p>DIM 5.XSL.A30</p>	
<p>ABUTMENT GRIP for NT-Holder fixture</p>	<p>W 12.300.000</p>	
<p>I-BUTMENT individual milled by NT-iDent</p> <p>Not available in the USA</p>	<p>TZ15.TiO.Si1</p>	



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm UTG-2CON 8050001S

GH 2,5 mm UTG-2CON 8050002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm UTG-2CON 8050001M

GH 2,5 mm UTG-2CON 8050002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 20 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800






2-CONNECT®-ABUTMENT SET (ANGLED)

incl. 2-CONnect®-Coping, mounting aid,
Coping-screw and Male Part Screw, angled 20°,
Gingiva Height 2,5 mm

UTG-2CON 80500A3S



2-CONNECT®-ABUTMENT

angled 20°, Gingiva Height 2,5 mm
Male Part inclusive Screw,
mounting aid

UTG-2CON 80500A3M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning
KS 62 = 1,6 mm

2CON 9.S3D1.600



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 15 Ncm

KS 62



MALE PART SCREW

for angled 20° 2-CONnect® Abutment
Recommended tightening torques: 20 Ncm

UTG 6.Ms5.000



2-CONNECT®-DIM ANALOG®

KS 62 = 1,6 mm

2CON 5.DIM.160



SLEEVE (belongs to DIM ANALOG®)
plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw,
KS 62 = 1,6 mm

2CON TR-00031.6



2-CONNECT®-LAB ANALOG

KS 62 = 1,6 mm

2CON 5.La1.600

