

# BEG-SERIES

compatible to: BEGO Implant Systems® / Semados® S, RS-RX, RI-Line



**PLATFORM**

3,25 / 3,75 mm

4,1 mm

4,5 mm

5,5 mm

**TITANIUM BASE**

for individual milled  
Zirconium Abutment  
incl. screw  
Titanium Grade 5 ELI



BEG 80323700

BEG 80410000

BEG 80450000

BEG 80550000

**SCAN BODY 3D-GUIDE®**

for intra + extra  
oral scanning



BEG 9.S3D3.237

BEG 9.S3D4.100

BEG 9.S3D4.500

BEG 9.S3D5.500

★

**NT-MILL LAB ANALOG**

Stainless Steel



BEG 5.La3.237

BEG 5.La4.100

BEG 5.La4.500

BEG 5.La5.500

★

**IMPLANT PICK-UP**

for open-tray with  
integrated screw  
Coming soon



BEG TR-00013.2

BEG TR-00024.1

BEG TR-00034.5

BEG TR-00045.5

**ABUTMENT SCREW**

Recommended tightening  
torques 30 Ncm  
Titanium Grade 5 ELI



BEG 603255

BEG 603255

BEG 603255

BEG 603255

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# BEG-SERIES

compatible to: BEGO Implant Systems® / Semados® S, RS-RX, RI-Line



**PLATFORM**

3,25 / 3,75 mm

4,1 mm

4,5 mm

5,5 mm

**2-CONNECT®-ABUTMENT SET**

incl. 2-CONnect®-Coping and Coping-screw  
Titanium Grade 5 ELI

Coming soon



Gingiva Height	0,5 mm	BEG-2CON 8032371S	BEG-2CON 8041001S	BEG-2CON 8045001S	BEG-2CON 8055001S
	2,5 mm	BEG-2CON 8032372S	BEG-2CON 8041002S	BEG-2CON 8045002S	BEG-2CON 8055002S

**2-CONNECT®-ABUTMENT**

Recommended tightening torques: 30 Ncm  
Titanium Grade 5 ELI

Coming soon

M (Male Part)



Gingiva Height	0,5 mm	BEG-2CON 8032371M	BEG-2CON 8041001M	BEG-2CON 8045001M	BEG-2CON 8055001M
	2,5 mm	BEG-2CON 8032372M	BEG-2CON 8041002M	BEG-2CON 8045002M	BEG-2CON 8055002M

**2-CONNECT®-COPING**

Titanium Grade 5 ELI



F (Female Part)

2CON F-U	2CON F-U	2CON F-U	2CON F-U
----------	----------	----------	----------

★

**COPING-SCREW**

Recommended tightening torques: 25 Ncm  
Titanium Grade 5 ELI



KS 61	KS 61	KS 61	KS 61
-------	-------	-------	-------

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures. Products indicated with ★ are available in the USA

# BEG-SERIES

compatible to: BEGO Implant Systems® / Semados® S, RS-RX, RI-Line



**PLATFORM**

3,25 / 3,75 mm

4,1 mm

4,5 mm

5,5 mm

**2-CONNECT®-ABUTMENT SET**

incl. 2-CONnect®-Coping,  
Coping-screw & Male Part Screw  
angled 20°  
Gingiva Height 2,5 mm  
Titanium Grade 5 ELI

Coming soon



BEG-2CON 80323A3S

BEG-2CON 80410A3S

BEG-2CON 80450A3S

BEG-2CON 80550A3S

**2-CONNECT®-ABUTMENT**

angled 20°  
Gingiva Height 2,5 mm  
Titanium Grade 5 ELI

Coming soon

M (Male Part)



BEG-2CON 80323A3M

BEG-2CON 80410A3M

BEG-2CON 80450A3M

BEG-2CON 80550A3M

**2-CONNECT®-COPING**

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

2CON F-U

2CON F-U

★

**COPING-SCREW**

Recommended tightening  
torques: 15 Ncm  
Titanium Grade 5 ELI

Coming soon



KS 62

KS 62

KS 62

KS 62

**MALE PART SCREW**

for angled 20° 2-CONnect®  
Abutment | Recommended  
tightening torques: 20 Ncm  
Titanium Grade 5 ELI

Coming soon



BEG 6.Ms3.255

BEG 6.Ms3.255

BEG 6.Ms3.255

BEG 6.Ms3.255

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# BEG-SERIES

compatible to: BEGO Implant Systems® / Semados® S, RS-RX, RI-Line



**PLATFORM**

3,25 / 3,75 mm

4,1 mm

4,5 mm

5,5 mm

**2-CONNECT®-  
IMPLANT PICK-UP**  
integrated screw



● Straight KS 61 =  $\varnothing$  1,8 mm

2CON TR-00021.8

2CON TR-00021.8

2CON TR-00021.8

2CON TR-00021.8

● Angled KS 62 =  $\varnothing$  1,6 mm

2CON TR-00031.6

2CON TR-00031.6

2CON TR-00031.6

2CON TR-00031.6

**2-CONNECT®-  
SCAN BODY 3D-GUIDE®**  
for intra + extra  
oral scanning



● Straight KS 61 =  $\varnothing$  1,8 mm

2CON 9.S3D1.800

2CON 9.S3D1.800

2CON 9.S3D1.800

2CON 9.S3D1.800

★

● Angled KS 62 =  $\varnothing$  1,6 mm

2CON 9.S3D1.600

2CON 9.S3D1.600

2CON 9.S3D1.600

2CON 9.S3D1.600

★

**2-CONNECT®-  
LAB ANALOG**



● Straight KS 61 =  $\varnothing$  1,8 mm

2CON 5.La1.800

2CON 5.La1.800

2CON 5.La1.800

2CON 5.La1.800

★

● Angled KS 62 =  $\varnothing$  1,6 mm

2CON 5.La1.600

2CON 5.La1.600

2CON 5.La1.600

2CON 5.La1.600

★

**2-CONNECT®-  
DIM ANALOG**



● Straight KS 61 =  $\varnothing$  1,8 mm

2CON 5.DIM.180

2CON 5.DIM.180

2CON 5.DIM.180

2CON 5.DIM.180

★

● Angled KS 62 =  $\varnothing$  1,6 mm

2CON 5.DIM.160

2CON 5.DIM.160


2CON 5.DIM.160

2CON 5.DIM.160

★

**SLEEVE**

for plaster model  
Stainless Steel



DIM 5.XSL.A30

DIM 5.XSL.A30

DIM 5.XSL.A30





























DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# BEG-SERIES

compatible to: BEGO Implant Systems® / Semados® S, RS-RX, RI-Line





				
<b>PLATFORM</b>	3,25 / 3,75 mm	4,1 mm	4,5 mm	5,5 mm
<b>NT-PREFORM®</b> incl. coloured milling screw Titanium Grade 5 ELI				
Gingiva Height 0,5 mm	BEG 10.Mi3.237.05	BEG 10.Mi4.100.05	BEG 10.Mi4.500.05	BEG 10.Mi5.500.05
<b>NT-MILL LAB ANALOG</b> for nt-Preform® Stainless Steel				
	BEG 5.La3.237	BEG 5.La4.100	BEG 5.La4.500	BEG 5.La5.500 ★
<b>ABUTMENT GRIP</b> for nt-Holder fixture Stainless Steel				
	W 12.100.000	W 12.200.000	W 12.300.000	W 12.300.000 ★
<b>SCAN BODY 3D-GUIDE®</b> for intra + extra oral scanning				
	BEG 9.S3D3.237	BEG 9.S3D4.100	BEG 9.S3D4.500	BEG 9.S3D5.500 ★
<b>COLOURED MILLING SCREW</b> for nt-Preform® Recommended tightening torques 10 Ncm				
	BEG 6.Sc3.255	BEG 6.Sc3.255	BEG 6.Sc3.255	BEG 6.Sc3.255 ★
<b>ABUTMENT SCREW</b> Recommended tightening torques 30 Ncm Titanium Grade 5 ELI				
	BEG 603255	BEG 603255	BEG 603255	BEG 603255

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures. Products indicated with ★ are available in the USA 

**NEW**

# BEG-SERIES DIM ANALOG

compatible to: BEGO Implant Systems® / Semados® S, RS-RX, RI-Line

					
<b>PLATFORM</b>	3,25 / 3,75 mm	4,1 mm	4,5 mm	5,5 mm	
<b>DIM ANALOG</b>					
	BEG 5.DIM.323	BEG 5.DIM.410	BEG 5.DIM.450	BEG 5.DIM.550	★
<b>SLEEVE</b> for plaster model Stainless Steel					
	DIM 5.XSL.A20	DIM 5.XSL.A20	DIM 5.XSL.A30	DIM 5.XSL.A30	★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA 