

# BIO-SERIES

compatible to: BioHorizons® / internal



**PLATFORM**

3,0 mm

3,5 mm

4,5 mm

5,7 mm

**TITANIUM BASE**

for individual milled  
Zirconium Abutment  
incl. screw  
Titanium Grade 5 ELI



BIO 80300000

BIO 80350000

BIO 80450000

BIO 80570000

**SCAN BODY 3D-GUIDE®**

for intra + extra  
oral scanning



BIO 9.S3D3.000

BIO 9.S3D3.500

BIO 9.S3D4.500

BIO 9.S3D5.700

★

**NT-MILL LAB ANALOG**

Stainless Steel



BIO 5.La3.000

BIO 5.La3.500

BIO 5.La4.500

BIO 5.La5.700

★

**IMPLANT PICK-UP**

for open-tray with  
integrated screw  
Coming soon



BIO TR-00013.0

BIO TR-00023.5

BIO TR-00034.5

BIO TR-00045.7

**ABUTMENT SCREW**

Recommended tightening  
torques 30 Ncm  
Titanium Grade 5 ELI



BIO 603057

BIO 603057

BIO 603057

BIO 603057

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# BIO-SERIES

compatible to: BioHorizons® / internal



**PLATFORM**

3,0 mm

3,5 mm

4,5 mm

5,7 mm

**2-CONNECT®-ABUTMENT SET**

incl. 2-CONNECT®-Coping and Coping-screw  
Titanium Grade 5 ELI  
Coming soon



Gingiva Height

0,5 mm  
2,5 mm

BIO-2CON 8030001S	BIO-2CON 8035001S	BIO-2CON 8045001S	BIO-2CON 8057001S
BIO-2CON 8030002S	BIO-2CON 8035002S	BIO-2CON 8045002S	BIO-2CON 8057002S

**2-CONNECT®-ABUTMENT**

Recommended tightening torques: 30 Ncm  
Titanium Grade 5 ELI  
Coming soon



M (Male Part)

Gingiva Height

0,5 mm  
2,5 mm

BIO-2CON 8030001M	BIO-2CON 8035001M	BIO-2CON 8045001M	BIO-2CON 8057001M
BIO-2CON 8030002M	BIO-2CON 8035002M	BIO-2CON 8045002M	BIO-2CON 8057002M

**2-CONNECT®-COPING**

Titanium Grade 5 ELI

F (Female Part)



2CON F-U	2CON F-U	2CON F-U	2CON F-U	★
----------	----------	----------	----------	---

**COPING-SCREW**

Recommended tightening torques: 25 Ncm  
Titanium Grade 5 ELI



KS 61	KS 61	KS 61	KS 61	★
-------	-------	-------	-------	---

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures. Products indicated with ★ are available in the USA

# BIO-SERIES

compatible to: BioHorizons® / internal



**PLATFORM**

3,0 mm

3,5 mm

4,5 mm

5,7 mm

**2-CONNECT®-ABUTMENT SET**

incl. 2-CONnect®-Coping,  
Coping-screw & Male Part Screw  
angled 20°  
Gingiva Height 2,5 mm  
Titanium Grade 5 ELI

Coming soon



BIO-2CON 80300A3S

BIO-2CON 80350A3S

BIO-2CON 80450A3S

BIO-2CON 80570A3S

**2-CONNECT®-ABUTMENT**

angled 20°  
Gingiva Height 2,5 mm  
Titanium Grade 5 ELI

Coming soon

M (Male Part)



BIO-2CON 80300A3M

BIO-2CON 80350A3M

BIO-2CON 80450A3M

BIO-2CON 80570A3M

**2-CONNECT®-COPING**

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

2CON F-U

2CON F-U

★

**COPING-SCREW**

Recommended tightening  
torques: 15 Ncm  
Titanium Grade 5 ELI

Coming soon



KS 62

KS 62

KS 62

KS 62

**MALE PART SCREW**

for angled 20° 2-CONnect®  
Abutment | Recommended  
tightening torques: 20 Ncm  
Titanium Grade 5 ELI

Coming soon



BIO 6.Ms3.057

BIO 6.Ms3.057

BIO 6.Ms3.057

BIO 6.Ms3.057

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# BIO-SERIES

compatible to: BioHorizons® / internal



**PLATFORM**

3,0 mm

3,5 mm

4,5 mm

5,7 mm

**2-CONNECT®-  
IMPLANT PICK-UP**  
integrated screw



● Straight KS 61 = $\varnothing$ 1,8 mm	2CON TR-00021.8	2CON TR-00021.8	2CON TR-00021.8	2CON TR-00021.8
● Angled KS 62 = $\varnothing$ 1,6 mm	2CON TR-00031.6	2CON TR-00031.6	2CON TR-00031.6	2CON TR-00031.6

**2-CONNECT®-  
SCAN BODY 3D-GUIDE®**  
for intra + extra  
oral scanning



● Straight KS 61 = $\varnothing$ 1,8 mm	2CON 9.S3D1.800	2CON 9.S3D1.800	2CON 9.S3D1.800	2CON 9.S3D1.800	★
● Angled KS 62 = $\varnothing$ 1,6 mm	2CON 9.S3D1.600	2CON 9.S3D1.600	2CON 9.S3D1.600	2CON 9.S3D1.600	★

**2-CONNECT®-  
LAB ANALOG**



● Straight KS 61 = $\varnothing$ 1,8 mm	2CON 5.La1.800	2CON 5.La1.800	2CON 5.La1.800	2CON 5.La1.800	★
● Angled KS 62 = $\varnothing$ 1,6 mm	2CON 5.La1.600	2CON 5.La1.600	2CON 5.La1.600	2CON 5.La1.600	★

**2-CONNECT®-  
DIM ANALOG**



● Straight KS 61 = $\varnothing$ 1,8 mm	2CON 5.DIM.180	2CON 5.DIM.180	2CON 5.DIM.180	2CON 5.DIM.180	★
● Angled KS 62 = $\varnothing$ 1,6 mm	2CON 5.DIM.160	2CON 5.DIM.160	2CON 5.DIM.160	2CON 5.DIM.160	★

**SLEEVE**  
for plaster model  
Stainless Steel



DIM 5.XSL.A30	DIM 5.XSL.A30	DIM 5.XSL.A30	DIM 5.XSL.A30	★
---------------	---------------	---------------	---------------	---

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# BIO-SERIES

compatible to: BioHorizons® / internal



**PLATFORM**

3,0 mm

3,5 mm

4,5 mm

5,7 mm

**NT-PREFORM®**

incl. coloured milling screw  
Titanium Grade 5 ELI



Gingiva Height 0,5 mm

BIO 10.Mi3.000.05

BIO 10.Mi3.500.05

BIO 10.Mi4.500.05

BIO 10.Mi5.700.05

**NT-MILL LAB ANALOG**

for nt-Preform®  
Stainless Steel



BIO 5.La3.000

BIO 5.La3.500

BIO 5.La4.500

BIO 5.La5.700

★

**ABUTMENT GRIP**

for nt-Holder fixture  
Stainless Steel



W 12.100.000

W 12.100.000

W 12.200.000

W 12.300.000

★

**SCAN BODY 3D-GUIDE®**

for intra + extra  
oral scanning



BIO 9.S3D3.000

BIO 9.S3D3.500

BIO 9.S3D4.500

BIO 9.S3D5.700

★

**COLOURED MILLING SCREW**

for nt-Preform®  
Recommended tightening  
torques 10 Ncm



BIO 6.Sc3.057

BIO 6.Sc3.057

BIO 6.Sc3.057

BIO 6.Sc3.057

★

**ABUTMENT SCREW**

Recommended  
tightening torques  
30 Ncm  
Titanium Grade 5 ELI



BIO 603057

BIO 603057

BIO 603057

BIO 603057

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# BIO-SERIES DIM ANALOG

compatible to: BioHorizons® / internal

**NEW**



**PLATFORM**

3,0 mm

3,5 mm

4,5 mm

5,7 mm

**DIM ANALOG**



BIO 5.DIM.300

BIO 5.DIM.350

BIO 5.DIM.450

BIO 5.DIM.570

★

**SLEEVE**

for plaster model  
Stainless Steel



DIM 5.XSL.A10

DIM 5.XSL.A10

DIM 5.XSL.A30

DIM 5.XSL.A60

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA 