

BRE-SERIES

compatible to: Bredent Medical® / Sky®



PLATFORM

3,5 mm

4,0 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screw
Titanium Grade 5 ELI



BRE 80350000

BRE 80400000

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning




BRE 9.S3D3.500

BRE 9.S3D4.000

★

NT-MILL LAB ANALOG

Stainless Steel



BRE 5.La3.500

BRE 5.La4.000

★

IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



BRE TR-00013.5

BRE TR-00024.0

ABUTMENT SCREW

Recommended tightening
torques 25 Ncm
Titanium Grade 5 ELI



BRE 603540

BRE 603540

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

BRE-SERIES

compatible to: Bredent Medical® / Sky®



PLATFORM

3,5 mm

4,0 mm

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping and Coping-screw

Titanium Grade 5 ELI

Coming soon



Gingiva Height : 0,5 mm
: 2,5 mm

BRE-2CON 8035001S

BRE-2CON 8040001S

BRE-2CON 8035002S

BRE-2CON 8040002S

2-CONNECT®-ABUTMENT

Recommended tightening torques: 25 Ncm

Titanium Grade 5 ELI

Coming soon

M (Male Part)



Gingiva Height : 0,5 mm
: 2,5 mm

BRE-2CON 8035001M

BRE-2CON 8040001M

BRE-2CON 8035002M

BRE-2CON 8040002M

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

★

COPING-SCREW

Recommended tightening torques: 20 Ncm

Titanium Grade 5 ELI



KS 61

KS 61

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures. Products indicated with ★ are available in the USA

BRE-SERIES

compatible to: *Bredent Medical® / Sky®*



PLATFORM

3,5 mm

4,0 mm

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping,
Coping-screw & Male Part Screw
angled 20°

Gingiva Height 2,5 mm

Titanium Grade 5 ELI

Coming soon



BRE-2CON 80350A3S

BRE-2CON 80400A3S

2-CONNECT®-ABUTMENT

angled 20°

Gingiva Height 2,5 mm

Titanium Grade 5 ELI

Coming soon

M (Male Part)



BRE-2CON 80350A3M

BRE-2CON 80400A3M

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

★

COPING-SCREW

Recommended tightening
torques: 15 Ncm

Titanium Grade 5 ELI

Coming soon



KS 62

KS 62

MALE PART SCREW

for angled 20° 2-CONnect®
Abutment | Recommended
tightening torques: 20 Ncm

Titanium Grade 5 ELI

Coming soon



BRE 6.Ms3.540

BRE 6.Ms3.540

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 🇺🇸

BRE-SERIES

compatible to: Bredent Medical® / Sky®



PLATFORM

3,5 mm

4,0 mm

**2-CONNECT®-
IMPLANT PICK-UP**
integrated screw



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON TR-00021.8
2CON TR-00031.6

2CON TR-00021.8
2CON TR-00031.6

**2-CONNECT®-
SCAN BODY 3D-GUIDE®**
for intra + extra
oral scanning



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON 9.S3D1.800
2CON 9.S3D1.600

2CON 9.S3D1.800
2CON 9.S3D1.600

★
★

**2-CONNECT®-
LAB ANALOG**



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON 5.La1.800
2CON 5.La1.600

2CON 5.La1.800
2CON 5.La1.600

★
★

**2-CONNECT®-
DIM ANALOG**



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON 5.DIM.180
2CON 5.DIM.160

2CON 5.DIM.180
2CON 5.DIM.160

★
★

SLEEVE
for plaster model
Stainless Steel



DIM 5.XSL.A30

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

BRE-SERIES

compatible to: Bredent Medical® / Sky®



PLATFORM

3,5 mm

4,0 mm

NT-PREFORM®

incl. coloured milling screw
Titanium Grade 5 ELI




Gingiva Height 0,5 mm

BRE 10.Mi3.500.05

BRE 10.Mi4.000.05

NT-MILL LAB ANALOG

for nt-Preform®
Stainless Steel



BRE 5.La3.500

BRE 5.La4.000

★

ABUTMENT GRIP

for nt-Holder fixture
Stainless Steel



W 12.200.000

W 12.200.000

★

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



BRE 9.S3D3.500

BRE 9.S3D4.000

★

COLOURED MILLING SCREW

for nt-Preform®
Recommended tightening
torques 10 Ncm



BRE 6.Sc3.540

BRE 6.Sc3.540

★

ABUTMENT SCREW

Recommended
tightening torques
25 Ncm
Titanium Grade 5 ELI



BRE 603540

BRE 603540

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

BRE-SERIES DIM ANALOG

compatible to: *Bredent Medical*[®] / *Sky*[®]



PLATFORM

3,5 mm

4,0 mm

DIM ANALOG



BRE 5.DIM.350

BRE 5.DIM.400

★

SLEEVE

for plaster model
Stainless Steel



DIM 5.XSL.A20

DIM 5.XSL.A20

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA