

CAM-SERIES

compatible to: Altatec® / Camlog®



PLATFORM

3,3 mm

3,8 mm

4,3 mm

5,0 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screws
Titanium Grade 5 ELI



CAM 80330000

CAM 80380000

CAM 80430000

CAM 80500000

★

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



CAM 9.S3D3.300

CAM 9.S3D3.800

CAM 9.S3D4.300

CAM 9.S3D5.060

★

NT-MILL LAB ANALOG

Stainless Steel



CAM 5.La3.300

CAM 5.La3.800

CAM 5.La4.300

CAM 5.La5.060

★

IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



CAM TR-00013.3

CAM TR-00023.8

CAM TR-00034.3

—

LAB SCREW

coloured screw



CAM 6233843

CAM 6233843

CAM 6233843

CAM 6250600

★

ABUTMENT SCREW

Recommended tightening
torques 20 Ncm
Titanium Grade 5 ELI



CAM 6033843

CAM 6033843

CAM 6033843

CAM 6150600

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

CAM-SERIES

compatible to: Altatec® / Camlog®



PLATFORM

3,3 mm

3,8 mm

4,3 mm

5,0 mm

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping
and Coping-screw
Titanium Grade 5 ELI



Gingiva Height : 0,5 mm
: 2,5 mm

CAM-2CON 8033001S	CAM-2CON 8038001S	CAM-2CON 8043001S	CAM-2CON 8050001S
CAM-2CON 8033002S	CAM-2CON 8038002S	CAM-2CON 8043002S	CAM-2CON 8050002S

soon

2-CONNECT®-ABUTMENT

Recommended tightening
torques:
Titanium Grade 5 ELI
M (Male Part)



Gingiva Height : 0,5 mm
: 2,5 mm

CAM-2CON 8033001M	CAM-2CON 8038001M	CAM-2CON 8043001M	CAM-2CON 8050001M
CAM-2CON 8033002M	CAM-2CON 8038002M	CAM-2CON 8043002M	CAM-2CON 8050002M

soon

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)

2CON F-U	2CON F-U	2CON F-U	2CON F-U
----------	----------	----------	----------

★

COPING-SCREW

Recommended tightening
torques:
Titanium Grade 5 ELI



KS 61	KS 61	KS 61	KS 61
-------	-------	-------	-------

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

CAM-SERIES

compatible to: Altatec® / Camlog®



PLATFORM

3,3 mm

3,8 mm

4,3 mm

5,0 mm

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping,
Coping-screw & Male Part Screw
angled 20°
Gingiva Height 2,5 mm
Titanium Grade 5 ELI

Coming soon



over edge	CAM-2CON 80330A3S	CAM-2CON 80380A3S	CAM-2CON 80430A3S	CAM-2CON 80500A3S
over surface	CAM-2CON 80330C3S	CAM-2CON 80380C3S	CAM-2CON 80430C3S	CAM-2CON 80500C3S

2-CONNECT®-ABUTMENT

angled 20°
Gingiva Height 2,5 mm
Titanium Grade 5 ELI

Coming soon

M (Male Part)



over edge	CAM-2CON 80330A3M	CAM-2CON 80380A3M	CAM-2CON 80430A3M	CAM-2CON 80500A3M
over surface	CAM-2CON 80330C3M	CAM-2CON 80380C3M	CAM-2CON 80430C3M	CAM-2CON 80500C3M

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

2CON F-U

2CON F-U

★

COPING-SCREW

Recommended tightening
torques: 15 Ncm
Titanium Grade 5 ELI



KS 62

KS 62

KS 62

KS 62

★

MALE PART SCREW

for angled 20° 2-CONnect®
Abutment | Recommended
tightening torques: 20 Ncm
Titanium Grade 5 ELI

Coming soon



CAM 6.Ms3.384

CAM 6.Ms3.384

CAM 6.Ms3.384

CAM 6.Ms5.060

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

CAM-SERIES

compatible to: Altatec® / Camlog®



PLATFORM

3,3 mm

3,8 mm

4,3 mm

5,0 mm

2-CONNECT®- IMPLANT PICK-UP

integrated screw



● Straight KS 61 = \varnothing 1,8 mm

2CON TR-00021.8

2CON TR-00021.8

2CON TR-00021.8

2CON TR-00021.8

● Angled KS 62 = \varnothing 1,6 mm

2CON TR-00031.6

2CON TR-00031.6

2CON TR-00031.6

2CON TR-00031.6

2-CONNECT®- SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



● Straight KS 61 = \varnothing 1,8 mm

2CON 9.S3D1.800

2CON 9.S3D1.800

2CON 9.S3D1.800

2CON 9.S3D1.800

★

● Angled KS 62 = \varnothing 1,6 mm

2CON 9.S3D1.600

2CON 9.S3D1.600

2CON 9.S3D1.600

2CON 9.S3D1.600

★

2-CONNECT®- LAB ANALOG

Stainless Steel



● Straight KS 61 = \varnothing 1,8 mm

2CON 5.La1.800

2CON 5.La1.800

2CON 5.La1.800

2CON 5.La1.800

★

● Angled KS 62 = \varnothing 1,6 mm

2CON 5.La1.600

2CON 5.La1.600

2CON 5.La1.600

2CON 5.La1.600

★

2-CONNECT®- DIM ANALOG



● Straight KS 61 = \varnothing 1,8 mm

2CON 5.DIM.180

2CON 5.DIM.180

2CON 5.DIM.180

2CON 5.DIM.180

★

● Angled KS 62 = \varnothing 1,6 mm

2CON 5.DIM.160

2CON 5.DIM.160

2CON 5.DIM.160

2CON 5.DIM.160

★

SLEEVE

plaster model
Stainless Steel



DIM 5.XSL.A30

DIM 5.XSL.A30

DIM 5.XSL.A30

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 

CAM-SERIES

compatible to: Altatec® / Camlog®



PLATFORM

3,3 mm

3,8 mm

4,3 mm

5,0 mm

NT-PREFORM®

incl. coloured milling screw
Titanium Grade 5 ELI



Gingiva Height 0,5 mm

CAM 10.Mi3.300.05

CAM 10.Mi3.800.05

CAM 10.Mi4.300.05

CAM 10.Mi5.000.05

★

NT-MILL LAB ANALOG

for nt-Preform®
Stainless Steel



CAM 5.La3.300

CAM 5.La3.800

CAM 5.La4.300

CAM 5.La5.060

★

ABUTMENT GRIP

for nt-Holder fixture
Stainless Steel



W 12.100.000

W 12.100.000

W 12.200.000

W 12.300.000

★

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



CAM 9.S3D3.300

CAM 9.S3D3.800

CAM 9.S3D4.300

CAM 9.S3D5.060

★

COLOURED MILLING SCREW

for nt-Preform®
Recommended tightening
torques 10 Ncm



CAM 6.Sc3.384

CAM 6.Sc3.384

CAM 6.Sc3.384

CAM 6.Sc5.060

★

ABUTMENT SCREW

Recommended tightening
torques 20 Ncm
Titanium Grade 5 ELI



CAM 6033843

CAM 6033843

CAM 6033843

CAM 6150600

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

CAM-SERIES DIM ANALOG

compatible to: Altatec® / Camlog®



PLATFORM

3,3 mm

3,8 mm

4,3 mm

5,0/6,0 mm

DIM ANALOG



CAM 5.DIM.330

CAM 5.DIM.380

CAM 5.DIM.430

CAM 5.DIM.506



SLEEVE

for plaster model
Stainless Steel



DIM 5.XSL.A10

DIM 5.XSL.A40

DIM 5.XSL.A40

DIM 5.XSL.A50



Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA