

COL-SERIES

compatible to: Altatec® / Conelog®



PLATFORM

3,3 mm

3,8 / 4,3 mm

5,0 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screws
Titanium Grade 5 ELI



Gingiva Height
0,8 mm
2,0 mm

COL 80330008

COL 80384308

COL 80500008

COL 80330020

COL 80384320

COL 80500020

SCAN BODY

for Titanium Base
incl. coloured lab screw
PEEK



COL 703300

COL 703843

COL 705000

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



COL 9.S3D3.300

COL 9.S3D3.843

COL 9.S3D5.000

★

NT-MILL LAB ANALOG

Stainless Steel



COL 5.La3.300

COL 5.La3.843

COL 5.La5.000

★

IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



COL TR-00013.3

COL TR-00023.8

—

LAB SCREW

coloured screw



COL 623343

COL 623343

COL 625000

ABUTMENT SCREW

Recommended tightening
torques 20 Ncm
Titanium Grade 5 ELI



COL 603343

COL 603343

COL 605000

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

COL-SERIES

compatible to: Altatec® / Conelog®



PLATFORM

3,3 mm

3,8 / 4,3 mm

5,0 mm

2-CONNECT®-ABUTMENT SET

incl. 2-CONNECT®-Coping and Coping-screw

Titanium Grade 5 ELI

Coming soon



Gingiva Height : 0,5 mm
: 2,5 mm

COL-2CON 8033001S

COL-2CON 8038431S

COL-2CON 8050001S

COL-2CON 8033002S

COL-2CON 8038432S

COL-2CON 8050002S

2-CONNECT®-ABUTMENT

Recommended tightening torques:

Titanium Grade 5 ELI

M (Male Part)



20 Ncm



30 Ncm



30 Ncm

Gingiva Height : 0,5 mm
: 2,5 mm

COL-2CON 8033001M

COL-2CON 8038431M

COL-2CON 8050001M

COL-2CON 8033002M

COL-2CON 8038432M

COL-2CON 8050002M

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

2CON F-U

★

COPING-SCREW

Recommended tightening torques:

Titanium Grade 5 ELI



15 Ncm



20 Ncm



20 Ncm

KS 61

KS 61

KS 61

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 

COL-SERIES

compatible to: Altatec® / Conelog®



PLATFORM

3,3 mm

3,8 / 4,3 mm

5,0 mm

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping,
Coping-screw & Male Part Screw
angled 20°
Gingiva Height 2,5 mm
Titanium Grade 5 ELI
Coming soon



over edge	COL-2CON 80330A3S	COL-2CON 80384A3S	COL-2CON 80500A3S
over surface	COL-2CON 80330C3S	COL-2CON 80384C3S	COL-2CON 80500C3S

2-CONNECT®-ABUTMENT

angled 20°
Gingiva Height 2,5 mm
Titanium Grade 5 ELI
Coming soon
M (Male Part)



over edge	COL-2CON 80330A3M	COL-2CON 80384A3M	COL-2CON 80500A3M
over surface	COL-2CON 80330C3M	COL-2CON 80384C3M	COL-2CON 80500C3M

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)	2CON F-U	2CON F-U	2CON F-U	★
-----------------	----------	----------	----------	---

COPING-SCREW

Recommended tightening
torques: 15 Ncm
Titanium Grade 5 ELI
Coming soon



KS 62	KS 62	KS 62
-------	-------	-------

MALE PART SCREW

for angled 20° 2-CONnect®
Abutment | Recommended
tightening torques: 20 Ncm
Titanium Grade 5 ELI
Coming soon



COL 6.Ms3.343	COL 6.Ms3.343	COL 6.Ms5.000
---------------	---------------	---------------

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 🇺🇸

COL-SERIES

compatible to: Altatec® / Conelog®



PLATFORM

3,3 mm

3,8 / 4,3 mm

5,0 mm

2-CONNECT®- IMPLANT PICK-UP

integrated screw



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON TR-00021.8

2CON TR-00021.8

2CON TR-00021.8

2CON TR-00031.6

2CON TR-00031.6

2CON TR-00031.6

2-CONNECT®- SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON 9.S3D1.800

2CON 9.S3D1.800

2CON 9.S3D1.800

★

2CON 9.S3D1.600

2CON 9.S3D1.600

2CON 9.S3D1.600

★

2-CONNECT®- LAB ANALOG

Stainless Steel



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON 5.La1.800

2CON 5.La1.800

2CON 5.La1.800

★

2CON 5.La1.600

2CON 5.La1.600

2CON 5.La1.600

★

2-CONNECT®- DIM ANALOG



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON 5.DIM.180

2CON 5.DIM.180

2CON 5.DIM.180

★

2CON 5.DIM.160

2CON 5.DIM.160

2CON 5.DIM.160

★

SLEEVE

plaster model
Stainless Steel



DIM 5.XSL.A30

DIM 5.XSL.A30

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

COL-SERIES

compatible to: Altatec® / Conelog®



3,3 mm



3,8 / 4,3 mm



5,0 mm

NT-PREFORM®

incl. coloured milling screw
Titanium Grade 5 ELI



Gingiva Height : 1,3 mm
: 0,5 mm

COL 10.Mi3.300

COL 10.Mi3.843

COL 10.Mi5.000

COL 10.Mi3.300.05

COL 10.Mi3.843.05

COL 10.Mi5.000.05

NT-MILL LAB ANALOG

for nt-Preform®
Stainless Steel



COL 5.La3.300

COL 5.La3.843

COL 5.La5.000

★

ABUTMENT GRIP

for nt-Holder fixture
Stainless Steel



W 12.100.000

W 12.200.000

W 12.300.000

★

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



COL 9.S3D3.300

COL 9.S3D3.843

COL 9.S3D5.000

★

COLOURED MILLING SCREW

for nt-Preform®
Recommended tightening
torques 10 Ncm



COL 6.Sc3.343

COL 6.Sc3.343

COL 6.Sc5.000

★

ABUTMENT SCREW

Recommended
tightening torques 20 Ncm
Titanium Grade 5 ELI



COL 603343

COL 603343

COL 605000

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

COL-SERIES DIM ANALOG

compatible to: *Altatec® / Conelog®*



PLATFORM

3,3 mm

3,8 / 4,3 mm

5,0 mm

DIM ANALOG



COL 5.DIM.330

COL 5.DIM.384

COL 5.DIM.500

★

SLEEVE

for plaster model
Stainless Steel



DIM 5.XSL.A10

DIM 5.XSL.A20

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 