

H-SERIES

compatible to: Biomet 3i® / Certain®



PLATFORM

3,4 mm NP

4,1 mm RP

5,0 mm WP

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screw
Titanium Grade 5 ELI



H 800

H 810

H 820

★

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



H 9.S3D3.400

H 9.S3D4.150

H 9.S3D4.150

★

LAB ANALOG

Stainless Steel



H 50

H 51

H 52

★

IMPLANT PICK-UP

for open-tray with
integrated screw

Coming soon



H TR-NP013.4

H TR-RP024.1

H TR-WP035.0

ABUTMENT SCREW

Recommended tightening
torques 20 Ncm
Titanium Grade 5 ELI



H 60

H 60

H 60

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

H-SERIES

compatible to: Biomet 3i® / Certain®



PLATFORM

3,4 mm NP

4,1 mm RP

2-CONNECT®-ABUTMENT SET

incl. 2-CONNECT®-Coping
and Coping-screw
Titanium Grade 5 ELI



Gingiva Height	0,5 mm	H 800 S	H 810 S	★
	2,5 mm	H-2CON 8034002S	H-2CON 8041002S	★

2-CONNECT®-ABUTMENT

Recommended tightening
torques 20 Ncm
Titanium Grade 5 ELI

M (Male Part)



Gingiva Height	0,5 mm	H 800 M	H 810 M	★
	2,5 mm	H-2CON 8034002M	H-2CON 8041002M	★

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)	2CON F-U	2CON F-U	★
-----------------	----------	----------	---

COPING-SCREW

Recommended tightening
torques 15 Ncm
Titanium Grade 5 ELI



KS 60	KS 60	★
-------	-------	---

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

H-SERIES

compatible to: Biomet 3i® / Certain®



PLATFORM

3,4 mm NP

4,1 mm RP

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping,
Coping-screw & Male Part Screw
angled 20°
Gingiva Height 2,5 mm
Titanium Grade 5 ELI

Coming soon



H-2CON 80340A3S

H-2CON 80410A3S

2-CONNECT®-ABUTMENT

angled 20°
Gingiva Height 2,5 mm
Titanium Grade 5 ELI

Coming soon

M (Male Part)



H-2CON 80410A3M

H-2CON 80410A3M

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

★

COPING-SCREW

Recommended tightening
torques 15 Ncm
Titanium Grade 5 ELI

Coming soon



KS 62

KS 62

MALE PART SCREW

for angled 20° 2-CONnect®
Abutment | Recommended
tightening torques: 20 Ncm
Titanium Grade 5 ELI

Coming soon



H 6.Ms3.450

H 6.Ms3.450

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

H-SERIES

compatible to: Biomet 3i® / Certain®



PLATFORM

3,4 mm NP

4,1 mm RP

2-CONNECT®- IMPLANT PICK-UP

integrated screw



● Straight KS 60 = \varnothing 1,4 mm

2CON TR-00011.4

2CON TR-00011.4

● Angled KS 62 = \varnothing 1,6 mm

2CON TR-00031.6

2CON TR-00031.6

2-CONNECT®- SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



● Straight KS 60 = \varnothing 1,4 mm

2CON 9.S3D1.400

2CON 9.S3D1.400

★

● Angled KS 62 = \varnothing 1,6 mm

2CON 9.S3D1.600

2CON 9.S3D1.600

★

2-CONNECT®- LAB ANALOG

Stainless Steel



● Straight KS 60 = \varnothing 1,4 mm

2CON 5.La1.400

2CON 5.La1.400

★

● Angled KS 62 = \varnothing 1,6 mm

2CON 5.La1.600

2CON 5.La1.600

★

2-CONNECT®- DIM ANALOG



● Straight KS 60 = \varnothing 1,4 mm

2CON 5.DIM.140

2CON 5.DIM.140

★

● Angled KS 62 = \varnothing 1,6 mm

2CON 5.DIM.160

2CON 5.DIM.160

★

SLEEVE

plaster model
Stainless Steel



DIM 5.XSL.A30

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 

H-SERIES

compatible to: Biomet 3i® / Certain®



PLATFORM

3,4 mm NP

4,1 mm RP

5,0 mm WP

NT-PREFORM®

incl. coloured Milling Screw
Titanium Grade 5 ELI



Gingiva Height

1,3 mm

0,5 mm

H 10.Mi3.400

H 10.Mi4.100

H 10.Mi5.000

★

H 10.Mi3.400.05

H 10.Mi4.100.05

H 10.Mi5.000.05

★

NT-MILL LAB ANALOG

for nt-Preform®
Stainless Steel



H 5.La3.400

H 5.La4.100

H 5.La5.000

★

ABUTMENT GRIP

for nt-Holder fixture
Stainless Steel



W 12.100.000

W 12.200.000

W 12.300.000

★

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



H 9.S3D3.400

H 9.S3D4.150

H 9.S3D4.150

★

COLOURED MILLING SCREW

for nt-Preform®
Recommended tightening
torques 10 Ncm



H 6.Sc3.450

H 6.Sc3.450

H 6.Sc3.450

★

ABUTMENT SCREW

Recommended
tightening torques
20 Ncm
Titanium Grade 5 ELI



H 60

H 60

H 60

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

H-SERIES DIM ANALOG

compatible to: Biomet 3i® / Certain®



PLATFORM

3,4 mm NP

4,1 mm RP

5,0 mm WP

DIM ANALOG



H 5.DIM.340

H 5.DIM.410

H 5.DIM.500

★

SLEEVE

for plaster model
Stainless Steel



DIM 5.XSL.A10

DIM 5.XSL.A20

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 

NEXT GENERATION DESIGN – MADE IN GERMANY ///