

# HIO-SERIES

compatible to: HIOSSEN® / OSSTEM TS®



**PLATFORM**

3,5 mm

4,0 - 6,0 mm

**TITANIUM BASE**

for individual milled  
Zirconium Abutment  
incl. screw  
Titanium Grade 5 ELI



HIO 80350000

HIO 80406000

**SCAN BODY 3D-GUIDE®**

for intra + extra  
oral scanning



HIO 9.S3D3.500

HIO 9.S3D4.060

★

**NT-MILL LAB ANALOG**

Stainless Steel



HIO 5.La3.500

HIO 5.La4.060

★

**IMPLANT PICK-UP**

for open-tray with  
integrated screw

Coming soon



HIO TR-00013.5

HIO TR-00024.0

**ABUTMENT SCREW**

Recommended tightening  
torques 3,5 mm = 20 Ncm  
4,0 - 6,0 mm = 30 Ncm  
Titanium Grade 5 ELI



HIO 603500

HIO 604060

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA 🇺🇸

# HIO-SERIES

compatible to: HIOSSEN® / OSSTEM TS®



**PLATFORM**

3,5 mm

4,0 - 6,0 mm

**2-CONNECT®-ABUTMENT SET**

incl. 2-CONnect®-Coping  
and Coping-screw  
Titanium Grade 5 ELI

Coming soon



Gingiva Height : 0,5 mm  
: 2,5 mm

HIO-2CON 8035001S

HIO-2CON 8040601S

HIO-2CON 8035002S

HIO-2CON 8040602S

**2-CONNECT®-ABUTMENT**

Recommended tightening  
torques: 3,5 mm = 20 Ncm  
4,0 - 6,0 mm = 30 Ncm  
Titanium Grade 5 ELI

Coming soon

M (Male Part)



Gingiva Height : 0,5 mm  
: 2,5 mm

HIO-2CON 8035001M

HIO-2CON 8040601M

HIO-2CON 8035002M

HIO-2CON 8040602M

**2-CONNECT®-COPING**

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

★

**COPING-SCREW**

Recommended tightening  
torques:  
Titanium Grade 5 ELI



15 Ncm



25 Ncm

KS 61

KS 61

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# HIO-SERIES

compatible to: HIOSSEN® / OSSTEM TS®



**PLATFORM**

3,5 mm

4,0 - 6,0 mm

**2-CONNECT®-ABUTMENT SET**

incl. 2-CONnect®-Coping,  
Coping-screw & Male Part Screw  
angled 20°

Gingiva Height 2,5 mm  
Titanium Grade 5 ELI

Coming soon



HIO-2CON 80350A3S

HIO-2CON 80406A3S

**2-CONNECT®-ABUTMENT**

angled 20°

Gingiva Height 2,5 mm  
Titanium Grade 5 ELI

Coming soon

M (Male Part)



HIO-2CON 80350A3M

HIO-2CON 80406A3M

**2-CONNECT®-COPING**

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

★

**COPING-SCREW**

Recommended tightening  
torques: 15 Ncm

Titanium Grade 5 ELI

Coming soon



KS 62

KS 62

**MALE PART SCREW**

for angled 20° 2-CONnect®  
Abutment | Recommended  
tightening torques: 20 Ncm

Titanium Grade 5 ELI

Coming soon



HIO 6.Ms3.500

HIO 6.Ms4.060

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# HIO-SERIES

compatible to: HIOSSEN® / OSSTEM TS®



**PLATFORM**

3,5 mm

4,0 - 6,0 mm

**2-CONNECT®-  
IMPLANT PICK-UP**  
integrated screw



- Straight KS 61 =  $\varnothing$  1,8 mm
- Angled KS 62 =  $\varnothing$  1,6 mm

2CON TR-00021.8
2CON TR-00031.6

2CON TR-00021.8
2CON TR-00031.6

**2-CONNECT®-  
SCAN BODY 3D-GUIDE®**  
for intra + extra  
oral scanning



- Straight KS 61 =  $\varnothing$  1,8 mm
- Angled KS 62 =  $\varnothing$  1,6 mm

2CON 9.S3D1.800
2CON 9.S3D1.600

2CON 9.S3D1.800
2CON 9.S3D1.600

★  
★

**2-CONNECT®-  
LAB ANALOG**



- Straight KS 61 =  $\varnothing$  1,8 mm
- Angled KS 62 =  $\varnothing$  1,6 mm

2CON 5.La1.800
2CON 5.La1.600

2CON 5.La1.800
2CON 5.La1.600

★  
★

**2-CONNECT®-  
DIM ANALOG**



- Straight KS 61 =  $\varnothing$  1,8 mm
- Angled KS 62 =  $\varnothing$  1,6 mm

2CON 5.DIM.180
2CON 5.DIM.160

2CON 5.DIM.180
2CON 5.DIM.160

★  
★

**SLEEVE**  
for plaster model  
Stainless Steel



DIM 5.XSL.A30

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# HIO-SERIES

compatible to: HIOSSEN® / OSSTEM TS®



**PLATFORM**

3,5 mm

4,0 - 6,0 mm

**NT-PREFORM®**

incl. coloured milling screw  
Titanium Grade 5 ELI



Gingiva Height 0,5 mm

HIO 10.Mi3.500.05

HIO 10.Mi4.060.05

**NT-MILL LAB ANALOG**

for nt-Preform®  
Stainless Steel



HIO 5.La3.500

HIO 5.La4.060

★

**ABUTMENT GRIP**

for nt-Holder fixture  
Stainless Steel



W 12.100.000

W 12.200.000

★

**SCAN BODY 3D-GUIDE®**

for intra + extra  
oral scanning



HIO 9.S3D3.500

HIO 9.S3D4.060

★

**COLOURED MILLING SCREW**

for nt-Preform®  
Recommended tightening  
torques 10 Ncm



HIO 6.Sc3.500

HIO 6.Sc4.060

★

**ABUTMENT SCREW**

Recommended tightening  
torques 3,5 mm = 20 Ncm  
4,0 - 6,0 mm = 30 Ncm  
Titanium Grade 5 ELI



HIO 603500

HIO 604060

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

NEW

# HIO-SERIES DIM ANALOG

compatible to: HIOSSEN® / OSSTEM TS®



**PLATFORM**

3,5 mm

4,0 - 6,0 mm

**DIM ANALOG**



HIO 5.DIM.350

HIO 5.DIM.406

★

**SLEEVE**

for plaster model  
Stainless Steel



DIM 5.XSL.A10

DIM 5.XSL.A20

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA 