

K-SERIES

compatible to: Nobel Biocare™ / Brånemark®



PLATFORM

3,5 mm

4,1 mm

5,1 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screw
Titanium Grade 5 ELI



K 800

K 810

K 820

★

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



K 9.S3D3.500

K 9.S3D4.100

K 9.S3D5.100

★

LAB ANALOG

Stainless Steel



K 50

K 51

K 52

★

STRAIGHT ABUTMENT

Gingiva Height 2,5 mm
incl. screw
Titanium Grade 5 ELI



—

K 110

—

★

IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



K TR-00013.5

K TR-00024.1

—

ABUTMENT SCREW

Recommended tightening
torques 35 Ncm
Titanium Grade 5 ELI



K 60

K 61

K 62

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 

K-SERIES

compatible to: Nobel Biocare™ / Brånemark®



PLATFORM

3,5 mm

4,1 mm

5,1 mm

NT-PREFORM®

incl. coloured milling screw
Titanium Grade 5 ELI



Gingiva Height

1,3 mm

0,5 mm

K 10.Mi3.500

K 10.Mi4.100

K 10.Mi5.100



K 10.Mi3.500.05

K 10.Mi4.100.05

K 10.Mi5.100.05



NT-MILL LAB ANALOG

for nt-Preform®
Stainless Steel



K 5.La3.500

K 5.La4.100

K 5.La5.100



ABUTMENT GRIP

for nt-Holder fixture
Stainless Steel



W 12.100.000

W 12.200.000

W 12.300.000



SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



K 9.S3D3.500

K 9.S3D4.100

K 9.S3D5.100



COLOURED MILLING SCREW

for nt-Preform®
Recommended tightening
torques 10 Ncm



K 6.Sc3.500

K 6.Sc4.100

K 6.Sc5.100



ABUTMENT SCREW

Recommended
tightening torques 35 Ncm
Titanium Grade 5 ELI



K 60

K 61

K 62



Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

K-SERIES DIM ANALOG

compatible to: Nobel Biocare™ / Brånemark®



PLATFORM

3,5 mm

4,1 mm

5,1 mm

DIM ANALOG



K 5.DIM.350

K 5.DIM.410

K 5.DIM.510

★

SLEEVE

for plaster model
Stainless Steel



DIM 5.XSL.A10

DIM 5.XSL.A20

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 