

# MIC-SERIES

compatible to: MIS® / C1®



**PLATFORM**

3,3 mm NP

3,75 / 4,2 mm SP

5,0 mm WP

**TITANIUM BASE**

for individual milled Zirconium Abutment incl. screw  
Titanium Grade 5 ELI



MIC 80330000

MIC 80374200

MIC 80500000

**SCAN BODY 3D-GUIDE®**

for intra + extra oral scanning



MIC 9.S3D3.300

MIC 9.S3D3.742

MIC 9.S3D5.000

★

**NT-MILL LAB ANALOG**

Stainless Steel



MIC 5.La3.300

MIC 5.La3.742

MIC 5.La5.000

★

**IMPLANT PICK-UP**

for open-tray with integrated screw  
Coming soon



MIC TR-00013.3

MIC TR-00023.7

MIC TR-00035.0

**ABUTMENT SCREW**

Recommended tightening torques 30 Ncm  
Titanium Grade 5 ELI



MIC 603300

MIC 603750

MIC 603750

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures. Products indicated with ★ are available in the USA

# MIC-SERIES

compatible to: MIS® / C1®



**PLATFORM**

3,3 mm NP

3,75 / 4,2 mm SP

5,0 mm WP

**2-CONNECT®-ABUTMENT SET**

incl. 2-CONnect®-Coping  
and Coping-screw  
Titanium Grade 5 ELI



Gingiva Height : 0,5 mm  
: 2,5 mm

MIC-2CON 8033001S

MIC-2CON 8037421S

MIC-2CON 8050001S

MIC-2CON 8033002S

MIC-2CON 8037422S

MIC-2CON 8050002S **soon**

**2-CONNECT®-ABUTMENT**

Recommended tightening  
torques: 30 Ncm  
Titanium Grade 5 ELI  
**M** (Male Part)



Gingiva Height : 0,5 mm  
: 2,5 mm

MIC-2CON 8033001M

MIC-2CON 8037421M

MIC-2CON 8050001M

MIC-2CON 8033002M

MIC-2CON 8037422M

MIC-2CON 8050002M **soon**

**2-CONNECT®-COPING**

Titanium Grade 5 ELI



**F** (Female Part)

2CON F-U

2CON F-U

2CON F-U

★

**COPING-SCREW**

Recommended tightening  
torques: 25 Ncm  
Titanium Grade 5 ELI



KS 61

KS 61

KS 61

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# MIC-SERIES

compatible to: MIS® / C1®



**PLATFORM**

3,3 mm NP

3,75 / 4,2 mm SP

5,0 mm WP

**2-CONNECT®-ABUTMENT SET**

incl. 2-CONnect®-Coping,  
Coping-screw & Male Part Screw  
angled 20°  
Gingiva Height 2,5 mm  
Titanium Grade 5 ELI

Coming soon



over edge

MIC-2CON 80330A3S

MIC-2CON 80374A3S

MIC-2CON 80500A3S

over surface

MIC-2CON 80330A3S

—

—

**2-CONNECT®-ABUTMENT**

angled 20°  
Gingiva Height 2,5 mm  
Titanium Grade 5 ELI

**M (Male Part)**

Coming soon



over edge

MIC-2CON 80330A3M

MIC-2CON 80374A3M

MIC-2CON 80500A3M

over surface

MIC-2CON 80330A3M

—

—

**2-CONNECT®-COPING**

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

2CON F-U

★

**COPING-SCREW**

Recommended tightening  
torques: 15 Ncm

Titanium Grade 5 ELI

Coming soon



KS 62

KS 62

KS 62

**MALE PART SCREW**

for angled 20° 2-CONnect®  
Abutment | Recommended  
tightening torques: 20 Ncm

Titanium Grade 5 ELI

Coming soon



MIC 6.Ms3.300

MIC 6.Ms3.750

MIC 6.Ms3.750

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# MIC-SERIES

compatible to: MIS® / C1®



**PLATFORM**

3,3 mm NP

3,75 / 4,2 mm SP

5,0 mm WP

**2-CONNECT®-  
IMPLANT PICK-UP**  
integrated screw



- Straight KS 61 =  $\varnothing$  1,8 mm
- Angled KS 62 =  $\varnothing$  1,6 mm

2CON TR-00021.8
2CON TR-00031.6

2CON TR-00021.8
2CON TR-00031.6

2CON TR-00021.8
2CON TR-00031.6

**2-CONNECT®-  
SCAN BODY 3D-GUIDE®**  
for intra + extra  
oral scanning



- Straight KS 61 =  $\varnothing$  1,8 mm
- Angled KS 62 =  $\varnothing$  1,6 mm

2CON 9.S3D1.800
2CON 9.S3D1.600

2CON 9.S3D1.800
2CON 9.S3D1.600

2CON 9.S3D1.800
2CON 9.S3D1.600

★  
★

**2-CONNECT®-  
LAB ANALOG**



- Straight KS 61 =  $\varnothing$  1,8 mm
- Angled KS 62 =  $\varnothing$  1,6 mm

2CON 5.La1.800
2CON 5.La1.600

2CON 5.La1.800
2CON 5.La1.600

2CON 5.La1.800
2CON 5.La1.600

★  
★

**2-CONNECT®-  
DIM ANALOG**



- Straight KS 61 =  $\varnothing$  1,8 mm
- Angled KS 62 =  $\varnothing$  1,6 mm

2CON 5.DIM.180
2CON 5.DIM.160

2CON 5.DIM.180
2CON 5.DIM.160

2CON 5.DIM.180
2CON 5.DIM.160

★  
★

**SLEEVE**  
for plaster model  
Stainless Steel



DIM 5.XSL.A30

DIM 5.XSL.A30

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# MIC-SERIES

compatible to: MIS® / C1®



3,3 mm NP



3,75 / 4,2 mm SP



5,0 mm WP

## NT-PREFORM®

incl. coloured milling screw  
Titanium Grade 5 ELI



MIC 10.Mi3.300.05



MIC 10.Mi3.742.05




MIC 10.Mi5.000.05


Gingiva Height 0,5 mm

## NT-MILL LAB ANALOG

for nt-Preform®  
Stainless Steel



MIC 5.La3.300



MIC 5.La3.742



MIC 5.La5.000

★

## ABUTMENT GRIP

for nt-Holder fixture  
Stainless Steel



W 12.100.000



W 12.200.000



W 12.300.000

★

## SCAN BODY 3D-GUIDE®

for intra + extra  
oral scanning



MIC 9.S3D3.300



MIC 9.S3D3.742



MIC 9.S3D5.000

★

## COLOURED MILLING SCREW

for nt-Preform®  
Recommended tightening  
torques 10 Ncm



MIC 6.Sc3.750



MIC 6.Sc3.750



MIC 6.Sc3.750

★

## ABUTMENT SCREW

Recommended  
tightening torques 30 Ncm  
Titanium Grade 5 ELI



MIC 603300



MIC 603750



MIC 603750

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# MIC-SERIES DIM ANALOG

compatible to: MIS® / C1®

**NEW**



**PLATFORM**

3,3 mm NP

3,75 / 4,2 mm SP

5,0 mm WP

**DIM ANALOG**



MIC 5.DIM.330

MIC 5.DIM.374

MIC 5.DIM.500

★

**SLEEVE**

for plaster model  
Stainless Steel



DIM 5.XSL.A10

DIM 5.XSL.A20

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA 