

MIS-SERIES

compatible to: MIS® / SEVEN®



PLATFORM

3,3 mm NP

3,75 / 4,2 mm SP

5,0 mm WP

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screw
Titanium Grade 5 ELI



MIS 80330000

MIS 80374200

MIS 80500000

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



MIS 9.S3D3.300

MIS 9.S3D3.742

MIS 9.S3D5.000 ★

NT-MILL LAB ANALOG

Stainless Steel



MIS 5.La3.300

MIS 5.La3.742

MIS 5.La5.000 ★

IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



MIS TR-00013.3

MIS TR-00023.7

MIS TR-00035.0

ABUTMENT SCREW

Recommended tightening
torques 30 Ncm
Titanium Grade 5 ELI



MIS 603300

MIS 603750

MIS 603750

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

MIS-SERIES

compatible to: MIS® / SEVEN®



PLATFORM

3,3 mm NP

3,75 / 4,2 mm SP

5,0 mm WP

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping
and Coping-screw
Titanium Grade 5 ELI



Gingiva Height : 0,5 mm
: 2,5 mm

MIS-2CON 8033001S

MIS-2CON 8037421S

MIS-2CON 8050001S

MIS-2CON 8033002S

MIS-2CON 8037422S

MIS-2CON 8050002S *soon*

2-CONNECT®-ABUTMENT

Recommended tightening
torques: 30 Ncm
Titanium Grade 5 ELI
M (Male Part)



Gingiva Height : 0,5 mm
: 2,5 mm

MIS-2CON 8033001M

MIS-2CON 8037421M

MIS-2CON 8050001M

MIS-2CON 8033002M

MIS-2CON 8037422M

MIS-2CON 8050002M *soon*

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

2CON F-U

★

COPING-SCREW

Recommended tightening
torques: 25 Ncm
Titanium Grade 5 ELI



KS 61

KS 61

KS 61

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

MIS-SERIES

compatible to: MIS® / SEVEN®



PLATFORM

3,3 mm NP

3,75 / 4,2 mm SP

5,0 mm WP

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping,
Coping-screw & Male Part Screw
angled 20°
Gingiva Height 2,5 mm
Titanium Grade 5 ELI
Coming soon



MIS-2CON 8033002S

MIS-2CON 8037422S

MIS-2CON 8050002S

2-CONNECT®-ABUTMENT

angled 20°
Gingiva Height 2,5 mm
Titanium Grade 5 ELI
Coming soon
M (Male Part)



MIS-2CON 8033002M

MIS-2CON 8037422M

MIS-2CON 8050002M

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

2CON F-U

★

COPING-SCREW

Recommended tightening
torques: 15 Ncm
Titanium Grade 5 ELI
Coming soon



KS 62

KS 62

KS 62

MALE PART SCREW

for angled 20° 2-CONnect®
Abutment | Recommended
tightening torques: 20 Ncm
Titanium Grade 5 ELI
Coming soon



MIS 6.Ms3.300

MIS 6.MS3.750

MIS 6.MS3.750

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

MIS-SERIES

compatible to: MIS® / SEVEN®



PLATFORM

3,3 mm NP

3,75 / 4,2 mm SP

5,0 mm WP

**2-CONNECT®-
IMPLANT PICK-UP**
integrated screw



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON TR-00021.8
2CON TR-00031.6

2CON TR-00021.8
2CON TR-00031.6

2CON TR-00021.8
2CON TR-00031.6

**2-CONNECT®-
SCAN BODY 3D-GUIDE®**
for intra + extra
oral scanning



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON 9.S3D1.800
2CON 9.S3D1.600

2CON 9.S3D1.800
2CON 9.S3D1.600

2CON 9.S3D1.800
2CON 9.S3D1.600

★
★

**2-CONNECT®-
LAB ANALOG**



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON 5.La1.800
2CON 5.La1.600

2CON 5.La1.800
2CON 5.La1.600

2CON 5.La1.800
2CON 5.La1.600

★
★

**2-CONNECT®-
DIM ANALOG**



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON 5.DIM.180
2CON 5.DIM.160

2CON 5.DIM.180
2CON 5.DIM.160

2CON 5.DIM.180
2CON 5.DIM.160

★
★

SLEEVE
for plaster model
Stainless Steel



DIM 5.XSL.A30

DIM 5.XSL.A30

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

MIS-SERIES

compatible to: MIS® / SEVEN®



PLATFORM

3,3 mm NP

3,75 / 4,2 mm SP

5,0 mm WP

NT-PREFORM®

incl. coloured milling screw
Titanium Grade 5 ELI



Gingiva Height 0,5 mm



MIS 10.Mi3.300.05

MIS 10.Mi3.742.05

MIS 10.Mi5.000.05

NT-MILL LAB ANALOG

for nt-Preform®
Stainless Steel



MIS 5.La3.300

MIS 5.La3.742

MIS 5.La5.000

★

ABUTMENT GRIP

for nt-Holder fixture
Stainless Steel



W 12.100.000

W 12.200.000

W 12.300.000

★

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



MIS 9.S3D3.300

MIS 9.S3D3.742

MIS 9.S3D5.000

★

COLOURED MILLING SCREW

for nt-Preform®
Recommended tightening
torques 10 Ncm



MIS 6.Sc3.300

MIS 6.Sc3.750

MIS 6.Sc3.750

★

ABUTMENT SCREW

Recommended
tightening torques 30 Ncm
Titanium Grade 5 ELI



MIS 603300

MIS 603750

MIS 603750

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

MIS-SERIES DIM ANALOG

compatible to: MIS® / SEVEN®



PLATFORM

3,3 mm NP

3,75 / 4,2 mm SP

5,0 mm WP

DIM ANALOG



MIS 5.DIM.330

MIS 5.DIM.374

MIS 5.DIM.500

★

SLEEVE

for plaster model
Stainless Steel



DIM 5.XSL.A10

DIM 5.XSL.A20

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA