

UPR-SERIES

compatible to: Sweden & Martina Implantology™ / Premium Kohno®



PLATFORM

3,3 mm

3,8 mm

4,25 mm

5,0 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screw
Titanium Grade 5 ELI



UPR 80330000

UPR 80380000

UPR 80420000

UPR 80500000

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



UPR 9.S3D3.300

UPR 9.S3D3.800


UPR 9.S3D4.200

UPR 9.S3D5.000

★

NT-MILL LAB ANALOG

Stainless Steel



UPR 5.La3.300

UPR 5.La3.800

UPR 5.La4.200

UPR 5.La5.000

★

IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



UPR TR-00013.3




UPR TR-00023.8

UPR TR-00034.2

UPR TR-00045.0

ABUTMENT SCREW

Recommended tightening
torques 25 Ncm
Titanium Grade 5 ELI



UPR 603338

UPR 603338

UPR 604250

UPR 604250

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures. Products indicated with ★ are available in the USA 🇺🇸

UPR-SERIES

compatible to: Sweden & Martina Implantology™ / Premium Kohno®



PLATFORM

3,3 mm

3,8 mm

4,25 mm

5,0 mm

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping and Coping-screw
Titanium Grade 5 ELI
Coming soon



Gingiva Height	0,5 mm	UPR-2CON 8033001S	UPR-2CON 8038001S	UPR-2CON 8042001S	UPR-2CON 8050001S
	2,5 mm	UPR-2CON 8033002S	UPR-2CON 8038002S	UPR-2CON 8042002S	UPR-2CON 8050002S

2-CONNECT®-ABUTMENT

Recommended tightening torques: 25 Ncm
Titanium Grade 5 ELI
Coming soon
M (Male Part)



Gingiva Height	0,5 mm	UPR-2CON 8033001M	UPR-2CON 8038001M	UPR-2CON 8042001M	UPR-2CON 8050001M
	2,5 mm	UPR-2CON 8033002M	UPR-2CON 8038002M	UPR-2CON 8042002M	UPR-2CON 8050002M

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)	2CON F-U	2CON F-U	2CON F-U	2CON F-U	★
-----------------	----------	----------	----------	----------	---

COPING-SCREW

Recommended tightening torques: 20 Ncm
Titanium Grade 5 ELI



KS 61	KS 61	KS 61	KS 61
-------	-------	-------	-------

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures. Products indicated with ★ are available in the USA

UPR-SERIES

compatible to: Sweden & Martina Implantology™ / Premium Kohno®



PLATFORM

3,3 mm

3,8 mm

4,25 mm

5,0 mm

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping,
Coping-screw & Male Part Screw
angled 20°
Gingiva Height 2,5 mm
Titanium Grade 5 ELI

Coming soon



UPR-2CON 80330A3S

UPR-2CON 80380A3S

UPR-2CON 80420A3S

UPR-2CON 80500A3S

2-CONNECT®-ABUTMENT

angled 20°
Gingiva Height 2,5 mm
Titanium Grade 5 ELI

Coming soon

M (Male Part)



UPR-2CON 80330A3M

UPR-2CON 80380A3M

UPR-2CON 80420A3M

UPR-2CON 80500A3M

2-CONNECT®-COPING

Titanium Grade 5 ELI

F (Female Part)



2CON F-U

2CON F-U

2CON F-U


2CON F-U

★

COPING-SCREW

Recommended tightening
torques: 15 Ncm
Titanium Grade 5 ELI

Coming soon



KS 62

KS 62

KS 62

KS 62

MALE PART SCREW

for angled 20° 2-CONnect®
Abutment | Recommended
tightening torques: 20 Ncm
Titanium Grade 5 ELI

Coming soon



UPR 6.Ms.3.338

UPR 6.Ms.3.338

UPR 6.Ms4.250

UPR 6.Ms4.250

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

UPR-SERIES

compatible to: Sweden & Martina Implantology™ / Premium Kohno®



PLATFORM

3,3 mm

3,8 mm

4,25 mm

5,0 mm

**2-CONNECT®-
IMPLANT PICK-UP**
integrated screw



● Straight KS 61 = \varnothing 1,8 mm	2CON TR-00021.8	2CON TR-00021.8	2CON TR-00021.8	2CON TR-00021.8
● Angled KS 62 = \varnothing 1,6 mm	2CON TR-00031.6	2CON TR-00031.6	2CON TR-00031.6	2CON TR-00031.6

**2-CONNECT®-
SCAN BODY 3D-GUIDE®**
for intra + extra
oral scanning



● Straight KS 61 = \varnothing 1,8 mm	2CON 9.S3D1.800	2CON 9.S3D1.800	2CON 9.S3D1.800	2CON 9.S3D1.800	★
● Angled KS 62 = \varnothing 1,6 mm	2CON 9.S3D1.600	2CON 9.S3D1.600	2CON 9.S3D1.600	2CON 9.S3D1.600	★

**2-CONNECT®-
LAB ANALOG**



● Straight KS 61 = \varnothing 1,8 mm	2CON 5.La1.800	2CON 5.La1.800	2CON 5.La1.800	2CON 5.La1.800	★
● Angled KS 62 = \varnothing 1,6 mm	2CON 5.La1.600	2CON 5.La1.600	2CON 5.La1.600	2CON 5.La1.600	★

**2-CONNECT®-
DIM ANALOG**



● Straight KS 61 = \varnothing 1,8 mm	2CON 5.DIM.180	2CON 5.DIM.180	2CON 5.DIM.180	2CON 5.DIM.180	★
● Angled KS 62 = \varnothing 1,6 mm	2CON 5.DIM.160	2CON 5.DIM.160	2CON 5.DIM.160	2CON 5.DIM.160	★

SLEEVE
for plaster model
Stainless Steel



DIM 5.XSL.A30	DIM 5.XSL.A30	DIM 5.XSL.A30	DIM 5.XSL.A30	★
---------------	---------------	---------------	---------------	---

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures. Products indicated with ★ are available in the USA

UPR-SERIES

compatible to: Sweden & Martina Implantology™ / Premium Kohno®



PLATFORM

3,3 mm

3,8 mm

4,25 mm

5,0 mm

NT-PREFORM®

incl. coloured milling screw
Titanium Grade 5 ELI



Gingiva Height 0,5 mm

UPR 10.Mi3.300.05

UPR 10.Mi3.800.05

UPR 10.Mi4.200.05

UPR 10.Mi5.000.05

NT-MILL LAB ANALOG

for nt-Preform®
Stainless Steel



UPR 5.La3.300

UPR 5.La3.800

UPR 5.La4.200

UPR 5.La5.000

★

ABUTMENT GRIP

for nt-Holder fixture
Stainless Steel



W 12.100.000

W 12.200.000

W 12.200.000

W 12.300.000

★

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



UPR 9.S3D3.300

UPR 9.S3D3.800

UPR 9.S3D4.200

UPR 9.S3D5.000

★

COLOURED MILLING SCREW

for nt-Preform®
Recommended tightening
torques 10 Ncm



UPR 6.Sc3.338

UPR 6.Sc3.338

UPR 6.Sc4.250

UPR 6.Sc4.250

★

ABUTMENT SCREW

Recommended tightening
torques 25 Ncm
Titanium Grade 5 ELI



UPR 603338

UPR 603338


UPR 604250

UPR 604250

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

UPR-SERIES DIM ANALOG

compatible to: Sweden & Martina Implantology™ / Premium Kohno®

					
PLATFORM	3,3 mm	3,8 mm	4,25 mm	5,0 mm	
DIM ANALOG					
	UPR 5.DIM.330	UPR 5.DIM.380	UPR 5.DIM.420	UPR 5.DIM.500	★
SLEEVE for plaster model Stainless Steel					
	DIM 5.XSL.A10	DIM 5.XSL.A20	DIM 5.XSL.A20	DIM 5.XSL.A30	★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
 Products indicated with ★ are available in the USA 