

UTG-SERIES

compatible to: Sweden & Martina Implantology™ / Premium TG®



PLATFORM

3,3 - 4,25 mm

5,0 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screw
Titanium Grade 5 ELI



UTG 80334200

UTG 80500000

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning




UTG 9.S3D3.342

UTG 9.S3D5.000

★

NT-MILL LAB ANALOG

Stainless Steel



UTG 5.La3.342

UTG 5.La5.000

★

IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



UTG TR-00013.3

UTG TR-00025.0

ABUTMENT SCREW

Recommended tightening
torques 25 Ncm
Titanium Grade 5 ELI



UTG 603342

UTG 605000

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

UTG-SERIES

compatible to: Sweden & Martina Implantology™ / Premium TG®



PLATFORM

3,3 - 4,25 mm

5,0 mm

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping
and Coping-screw
Titanium Grade 5 ELI

Coming soon



Gingiva Height

0,5 mm	
2,5 mm	

UTG-2CON 8033421S

UTG-2CON 8050001S

UTG-2CON 8033422S

UTG-2CON 8050002S

2-CONNECT®-ABUTMENT

Recommended tightening
torques 25 Ncm
Titanium Grade 5 ELI

Coming soon

M (Male Part)



Gingiva Height

0,5 mm	
2,5 mm	

UTG-2CON 8033421M

UTG-2CON 8050001M

UTG-2CON 8033422M

UTG-2CON 8050002M

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

★

COPING-SCREW

Recommended tightening
torques: 20 Ncm
Titanium Grade 5 ELI



KS 61

KS 61

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

UTG-SERIES

compatible to: Sweden & Martina Implantology™ / Premium TG®



PLATFORM

3,3 - 4,25 mm

5,0 mm

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping,
Coping-screw & Male Part Screw
angled 20°
Gingiva Height 2,5 mm
Titanium Grade 5 ELI

Coming soon



UTG-2CON 80334A3S

UTG-2CON 80500A3S

2-CONNECT®-ABUTMENT

angled 20°
Gingiva Height 2,5 mm
Titanium Grade 5 ELI

Coming soon

M (Male Part)



UTG-2CON 80334A3M

UTG-2CON 80500A3M

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U

★

COPING-SCREW

Recommended tightening
torques: 15 Ncm

Titanium Grade 5 ELI

Coming soon



KS 62

KS 62

MALE PART SCREW

for angled 20° 2-CONnect®
Abutment | Recommended
tightening torques: 20 Ncm
Titanium Grade 5 ELI

Coming soon



UTG 6.Ms3.342

UTG 6.Ms5.000

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

UTG-SERIES

compatible to: Sweden & Martina Implantology™ / Premium TG®



PLATFORM

3,3 - 4,25 mm

5,0 mm

2-CONNECT®- IMPLANT PICK-UP

integrated screw



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON TR-00021.8	2CON TR-00031.6
-----------------	-----------------

2CON TR-00021.8	2CON TR-00031.6
-----------------	-----------------

2-CONNECT®- SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON 9.S3D1.800	2CON 9.S3D1.600
-----------------	-----------------

2CON 9.S3D1.800	2CON 9.S3D1.600
-----------------	-----------------

★

★

2-CONNECT®- LAB ANALOG



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON 5.La1.800	2CON 5.La1.600
----------------	----------------

2CON 5.La1.800	2CON 5.La1.600
----------------	----------------

★

★

2-CONNECT®- DIM ANALOG



- Straight KS 61 = \varnothing 1,8 mm
- Angled KS 62 = \varnothing 1,6 mm

2CON 5.DIM.180	2CON 5.DIM.160
----------------	----------------

2CON 5.DIM.180	2CON 5.DIM.160
----------------	----------------

★

★

SLEEVE

for plaster model
Stainless Steel



DIM 5.XSL.A30	DIM 5.XSL.A30
---------------	---------------

DIM 5.XSL.A30	DIM 5.XSL.A30
---------------	---------------

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

UTG-SERIES

compatible to: Sweden & Martina Implantology™ / Premium TG®



PLATFORM

3,3 - 4,25 mm

5,0 mm

NT-PREFORM®

incl. coloured milling screw
Titanium Grade 5 ELI



Gingiva Height 0,5 mm

UTG 10.Mi3.342.05

UTG 10.Mi5.000.05

NT-MILL LAB ANALOG

for nt-Preform®
Stainless Steel



UTG 5.La3.342

UTG 5.La5.000

★

ABUTMENT GRIP

for nt-Holder fixture
Stainless Steel



W 12.200.000

W 12.300.000

★

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



UTG 9.S3D3.342

UTG 9.S3D5.000

★

COLOURED MILLING SCREW

for nt-Preform®
Recommended tightening
torques 10 Ncm



UTG 6.Sc3.342

UTG 6.Sc5.000

★

ABUTMENT SCREW

Recommended tightening
torques 25 Ncm
Titanium Grade 5 ELI



UTG 603342

UTG 605000

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

UTG-SERIES DIM ANALOG

compatible to: Sweden & Martina Implantology™ / Premium TG®

NEW



PLATFORM

3,3 - 4,25 mm

5,0 mm

DIM ANALOG



UTG 5.DIM.334

UTG 5.DIM.500

★

SLEEVE

for plaster model
Stainless Steel



DIM 5.XSL.A20

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 

NEXT GENERATION DESIGN – MADE IN GERMANY ///