


T-SERIES

compatible to:


DENTSPLY Implants®
Xive®



PLATFORMS:




3,4 mm



3,8 mm



4,5 mm



5,5 mm

T-SERIES

TITANIUM BASE

incl. screw

T 800



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

T 9.S3D3.400



LAB ANALOG

T 50



IMPLANT PICK-UP

for open-tray with integrated screw

T TR-00013.4



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

T 60



DIM ANALOG®

T 5.DIM.340



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A10



**NT-PREFORM®**

incl. coloured milling screw

GH 0,5 mm

T 10.Mi3.400.05

GH 1,3 mm

T 10.Mi3.400

**SCAN BODY 3D-GUIDE®**

for intra + extra oral scanning

T 9.S3D3.400

**NT-MILL LAB ANALOG**


for NT-Preform®

T 5.La3.400

**COLOURED MILLING SCREW**

for NT-Preform®

T 6.Sc3.455

**ABUTMENT SCREW**

Overview of all screwdrivers on pages 344 and 345

T 60

**DIM ANALOG®**

T 5.DIM.340

**SLEEVE** (belongs to DIM ANALOG®)
plaster model

DIM 5.XSL.A10

**ABUTMENT GRIP**

for NT-Holder fixture

W 12.100.000

**CUSTOMIZED ABUTMENT**

milled by NT-iDent

Not available in the USA

TZ15.TiO.Si1



T-SERIES

2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm

T 800 S

GH 2,5 mm

T-2CON 8034002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm

T 800 M

GH 2,5 mm

T-2CON 8034002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 60 = 1,4 mm

2CON 9.S3D1.400



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 15 Ncm

KS 60



2-CONNECT®-DIM ANALOG®

KS 60 = 1,4 mm

2CON 5.DIM.140



SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 60 = 1,4 mm

2CON TR-00011.4



2-CONNECT®-LAB ANALOG

KS 60 = 1,4 mm

2CON 5.La1.400



**2-CONNECT®-ABUTMENT SET (ANGLED)**

incl. 2-CONnect®-Coping, mounting aid,
Coping-screw and Male Part Screw, angled 20°,
Gingiva Height 2,5 mm

T-2CON 80340A3S

**2-CONNECT®-ABUTMENT**

angled 20°, Gingiva Height 2,5 mm
Male Part incl. screw,
mounting aid

T-2CON 80340A3M

**2-CONNECT®-SCAN BODY 3D-GUIDE®**

for intra + extra oral scanning
KS 62 = 1,6 mm

2CON 9.S3D1.600

**2-CONNECT®-COPING**

Female Part

2CON F-U

**COPING-SCREW**

Recommended tightening torques: 15 Ncm

KS 62

**MALE PART SCREW**

for angled 20° 2-CONnect® Abutment
Recommended tightening torques: 20 Ncm

T 6.Ms3.455

**2-CONNECT®-DIM ANALOG®**

KS 62 = 1,6 mm

2CON 5.DIM.160



SLEEVE (belongs to DIM ANALOG®)
plaster model

DIM 5.XSL.A30

**2-CONNECT®-IMPLANT PICK-UP**

integrated screw
KS 62 = 1,6 mm

2CON TR-00031.6

**2-CONNECT®-LAB ANALOG**

KS 62 = 1,6 mm

2CON 5.La1.600



T-SERIES

TITANIUM BASE

incl. screw

T 805



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

T 9.S3D3.800



LAB ANALOG

T 55



IMPLANT PICK-UP

for open-tray with integrated screw

T TR-00023.8



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

T 60



DIM ANALOG®

T 5.DIM.380



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A20



**NT-PREFORM®**

incl. coloured milling screw

GH 0,5 mm

T 10.Mi3.800.05

GH 1,3 mm

T 10.Mi3.800

**SCAN BODY 3D-GUIDE®**

for intra + extra oral scanning

T 9.S3D3.800

**NT-MILL LAB ANALOG**

for NT-Preform®

T 5.La3.800



groove

COLOURED MILLING SCREW

for NT-Preform®

T 6.Sc3.455

**ABUTMENT SCREW**

Overview of all screwdrivers on pages 344 and 345

T 60

**DIM ANALOG®**

T 5.DIM.380

**SLEEVE** (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A20

**ABUTMENT GRIP**

for NT-Holder fixture

W 12.200.000

**CUSTOMIZED ABUTMENT**

milled by NT-iDent

Not available in the USA

TZ15.TiO.Si1



T-SERIES

2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm

T 805 S

GH 2,5 mm

T-2CON 8038002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm

T 805 M

GH 2,5 mm

T-2CON 8038002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 20 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800



**2-CONNECT®-ABUTMENT SET (ANGLED)**

incl. 2-CONnect®-Coping, mounting aid,
Coping-screw and Male Part Screw, angled 20°,
Gingiva Height 2,5 mm

T-2CON 80380A3S

**2-CONNECT®-ABUTMENT**

angled 20°, Gingiva Height 2,5 mm
Male Part incl. screw,
mounting aid

T-2CON 80380A3M

**2-CONNECT®-SCAN BODY 3D-GUIDE®**

for intra + extra oral scanning
KS 62 = 1,6 mm

2CON 9.S3D1.600

**2-CONNECT®-COPING**

Female Part

2CON F-U

**COPING-SCREW**

Recommended tightening torques: 15 Ncm

KS 62

**MALE PART SCREW**

for angled 20° 2-CONnect® Abutment
Recommended tightening torques: 20 Ncm

T 6.Ms3.455

**2-CONNECT®-DIM ANALOG®**

KS 62 = 1,6 mm

2CON 5.DIM.160



SLEEVE (belongs to DIM ANALOG®)
plaster model

DIM 5.XSL.A30

**2-CONNECT®-IMPLANT PICK-UP**

integrated screw
KS 62 = 1,6 mm

2CON TR-00031.6

**2-CONNECT®-LAB ANALOG**

KS 62 = 1,6 mm

2CON 5.La1.600



T-SERIES

TITANIUM BASE

incl. screw

T 810



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

T 9.S3D4.555



LAB ANALOG

T 51



IMPLANT PICK-UP

for open-tray with integrated screw

T TR-00034.5



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

T 60



DIM ANALOG®

T 5.DIM.450



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A20



**NT-PREFORM®**

incl. coloured milling screw

GH 0,5 mm

T 10.Mi4.500.05

GH 1,3 mm

T 10.Mi4.500

**SCAN BODY 3D-GUIDE®**

for intra + extra oral scanning

T 9.S3D4.555

**NT-MILL LAB ANALOG**

for NT-Preform®

T 5.La4.500



groove

COLOURED MILLING SCREW

for NT-Preform®

T 6.Sc3.455

**ABUTMENT SCREW**

Overview of all screwdrivers on pages 344 and 345

T 60

**DIM ANALOG®**

T 5.DIM.450

**SLEEVE** (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A20

**ABUTMENT GRIP**

for NT-Holder fixture

W 12.200.000

**CUSTOMIZED ABUTMENT**

milled by NT-iDent

Not available in the USA

TZ15.TiO.Si1



T-SERIES

2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm

T 810 S

GH 2,5 mm

T-2CON 8045002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm

T 810 M

GH 2,5 mm

T-2CON 8045002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 20 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800



**2-CONNECT®-ABUTMENT SET (ANGLED)**

incl. 2-CONnect®-Coping, mounting aid,
Coping-screw and Male Part Screw, angled 20°,
Gingiva Height 2,5 mm

T-2CON 80450A3S

**2-CONNECT®-ABUTMENT**

angled 20°, Gingiva Height 2,5 mm
Male Part incl. screw,
mounting aid

T-2CON 80450A3M

**2-CONNECT®-SCAN BODY 3D-GUIDE®**

for intra + extra oral scanning
KS 62 = 1,6 mm

2CON 9.S3D1.600

**2-CONNECT®-COPING**

Female Part

2CON F-U

**COPING-SCREW**

Recommended tightening torques: 15 Ncm

KS 62

**MALE PART SCREW**

for angled 20° 2-CONnect® Abutment
Recommended tightening torques: 20 Ncm

T 6.Ms3.455

**2-CONNECT®-DIM ANALOG®**

KS 62 = 1,6 mm

2CON 5.DIM.160



SLEEVE (belongs to DIM ANALOG®)
plaster model

DIM 5.XSL.A30

**2-CONNECT®-IMPLANT PICK-UP**

integrated screw
KS 62 = 1,6 mm

2CON TR-00031.6

**2-CONNECT®-LAB ANALOG**

KS 62 = 1,6 mm

2CON 5.La1.600



T-SERIES

TITANIUM BASE

incl. screw

T 820



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

T 9.S3D4.555



LAB ANALOG

T 52



IMPLANT PICK-UP

for open-tray with integrated screw

T TR-00034.5



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

T 60



DIM ANALOG®

T 5.DIM.550



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



**NT-PREFORM®**

incl. coloured milling screw

GH 0,5 mm

T 10.Mi5.500.05

GH 1,3 mm

T 10.Mi5.500

**SCAN BODY 3D-GUIDE®**

for intra + extra oral scanning

T 9.S3D4.555

**NT-MILL LAB ANALOG**

for NT-Preform®

T 5.La5.500



groove

COLOURED MILLING SCREW

for NT-Preform®

T 6.Sc3.455

**ABUTMENT SCREW**

Overview of all screwdrivers on pages 344 and 345

T 60

**DIM ANALOG®**

T 5.DIM.550

**SLEEVE** (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30

**ABUTMENT GRIP**

for NT-Holder fixture

W 12.300.000

**CUSTOMIZED ABUTMENT**

milled by NT-iDent

Not available in the USA

TZ15.TiO.Si1

