

BIO -SERIES



compatible to:
BioHorizons® / internal



3,0 mm




4,5 mm



3,5 mm




5,7 mm



**PRECISION
MEETS QUALITY.**

**OVERVIEW
PLATFORMS:**



3,0 mm



3,5 mm




4,5 mm



5,7 mm

 **MADE IN
GERMANY**



3,0 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

BIO 80300000



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,5 mm BIO 82300005

GH 1,5 mm BIO 82300015



Scanbody 3D-Guide®

for intra + extra oral scanning

BIO 9.S3D3.000



Lab Analog

BIO 5.La3.000



Implant Pick.Up

for open-tray with integrated screw

BIO TR-00013.0



Abutment Screw

Recommended tightening torques 30 Ncm

BIO 603057



DIM-ANALOG®

BIO 5.DIM.300



Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSL.A10

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,0 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



nt-Preform®

incl. final Abutment Screw

GH 0,5 mm BIO 10.Mi3.000.05



Scanbody 3D-Guide®

for intra + extra oral scanning

BIO 9.S3D3.000



nt-Mill Lab Analog

for nt-Preform®

BIO 5.La3.000



colored milling screw

for nt-Preform®

BIO 6.Sc3.057



Abutment Screw

Recommended tightening torques 30 Ncm

BIO 603057



DIM-ANALOG®

BIO 5.DIM.300



Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSLA10



Abutment Grip

for nt-Preform® Holder

W12.100.00



Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,0 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm BIO-2CON 8030001S

GH 2,5 mm BIO-2CON 8030002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 30 Ncm

GH 0,5 mm BIO-2CON 8030001M

GH 2,5 mm BIO-2CON 8030002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,0 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm BIO-2CON 80300A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm BIO-2CON 80300A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

BIO 6.Ms3.057



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



3,5 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

BIO 80350000



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,5 mm BIO 82350005

GH 1,5 mm BIO 82350015



Scanbody 3D-Guide®

for intra + extra oral scanning

BIO 9.S3D3.500



Lab Analog

BIO 5.La3.500



Implant Pick.Up

for open-tray with integrated screw

BIO TR-00023.5



Abutment Screw

Recommended tightening torques 30 Ncm

BIO 603057



DIM-ANALOG®

BIO 5.DIM.350



Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSL.A10

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,5 mm
BioHorizons® / internal

ALL DETAILS
ONLINE



nt-Preform®
incl. final Abutment Screw

GH 0,5 mm BIO 10.Mi3.500.05



Scanbody 3D-Guide®
for intra + extra oral scanning

BIO 9.S3D3.500



nt-Mill Lab Analog
for nt-Preform®

BIO 5.La3.500



colored milling screw
for nt-Preform®

BIO 6.Sc3.057



Abutment Screw
Recommended tightening torques 30 Ncm

BIO 603057



DIM-ANALOG®

BIO 5.DIM.350



Sleeve (für DIM-ANALOG®)
for plaster model

DIM 5.XSL.A10



Abutment Grip
for nt-Preform® Holder

W12.100.00



Customized Abutment
milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,5 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm BIO-2CON 8035001S

GH 2,5 mm BIO-2CON 8035002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 30 Ncm

GH 0,5 mm BIO-2CON 8035001M

GH 2,5 mm BIO-2CON 8035002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,5 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm BIO-2CON 80350A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm BIO-2CON 80350A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

BIO 6.Ms3.057



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



4,5 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

BIO 80450000



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,5 mm BIO 82450005

GH 1,5 mm BIO 82450015



Scanbody 3D-Guide®

for intra + extra oral scanning

BIO 9.S3D4.500



Lab Analog

BIO 5.La4.500



Implant Pick.Up

for open-tray with integrated screw

BIO TR-00034.5



Abutment Screw

Recommended tightening torques 30 Ncm

BIO 603057



DIM-ANALOG®

BIO 5.DIM.450



Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSLA30

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,5 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



nt-Preform®

incl. final Abutment Screw

GH 0,5 mm BIO 10.Mi4.500.05



Scanbody 3D-Guide®

for intra + extra oral scanning

BIO 9.S3D4.500



nt-Mill Lab Analog

for nt-Preform®

BIO 5.La4.500



colored milling screw

for nt-Preform®

BIO 6.SC3.057



Abutment Screw

Recommended tightening torques 30 Ncm

BIO 603057



DIM-ANALOG®

BIO 5.DIM.450



Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSLA30



Abutment Grip

for nt-Preform® Holder

W12.200.00



Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,5 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm BIO-2CON 8045001S

GH 2,5 mm BIO-2CON 8045002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 30 Ncm

GH 0,5 mm BIO-2CON 8045001M

GH 2,5 mm BIO-2CON 8045002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,5 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm BIO-2CON 80450A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm BIO-2CON 80450A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

BIO 6.Ms3.057



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



5,7 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

BIO 80570000



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,5 mm BIO 82570005

GH 1,5 mm BIO 82570015



Scanbody 3D-Guide®

for intra + extra oral scanning

BIO 9.S3D5.700



Lab Analog

BIO 5.La5.700



Implant Pick.Up

for open-tray with integrated screw

BIO TR-00045.7



Abutment Screw

Recommended tightening torques 30 Ncm

BIO 603057



DIM-ANALOG®

BIO 5.DIM.570



Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSL.A60

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



5,7 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



nt-Preform®

incl. final Abutment Screw

GH 0,5 mm BIO 10.Mi5.700.05



Scanbody 3D-Guide®

for intra + extra oral scanning

BIO 9.S3D5.700



nt-Mill Lab Analog

for nt-Preform®

BIO 5.La5.700



colored milling screw

for nt-Preform®

BIO 6.SC3.057



Abutment Screw

Recommended tightening torques 30 Ncm

BIO 603057



DIM-ANALOG®

BIO 5.DIM.570



Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSL.A60



Abutment Grip

for nt-Preform® Holder

W12.300.00



Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



5,7 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm BIO-2CON 8057001S

GH 2,5 mm BIO-2CON 8057002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 30 Ncm

GH 0,5 mm BIO-2CON 8057001M

GH 2,5 mm BIO-2CON 8057002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



5,7 mm

BioHorizons® / internal

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm BIO-2CON 80570A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm BIO-2CON 80570A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

BIO 6.Ms3.057



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

BIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL

NOTES

[illegible]

EUROPE

NT-Trading GmbH & Co. KG

G.-Braun-Str. 18
76187 Karlsruhe
Germany

E info@nt.dental

W www.nt.dental

NORTH AMERICA

NT-US, Inc.

132 Pepes Farm Road,
Unit 3
Milford, CT 06460
USA

E info@nt.dental (general information)

E order-us@nt.dental (for orders)

W www.nt.dental

2023_NT_Serienbroschüre_BIO_EN_REV001