


CAM -SERIES



compatible to:
CAMLOG®



3,3 mm




4,3 mm



3,8 mm




5,0 mm



**PRECISION
MEETS QUALITY.**

OVERVIEW PLATFORMS:



A: over surface

B: over edge



3,3 mm



3,8 mm




4,3 mm



5,0 mm

 **MADE IN
GERMANY**



3,3 mm
CAMLOG®

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

CAM 80330000



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,5 mm CAM 82330005

GH 1,5 mm CAM 82330015



Scanbody 3D-Guide®

for intra + extra oral scanning

CAM 9.S3D3.300



Lab Analog

CAM 5.La3.300



Implant Pick.Up

for open-tray with integrated screw

CAM TR-00013.3



Abutment Screw

Recommended tightening torques: 20 Ncm

CAM 6033843



DIM-ANALOG®

CAM 5.DIM.330



Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSL.A10

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,3 mm
CAMLOG®

ALL DETAILS
ONLINE



nt-Preform®
incl. final Abutment Screw

CAM 10.Mi3.300.05



Scanbody 3D-Guide®
for intra + extra oral scanning

CAM 9.S3D3.300



nt-Mill Lab Analog
for nt-Preform®

CAM 5.La3.300



colored milling screw
for nt-Preform®

CAM 6.Sc3.843



Abutment Screw
Recommended tightening torques 20 Ncm

CAM 6033843



DIM-ANALOG®

CAM 5.DIM.330



Sleeve (für DIM-ANALOG®)
for plaster model

DIM 5.XSL.A10



Abutment Grip
for nt-Preform® Holder

W12.100.00



Customized Abutment
milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,3 mm
CAMLOG®

ALL DETAILS
ONLINE



2-CONNECT®-Abutment Set (straight)

incl. 2-CONNECT®-Coping Female Part and 2-CONNECT®-Coping-screw

GH 0,5 mm CAM-2CON 8033001S

GH 2,5 mm CAM-2CON 8033002S



2-CONNECT®-Abutment

Male Part, Recommended tightening torques 20 Ncm

GH 0,5 mm CAM-2CON 8033001M

GH 2,5 mm CAM-2CON 8033002M



2-CONNECT® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONNECT®-Coping

Female Part

2CON F-U



2-CONNECT®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 61



2-CONNECT®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONNECT®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONNECT®-lab Analog

for straight Male Parts


2CON 5.La1.800

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



3,3 mm
CAMLOG®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	CAM-2CON 80330C3S
B: over edge	GH 2,5 mm	CAM-2CON 80330A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	CAM-2CON 80330C3M
B: over edge	GH 2,5 mm	CAM-2CON 80330A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

CAM 6.Ms3.843



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY

WWW.NT.DENTAL



3,8 mm
CAMLOG®

ALL DETAILS
ONLINE



Titanium Base (engaged)
incl. Abutmentscrew

CAM 80380000



Titanium Base NX (non engaged)
incl. Abutmentscrew

GH 0,5 mm CAM 82380005

GH 1,5 mm CAM 82380015



Scanbody 3D-Guide®
for intra + extra oral scanning

CAM 9.S3D3.800



Lab Analog

CAM 5.La3.800



Implant Pick.Up
for open-tray with integrated screw

CAM TR-00023.8



Abutment Screw
Recommended tightening torques 20 Ncm

CAM 6033843



DIM-ANALOG®

CAM 5.DIM.380



Sleeve (für DIM-ANALOG®)
for plaster model


DIM 5.XSL.A40

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,8 mm
CAMLOG®

ALL DETAILS
ONLINE



nt-Preform®
incl. final Abutment Screw

GH 0,5 mm CAM 10.Mi3.800.05



Scanbody 3D-Guide®
for intra + extra oral scanning

CAM 9.S3D3.800



nt-Mill Lab Analog
for nt-Preform®

CAM 5.La3.800



colored milling screw
for nt-Preform®

CAM 6.Sc3.843



Abutment Screw
Recommended tightening torques 20 Ncm

CAM 6033843



DIM-ANALOG®

CAM 5.DIM.380



Sleeve (für DIM-ANALOG®)
for plaster model

DIM 5.XSL.A40



Abutment Grip
for nt-Preform® Holder

W12.100.00



Customized Abutment
milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,8 mm
CAMLOG®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm CAM-2CON 8038001S

GH 2,5 mm CAM-2CON 8038002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 30 Ncm

GH 0,5 mm CAM-2CON 8038001M

GH 2,5 mm CAM-2CON 8038002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY

WWW.NT.DENTAL



3,8 mm
CAMLOG®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	CAM-2CON 80380C3S
B: over edge	GH 2,5 mm	CAM-2CON 80380A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	CAM-2CON 80380C3M
B: over edge	GH 2,5 mm	CAM-2CON 80380A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

CAM 6.Ms3.843



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY

WWW.NT.DENTAL



4,3 mm
CAMLOG®

ALL DETAILS
ONLINE



Titanium Base (engaged)
incl. Abutmentscrew

CAM 80430000



Titanium Base NX (non engaged)
incl. Abutmentscrew

GH 0,5 mm CAM 82430005

GH 1,5 mm CAM 82430015



Scanbody 3D-Guide®
for intra + extra oral scanning

CAM 9.S3D4.300



Lab Analog

CAM 5.La4.300



Implant Pick.Up
for open-tray with integrated screw

CAM TR-00034.3



Abutment Screw
Recommended tightening torques 20 Ncm

CAM 6033843



DIM-ANALOG®

CAM 5.DIM.430



Sleeve (für DIM-ANALOG®)
for plaster model


DIM 5.XSL.A40

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,3 mm
CAMLOG®

ALL DETAILS
ONLINE



nt-Preform®
incl. final Abutment Screw

GH 0,5 mm CAM 10.Mi4.300.05



Scanbody 3D-Guide®
for intra + extra oral scanning

CAM 9.S3D4.300



nt-Mill Lab Analog
for nt-Preform®

CAM 5.La4.300



colored milling screw
for nt-Preform®

CAM 6.Sc3.843



Abutment Screw
Recommended tightening torques 20 Ncm

CAM 6033843



DIM-ANALOG®

CAM 5.DIM.430



Sleeve (für DIM-ANALOG®)
for plaster model

DIM 5.XSL.A40



Abutment Grip
for nt-Preform® Holder

W12.200.00



Customized Abutment
milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,3 mm
CAMLOG®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm CAM-2CON 8043001S

GH 2,5 mm CAM-2CON 8043002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 30 Ncm

GH 0,5 mm CAM-2CON 8043001M

GH 2,5 mm CAM-2CON 8043002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY

WWW.NT.DENTAL



4,3 mm
CAMLOG®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface GH 2,5 mm CAM-2CON 80430C3S
B: over edge GH 2,5 mm CAM-2CON 80430A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface GH 2,5 mm CAM-2CON 80430C3M
B: over edge GH 2,5 mm CAM-2CON 80430A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

CAM 6.Ms3.843



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



5,0 mm
CAMLOG®

ALL DETAILS
ONLINE



Titanium Base (engaged)
incl. Abutmentscrew

CAM 80500000



Titanium Base NX (non engaged)
incl. Abutmentscrew

GH 0,5 mm CAM 82500005

GH 1,5 mm CAM 82500015



Scanbody 3D-Guide®
for intra + extra oral scanning

CAM 9.S3D5.060



Lab Analog

CAM 5.La5.060



Implant Pick.Up
for open-tray with integrated screw

CAM TR-00045.0



Abutment Screw
Recommended tightening torques 20 Ncm

CAM 6150600



DIM-ANALOG®

CAM 5.DIM.506



Sleeve (für DIM-ANALOG®)
for plaster model


DIM 5.XSL.A50

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



5,0 mm
CAMLOG®

ALL DETAILS
ONLINE



nt-Preform®
incl. final Abutment Screw

GH 0,5 mm CAM 10.Mi5.000.05



Scanbody 3D-Guide®
for intra + extra oral scanning

CAM 9.S3D5.060



nt-Mill Lab Analog
for nt-Preform®

CAM 5.La5.060



colored milling screw
for nt-Preform®

CAM 6.Sc5.060



Abutment Screw
Recommended tightening torques 20 Ncm

CAM 6150600



DIM-ANALOG®

CAM 5.DIM.506



Sleeve (für DIM-ANALOG®)
for plaster model

DIM 5.XSLA50



Abutment Grip
for nt-Preform® Holder

W12.300.00



Customized Abutment
milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



5,0 mm
CAMLOG®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm CAM-2CON 8050001S

GH 2,5 mm CAM-2CON 8050002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 30 Ncm

GH 0,5 mm CAM-2CON 8050001M

GH 2,5 mm CAM-2CON 8050002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



5,0 mm
CAMLOG®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	CAM-2CON 80500C3S
B: over edge	GH 2,5 mm	CAM-2CON 80500A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	CAM-2CON 80500C3M
B: over edge	GH 2,5 mm	CAM-2CON 80500A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

CAM 6.Ms5.060



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

CAM-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY

WWW.NT.DENTAL

NOTES

[illegible]

EUROPE

NT-Trading GmbH & Co. KG

G.-Braun-Str. 18
76187 Karlsruhe
Germany

E info@nt.dental

W www.nt.dental

NORTH AMERICA

NT-US, Inc.

132 Pepes Farm Road,
Unit 3
Milford, CT 06460
USA

E info@nt.dental (general information)

E order-us@nt.dental (for orders)

W www.nt.dental

2023_NT_Serienbroschüre_CAM_EN_REV002