



E-SERIES




compatible to:
Noble Biocare™ / Replace Select®



3,5 mm NP




5,0 mm WP



4,3 mm RP




6,0 mm



**PRECISION
MEETS QUALITY.**

**OVERVIEW
PLATFORMS:**



A: over surface

B: over edge



3,5 mm NP



4,3 mm RP




5,0 mm WP



6,0 mm









 **MADE IN
GERMANY**




3,5 mm NP
Noble Biocare™ / Replace Select®

ALL DETAILS
ONLINE












	Titanium Base (engaged) incl. Abutmentscrew	E 800
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 0,4 mm E 82350004 GH 1,4 mm E 82350014
	Scanbody 3D-Guide® for intra + extra oral scanning	E 9.S3D3.500
	Lab Analog	E 50
	Implant Pick.Up for open-tray with integrated screw	E TR-NP013.5
	Abutment Screw Recommended tightening torques 35 Ncm	E 60
	DIM-ANALOG®	E 5.DIM.350
	Sleeve (für DIM-ANALOG®) for plaster model	DIM 5.XSL.A10




3,5 mm NP
Noble Biocare™ / Replace Select®

ALL DETAILS
ONLINE



	nt-Preform® incl. final Abutment Screw	E 10.Mi3.500.05
	Scanbody 3D-Guide® for intra + extra oral scanning	E 9.S3D3.500
	nt-Mill Lab Analog for nt-Preform®	E 5.La3.500
	colored milling screw for nt-Preform®	E 6.Sc3.500
	Abutment Screw Recommended tightening torques 35 Ncm	E 60
	DIM-ANALOG®	E 5.DIM.350
	Sleeve (für DIM-ANALOG®) for plaster model	DIM 5.XSLA10
	Abutment Grip for nt-Preform® Holder	W12.100.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1



3,5 mm NP
Noble Biocare™ / Replace Select®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)
incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm E 800 S
GH 2,5 mm E-2CON 8035002S



2-CONnect®-Abutment
Male Part, Recommended tightening torques 35 Ncm

GH 0,5 mm E 800 M
GH 2,5 mm E-2CON 8035002M



2-CONnect® Scanbody 3D-Guide®
for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING
For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping
Female Part

2CON F-U



2-CONnect®-Coping-Screw
Recommended tightening torques 30 Ncm

KS 61



2-CONnect®-DIM-ANALOG®
for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)
for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up
for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog
for straight Male Parts


2CON 5.La1.800

E-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,5 mm NP
Noble Biocare™ / Replace Select®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface GH 2,5 mm E-2CON 80350C3S
B: over edge GH 2,5 mm E-2CON 80350A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface GH 2,5 mm E-2CON 80350C3M
B: over edge GH 2,5 mm E-2CON 80350A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

E 6.Ms3.500



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw


2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts









2CON 5.La1.600




4,3 mm RP
Noble Biocare™ / Replace Select®










ALL DETAILS
ONLINE




	Titanium Base (engaged) incl. Abutmentscrew	E 810
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 0,4 mm E 82430004 GH 1,4 mm E 82430014
	Scanbody 3D-Guide® for intra + extra oral scanning	E 9.S3D4.300
	Lab Analog	E 51
	Implant Pick.Up for open-tray with integrated screw	E TR-RP024.3
	Abutment Screw Recommended tightening torques 35 Ncm	E 61
	DIM-ANALOG®	E 5.DIM.430
	Sleeve (für DIM-ANALOG®) for plaster model	DIM 5.XSL.A20



4,3 mm RP
Noble Biocare™ / Replace Select®

	nt-Preform® incl. final Abutment Screw	GH 0,5 mm E 10.Mi4.300.05
	Scanbody 3D-Guide® for intra + extra oral scanning	E 9.S3D4.300
	nt-Mill Lab Analog for nt-Preform®	E 5.La4.300
	colored milling screw for nt-Preform®	E 6.Sc4.350
	Abutment Screw Recommended tightening torques 35 Ncm	E 61
	DIM-ANALOG®	E 5.DIM.430
	Sleeve (für DIM-ANALOG®) for plaster model	DIM 5.XSLA20
	Abutment Grip for nt-Preform® Holder	W12.200.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1



4,3 mm RP

Noble Biocare™ / Replace Select®

**2-CONnect®-Abutment Set (straight)**

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm E 810 S

GH 2,5 mm E-2CON 8043002S

**2-CONnect®-Abutment**

Male Part, Recommended tightening torques 35 Ncm

GH 0,5 mm E 810 M

GH 2,5 mm E-2CON 8043002M

**2-CONnect® Scanbody 3D-Guide®**

for intra + extra oral scanning

2CON 9.S3D1.800

**LOcON COPING**

For more information, see our LOcON brochure

LOC 852001

**2-CONnect®-Coping**

Female Part

2CON F-U

**2-CONnect®-Coping-Screw**

Recommended tightening torques 30 Ncm

KS 61

**2-CONnect®-DIM-ANALOG®**

for straight Male Part, KS 61

2CON 5.DIM180

**Sleeve (for DIM-ANALOG®)**

for plaster model

DIM 5.XSL.A30

**2-CONnect®-Implant pick-up**


for open-tray with integrated screw

2CON TR-00021.8

**2-CONnect®-lab Analog**

for straight Male Parts

2CON 5.La1.800



4,3 mm RP

Noble Biocare™ / Replace Select®

**2-CONnect®-Abutment Set (angled)**incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	E-2CON 80430C3S
B: over edge	GH 2,5 mm	E-2CON 80430A3S

**2-CONnect®-Abutment (angled)**

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	E-2CON 80430C3M
B: over edge	GH 2,5 mm	E-2CON 80430A3M

**2-CONnect® Scanbody 3D-Guide®**

for intra + extra oral scanning

2CON 9.S3D1.600

**2-CONnect®-Coping**

Female Part

2CON F-U

**2-CONnect®-Coping-Screw**

Recommended tightening torques 15 Ncm

KS 62

**Male Part Screw**for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

E 6.Ms4.350

**2-CONnect®-DIM-ANALOG®**

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160

**Sleeve (for DIM-ANALOG®)**

for plaster model

DIM 5.XSL.A30

**2-CONnect®-Implant pick-up**


for open-tray with integrated screw

2CON TR-00031.6

**2-CONnect®-lab Analog**

for angled 2-CONnect® Abutment Male Parts









2CON 5.La1.600



5,0 mm WP
Noble Biocare™ / Replace Select®

ALL DETAILS
ONLINE




	Titanium Base (engaged) incl. Abutmentscrew	E 820
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 0,4 mm E 82500004 GH 1,4 mm E 82500014
	Scanbody 3D-Guide® for intra + extra oral scanning	E 9.S3D5.000
	Lab Analog	E 52
	Implant Pick.Up for open-tray with integrated screw	E TR-WP035.0
	Abutment Screw Recommended tightening torques 35 Ncm	E 61
	DIM-ANALOG®	E 5.DIM.500
	Sleeve (für DIM-ANALOG®) for plaster model	DIM 5.XSL.A30

E-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY










WWW.NT.DENTAL



5,0 mm WP
Noble Biocare™ / Replace Select®

ALL DETAILS
ONLINE




	nt-Preform® incl. final Abutment Screw	GH 0,5 mm E 10.Mi5.000.05
	Scanbody 3D-Guide® for intra + extra oral scanning	E 9.S3D5.000
	nt-Mill Lab Analog for nt-Preform®	E 5.La5.000
	colored milling screw for nt-Preform®	E 6.Sc4.350
	Abutment Screw Recommended tightening torques 35 Ncm	E 61
	DIM-ANALOG®	E 5.DIM.500
	Sleeve (für DIM-ANALOG®) for plaster model	DIM 5.XSLA30
	Abutment Grip for nt-Preform® Holder	W12.300.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

E-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



5,0 mm WP
Noble Biocare™ / Replace Select®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm E 820 S

GH 2,5 mm E-2CON 8050002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 35 Ncm

GH 0,5 mm E 820 M

GH 2,5 mm E-2CON 8050002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 30 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

E-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



5,0 mm WP
Noble Biocare™ / Replace Select®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface GH 2,5 mm E-2CON 80500C3S
B: over edge GH 2,5 mm E-2CON 80500A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw , mounting aid, angled 20°

A: over surface GH 2,5 mm E-2CON 80500C3M
B: over edge GH 2,5 mm E-2CON 80500A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

E 6.Ms4.350



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw


2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600










6,0 mm

Noble Biocare™ / Replace Select®

ALL DETAILS
ONLINE



	Titanium Base (engaged) incl. Abutmentscrew	E 830
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 0,4 mm E 82600004 GH 1,4 mm E 82600014
	Scanbody 3D-Guide® for intra + extra oral scanning	E 9.S3D6.000
	Lab Analog	E 53
	Abutment Screw Recommended tightening torques: 35 Ncm	E 61
	DIM-ANALOG®	E 5.DIM.600
	Sleeve (für DIM-ANALOG®) for plaster model	DIM 5.XSLA60

EUROPE

NT-Trading GmbH & Co. KG

G.-Braun-Str. 18
76187 Karlsruhe
Germany

E info@nt.dental

W www.nt.dental

NORTH AMERICA

NT-US, Inc.

132 Pepes Farm Road,
Unit 3
Milford, CT 06460
USA

E info@nt.dental (general information)

E order-us@nt.dental (for orders)

W www.nt.dental

2023_NT_Serienbroschüre_E_EN_REV001