

S-SERIES




compatible to:
Dentsply Implants® / Astra Tech
Implant System OsseoSpeed® TX



3,5 / 4,0 mm




4,5 / 5,0 mm



**PRECISION
MEETS QUALITY.**

**OVERVIEW
PLATFORMS:**




3,5 / 4,0 mm



4,5 / 5,0 mm

 **MADE IN
GERMANY**



3,5 / 4,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

S 800



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,8 mm S 82354008

GH 1,8 mm S 82354018



Scanbody 3D-Guide®

for intra + extra oral scanning

S 9.S3D3.540



Lab Analog

S 50



Implant Pick.Up

for open-tray with integrated screw

S TR-00013.5



Abutment Screw

Recommended tightening torques 25 Ncm

S 60



DIM-ANALOG®

S 5.DIM.354



Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSL.A10

S-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,5 / 4,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS
ONLINE



nt-Preform®

incl. final Abutment Screw

GH 0,5 mm S 10.Mi3.540.05



Scanbody 3D-Guide®

for intra + extra oral scanning

S 9.S3D3.540



nt-Mill Lab Analog

for nt-Preform®

S 5.La3.540



colored milling screw

for nt-Preform®

S 6.Sc3.540



Abutment Screw

Recommended tightening torques 25 Ncm

S 60



DIM-ANALOG®

S 5.DIM.354



Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSLA10



Abutment Grip

for nt-Preform® Holder

W12.100.00



Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

S-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,5 / 4,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm S 800 S

GH 2,5 mm S-2CON 8035402S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 25 Ncm

GH 0,5 mm S 800 M

GH 2,5 mm S-2CON 8035402M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

S-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



3,5 / 4,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm S-2CON 80354A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm S-2CON 80354A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

S 6.Ms3.540



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

S-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



4,5 / 5,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

S 820



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,8 mm S 82455008

GH 1,8 mm S 82455018



Scanbody 3D-Guide®

for intra + extra oral scanning

S 9.S3D4.550



Lab Analog

S 52



Implant Pick.Up

for open-tray with integrated screw

S TR-00024.5



Abutment Screw

Recommended tightening torques 25 Ncm

S 61



DIM-ANALOG®

S 5.DIM.455



Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSL.A20

S-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,5 / 5,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS
ONLINE



nt-Preform®

incl. final Abutment Screw

GH 0,5 mm S 10.Mi4.550.05



Scanbody 3D-Guide®

for intra + extra oral scanning

S 9.S3D4.550



nt-Mill Lab Analog

for nt-Preform®

S 5.La4.550



colored milling screw

for nt-Preform®

S 6.Sc4.550



Abutment Screw

Recommended tightening torques: 25 Ncm

S 61



DIM-ANALOG®

S 5.DIM.455



Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSL.A20



Abutment Grip

for nt-Preform® Holder

W12.200.00



Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

S-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,5 / 5,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm S 820 S

GH 2,5 mm S-2CON 8045502S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 25 Ncm

GH 0,5 mm S 820 M

GH 2,5 mm S-2CON8045502M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

S-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



4,5 / 5,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm S-2CON 80455A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm S-2CON 80455A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

S 6.Ms4550



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

S-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL

NOTES

[illegible]

EUROPE

NT-Trading GmbH & Co. KG

G.-Braun-Str. 18
76187 Karlsruhe
Germany

E info@nt.dental

W www.nt.dental

NORTH AMERICA

NT-US, Inc.

132 Pepes Farm Road,
Unit 3
Milford, CT 06460
USA

E info@nt.dental (general information)

E order-us@nt.dental (for orders)

W www.nt.dental

2023_NT_Serienbroschüre_S_EN_REV001