


SEV -SERIES



compatible to:
Dentsply Implants® / Astra Tech
Implant System™ EV



3,0 mm



4,8 mm




3,6 mm



5,4 mm




4,2 mm



**PRECISION
MEETS QUALITY.**

**OVERVIEW
PLATFORMS:**



3,0 mm



3,6 mm



4,2 mm




4,8 mm



5,4 mm

 **MADE IN
GERMANY**



3,0 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

SEV 80300000



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 1,0 mm SEV 82300010

GH 2,0 mm SEV 82300020



Scanbody 3D-Guide®

for intra + extra oral scanning

SEV 9.S3D3.000



Lab Analog

SEV 5.La3.000



Implant Pick.Up

for open-tray with integrated screw

SEV TR-00013.0



Abutment Screw

Recommended tightening torques 25 Ncm

SEV 603000



DIM-ANALOG®

SEV 5.DIM.300



Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSL.A10

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,0 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



nt-Preform®

incl. final Abutment Screw

GH 0,5 mm SEV 10.Mi3.000.05



Scanbody 3D-Guide®

for intra + extra oral scanning

SEV 9.S3D3.000



nt-Mill Lab Analog

for nt-Preform®

SEV 5.La3.000



colored milling screw

for nt-Preform®

SEV 6.Sc3.000



Abutment Screw

Recommended tightening torques 25 Ncm

SEV 603000



DIM-ANALOG®

SEV 5.DIM.300



Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSL.A10



Abutment Grip

for nt-Preform® Holder

W12.100.00



Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,0 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm SEV-2CON 8030001S

GH 2,5 mm SEV-2CON 8030002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 25 Ncm

GH 0,5 mm SEV-2CON 8030001M

GH 2,5 mm SEV-2CON 8030002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,0 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm SEV-2CON 80300A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm SEV-2CON 80300A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

SEV 6.Ms3.000



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



3,6 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

SEV 80360000



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 1,0 mm SEV 82360010

GH 2,0 mm SEV 82360020



Scanbody 3D-Guide®

for intra + extra oral scanning

SEV 9.S3D3.600



Lab Analog

SEV 5.La3.600



Implant Pick.Up

for open-tray with integrated screw

SEV TR-00023.6



Abutment Screw

Recommended tightening torques 25 Ncm

SEV 603600



DIM-ANALOG®

SEV 5.DIM.360



Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSL.A10

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,6 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



nt-Preform®

incl. final Abutment Screw

GH 0,5 mm SEV 10.Mi3.600.05



Scanbody 3D-Guide®

for intra + extra oral scanning

SEV 9.S3D3.600



nt-Mill Lab Analog

for nt-Preform®

SEV 5.La3.600



colored milling screw

for nt-Preform®

SEV 6.Sc3.600



Abutment Screw

Recommended tightening torques 25 Ncm

SEV 603600



DIM-ANALOG®

SEV 5.DIM.360



Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSL.A10



Abutment Grip

for nt-Preform® Holder

W12.100.00



Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



3,6 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm SEV-2CON 8036001S

GH 2,5 mm SEV-2CON 8036002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 25 Ncm

GH 0,5 mm SEV-2CON 8036001M

GH 2,5 mm SEV-2CON 8036002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



3,6 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm SEV-2CON 80360A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm SEV-2CON 80360A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

SEV 6.Ms3.600



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



4,2 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

SEV 80420000



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 1,0 mm SEV 82420010

GH 2,0 mm SEV 82420020



Scanbody 3D-Guide®

for intra + extra oral scanning

SEV 9.S3D4.200



Lab Analog

SEV 5.La4.200



Implant Pick.Up

for open-tray with integrated screw

SEV TR-00034.2



Abutment Screw

Recommended tightening torques 25 Ncm

SEV 604200



DIM-ANALOG®

SEV 5.DIM.420



Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSL.A20

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,2 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



nt-Preform®

incl. final Abutment Screw

GH 0,5 mm SEV 10.Mi4.200.05



Scanbody 3D-Guide®

for intra + extra oral scanning

SEV 9.S3D4.200



nt-Mill Lab Analog

for nt-Preform®

SEV 5.La4.200



colored milling screw

for nt-Preform®

SEV 6.Sc4.200



Abutment Screw

Recommended tightening torques 25 Ncm

SEV 604200



DIM-ANALOG®

SEV 5.DIM.420



Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSL.A20



Abutment Grip

for nt-Preform® Holder

W12.200.00



Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,2 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm SEV-2CON 8042001S

GH 2,5 mm SEV-2CON 8042002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 25 Ncm

GH 0,5 mm SEV-2CON 8042001M

GH 2,5 mm SEV-2CON 8042002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,2 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm SEV-2CON 80420A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm SEV-2CON 80420A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

SEV 6.Ms4.200



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



4,8 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

SEV 80480000



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 1,0 mm SEV 82480010

GH 2,0 mm SEV 82480020



Scanbody 3D-Guide®

for intra + extra oral scanning

SEV 9.S3D4.800



Lab Analog

SEV 5.La4.800



Implant Pick.Up

for open-tray with integrated screw

SEV TR-00044.8



Abutment Screw

Recommended tightening torques 25 Ncm

SEV 604854



DIM-ANALOG®

SEV 5.DIM.480



Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSL.A30

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,8 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



nt-Preform®

incl. final Abutment Screw

GH 0,5 mm SEV 10.Mi4.800.05



Scanbody 3D-Guide®

for intra + extra oral scanning

SEV 9.S3D4.800



nt-Mill Lab Analog

for nt-Preform®

SEV 5.La4.800



colored milling screw

for nt-Preform®

SEV 6.Sc4.854



Abutment Screw

Recommended tightening torques: 25 Ncm

SEV 604854



DIM-ANALOG®

SEV 5.DIM.480



Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSLA30



Abutment Grip

for nt-Preform® Holder

W12.300.00



Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,8 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm SEV-2CON 8048001S

GH 2,5 mm SEV-2CON 8048002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 25 Ncm

GH 0,5 mm SEV-2CON 8048001M

GH 2,5 mm SEV-2CON 8048002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



4,8 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm SEV-2CON 80480A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm SEV-2CON 80480A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

SEV 6.Ms4.854



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



5,4 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

SEV 80540000



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 1,0 mm SEV 82540010

GH 2,0 mm SEV 82540020



Scanbody 3D-Guide®

for intra + extra oral scanning

SEV 9.S3D5.400



Lab Analog

SEV 5.La5.400



Implant Pick.Up

for open-tray with integrated screw

SEV TR-00055.4



Abutment Screw

Recommended tightening torques 25 Ncm

SEV 604854



DIM-ANALOG®

SEV 5.DIM.540



Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSL.A30

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



5,4 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



nt-Preform®

incl. final Abutment Screw

GH 0,5 mm SEV 10.Mi5.400.05



Scanbody 3D-Guide®

for intra + extra oral scanning

SEV 9.S3D5.400



nt-Mill Lab Analog

for nt-Preform®

SEV 5.La5.400



colored milling screw

for nt-Preform®

SEV 6.Sc4.854



Abutment Screw

Recommended tightening torques 25 Ncm

SEV 604854



DIM-ANALOG®

SEV 5.DIM.540



Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



Abutment Grip

for nt-Preform® Holder

W12.300.00



Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



5,4 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm SEV-2CON 8054001S

GH 2,5 mm SEV-2CON 8054002S



2-CONnect®-Abutment

Male Part, Recommended tightening torques 25 Ncm

GH 0,5 mm SEV-2CON 8054001M

GH 2,5 mm SEV-2CON 8054002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



5,4 mm

Dentsply Implants® / Astra Tech Implant System™ EV

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm SEV-2CON 80540A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm SEV-2CON 80540A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20
Ncm

SEV 6.Ms4.854



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

SEV-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL

NOTES

[illegible]

EUROPE

NT-Trading GmbH & Co. KG

G.-Braun-Str. 18
76187 Karlsruhe
Germany

E info@nt.dental

W www.nt.dental

NORTH AMERICA

NT-US, Inc.

132 Pepes Farm Road,
Unit 3
Milford, CT 06460
USA

E info@nt.dental (general information)

E order-us@nt.dental (for orders)

W www.nt.dental

2023_NT_Serienbroschüre_SEV_EN_REV001