

# BEG -SERIES



compatible to:  
BEGO Implant Systems®  
Semados®, S, RS-RX, RI-Line



3,25 / 3,75 mm




4,5 mm



4,1 mm




5,5 mm



**PRECISION  
MEETS QUALITY.**

**OVERVIEW  
PLATFORMS:**



3,25 / 3,75 mm

---



4,1 mm

---



4,5 mm


---



5,5 mm

---

 **MADE IN  
GERMANY**



3,25 / 3,75 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### Titanium Base (engaged)

incl. Abutmentscrew

BEG 80323700



### Scanbody 3D-Guide®

for intra + extra oral scanning

BEG 9.S3D3.237



### Lab Analog

BEG 5.La3.237



### Implant Pick.Up

for open-tray with integrated screw

BEG TR-00013.2



### Abutment Screw

Recommended tightening torques 30 Ncm

BEG 603255



### DIM-ANALOG®

BEG 5.DIM.323



### Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSLA20

# BEG-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



3,25 / 3,75 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



**nt-Preform®**

incl. final Abutment Screw

GH 0,5 mm [BEG 10.MI3.237.05](#)



**Scanbody 3D-Guide®**

for intra + extra oral scanning

[BEG 9.S3D3.237](#)



**nt-Mill Lab Analog**

for nt-Preform®

[BEG 5.La3.237](#)



**colored milling screw**

for nt-Preform®

[BEG 6.Sc3.255](#)



**Abutment Screw**

Recommended tightening torques 30 Ncm

[BEG 603255](#)



**DIM-ANALOG®**

[BEG 5.DIM.323](#)



**Sleeve (für DIM-ANALOG®)**

for plaster model

[DIM 5.XSL.A20](#)



**Abutment Grip**

for nt-Preform® Holder

[W 12.100.000](#)



**Customized Abutment**

milled by NT-iDent / Not available in the USA


[TZ15.TiO.Si1](#)

**BEG-SERIES**

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



3,25 / 3,75 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm BEG-2CON 8032371S

GH 2,5 mm BEG-2CON 8032372S



### 2-CONnect®-Abutment

Male Part, Recommended tightening torques 30 Ncm

GH 0,5 mm BEG-2CON 8032371M

GH 2,5 mm BEG-2CON 8032372M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

# BEG-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



3,25 / 3,75 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm [BEG-2CON 80323A3S](#)



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm [BEG-2CON 80323A3M](#)



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

[2CON 9.S3D1.600](#)



### 2-CONnect®-Coping

Female Part

[2CON F-U](#)



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

[KS 62](#)



### Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20  
Ncm

[BEG 6.Ms3.255](#)



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

[2CON 5.DIM160](#)



### Sleeve (for DIM-ANALOG®)

for plaster model

[DIM 5.XSL.A30](#)



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

[2CON TR-00031.6](#)



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


[2CON 5.La1.600](#)

# BEG-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

[WWW.NT.DENTAL](http://WWW.NT.DENTAL)



4,1 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### **Titanium Base (engaged)**

incl. Abutmentscrew

BEG 80410000



### **Scanbody 3D-Guide®**

for intra + extra oral scanning

BEG 9.S3D4.100



### **Lab Analog**

BEG 5.La4.100



### **Implant Pick.Up**

for open-tray with integrated screw

BEG TR-00024.2



### **Abutment Screw**

Recommended tightening torques 30 Ncm

BEG 603255



### **DIM-ANALOG®**

BEG 5.DIM.410



### **Sleeve (für DIM-ANALOG®)**

for plaster model


DIM 5.XSLA20

# BEG-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



4,1 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### nt-Preform®

incl. final Abutment Screw

GH 0,5 mm BEG 10.Mi4.100.05



### Scanbody 3D-Guide®

for intra + extra oral scanning

BEG 9.S3D4.100



### nt-Mill Lab Analog

for nt-Preform®

BEG 5.La4.100



### colored milling screw

for nt-Preform®

BEG 6.Sc3.255



### Abutment Screw

Recommended tightening torques 30 Ncm

BEG 603255



### DIM-ANALOG®

BEG 5.DIM.410



### Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSL.A20



### Abutment Grip

for nt-Preform® Holder

W 12.200.000



### Customized Abutment

milled by NT-iDent / Not available in the USA

TZ15.TiO.Si1


# BEG-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**





4,1 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm BEG-2CON 8041001S

GH 2,5 mm BEG-2CON 8041002S



### 2-CONnect®-Abutment

Male Part, Recommended tightening torques 30 Ncm

GH 0,5 mm BEG-2CON 8041001M

GH 2,5 mm BEG-2CON 8041002M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

# BEG-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



4,1 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm [BEG-2CON 80410A3S](#)



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm [BEG-2CON 80410A3M](#)



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

[2CON 9.S3D1.600](#)



### 2-CONnect®-Coping

Female Part

[2CON F-U](#)



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

[KS 62](#)



### Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20  
Ncm

[BEG 6.Ms3.255](#)



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

[2CON 5.DIM160](#)



### Sleeve (for DIM-ANALOG®)

for plaster model

[DIM 5.XSL.A30](#)



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

[2CON TR-00031.6](#)



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


[2CON 5.La1.600](#)

# BEG-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

[WWW.NT.DENTAL](http://WWW.NT.DENTAL)



4,5 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### **Titanium Base (engaged)**

incl. Abutmentscrew

BEG 80450000



### **Scanbody 3D-Guide®**

for intra + extra oral scanning

BEG 9.S3D4.500



### **Lab Analog**

BEG 5.La4.500



### **Implant Pick.Up**

for open-tray with integrated screw

BEG TR-00034.5



### **Abutment Screw**

Recommended tightening torques 30 Ncm

BEG 603255



### **DIM-ANALOG®**

BEG 5.DIM.450



### **Sleeve (für DIM-ANALOG®)**

for plaster model


DIM 5.XSLA30

# BEG-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



4,5 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### nt-Preform®

incl. final Abutment Screw

GH 0,5 mm BEG 10.MI4.500.05



### Scanbody 3D-Guide®

for intra + extra oral scanning

BEG 9.S3D4.500



### nt-Mill Lab Analog

for nt-Preform®

BEG 5.La4.500



### colored milling screw

for nt-Preform®

BEG 6.Sc3.255



### Abutment Screw

Recommended tightening torques 30 Ncm

BEG 603255



### DIM-ANALOG®

BEG 5.DIM.450



### Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSLA30



### Abutment Grip

for nt-Preform® Holder

W 12.300.000



### Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

# BEG-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



4,5 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm BEG-2CON 8045001S

GH 2,5 mm BEG-2CON 8045002S



### 2-CONnect®-Abutment

Male Part, Recommended tightening torques 30 Ncm

GH 0,5 mm BEG-2CON 8045001M

GH 2,5 mm BEG-2CON 8045002M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

# BEG-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



4,5 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm [BEG-2CON 80450A3S](#)



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm [BEG-2CON 80450A3M](#)



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

[2CON 9.S3D1.600](#)



### 2-CONnect®-Coping

Female Part

[2CON F-U](#)



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

[KS 62](#)



### Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20  
Ncm

[BEG 6.Ms3.255](#)



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

[2CON 5.DIM160](#)



### Sleeve (for DIM-ANALOG®)

for plaster model

[DIM 5.XSL.A30](#)



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

[2CON TR-00031.6](#)



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


[2CON 5.La1.600](#)

# BEG-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

[WWW.NT.DENTAL](http://WWW.NT.DENTAL)



5,5 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### Titanium Base (engaged)

incl. Abutmentscrew

BEG 80550000



### Scanbody 3D-Guide®

for intra + extra oral scanning

BEG 9.S3D5.500



### Lab Analog

BEG 5.La5.500



### Implant Pick.Up

for open-tray with integrated screw

BEG TR-00045.5



### Abutment Screw

Recommended tightening torques 30 Ncm

BEG 603255



### DIM-ANALOG®

BEG 5.DIM.550



### Sleeve (für DIM-ANALOG®)

for plaster model


DIM 5.XSLA30

# BEG-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



5,5 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### nt-Preform®

incl. final Abutment Screw

GH 0,5 mm BEG 10.Mi5.500.05



### Scanbody 3D-Guide®

for intra + extra oral scanning

BEG 9.S3D5.500



### nt-Mill Lab Analog

for nt-Preform®

BEG 5.La5.500



### colored milling screw

for nt-Preform®

BEG 6.Sc3.255



### Abutment Screw

Recommended tightening torques 30 Ncm

BEG 603255



### DIM-ANALOG®

BEG 5.DIM.550



### Sleeve (für DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### Abutment Grip

for nt-Preform® Holder

W 12.300.000



### Customized Abutment

milled by NT-iDent / Not available in the USA

TZ15.TiO.Si1


# BEG-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**





5,5 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm BEG-2CON 8055001S

GH 2,5 mm BEG-2CON 8055002S



### 2-CONnect®-Abutment

Male Part, Recommended tightening torques 30 Ncm

GH 0,5 mm BEG-2CON 8055001M

GH 2,5 mm BEG-2CON 8055002M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

# BEG-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



5,5 mm

BEGO Implant Systems® Semados® S, RS-RX, RI-Line

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm BEG-2CON 80550A3S



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm BEG-2CON 80550A3M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



### Male Part Screw

for angled 2-CONnect® Abutment, Recommended tightening torques 20  
Ncm

BEG 6.Ms3.255



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

# BEG-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

MADE IN  
GERMANY

WWW.NT.DENTAL

## NOTIZEN

[illegible]

## EUROPE

### **NT-Trading GmbH & Co. KG**

G.-Braun-Str. 18  
76187 Karlsruhe  
Germany

**E** [info@nt.dental](mailto:info@nt.dental)

**W** [www.nt.dental](http://www.nt.dental)

## NORTH AMERICA

### **NT-US, Inc.**

132 Pepes Farm Road,  
Unit 3  
Milford, CT 06460  
USA

**E** [info@nt.dental](mailto:info@nt.dental) (general information)

**E** [order-us@nt.dental](mailto:order-us@nt.dental) (for orders)

**W** [www.nt.dental](http://www.nt.dental)

2023\_NT\_Serienbroschüre\_BEG\_EN\_REV001