

COL -SERIES



compatible to:
CONELOG®




3,3 mm



5,0 mm




3,8 / 4,3 mm



**PRECISION
MEETS QUALITY.**

OVERVIEW PLATFORMS:



A: over surface

B: over edge



3,3 mm




3,8 / 4,3 mm



5,0 mm

 **MADE IN
GERMANY**



3,3 mm
CONELOG®

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrews

GH 0,8 mm COL 80330008

GH 2,0 mm COL 80330020



Titanium Base NX (non engaged)

incl. Abutmentscrews

GH 1,0 mm COL 82330010

GH 2,0 mm COL 82330020



Scanbody 3D-Guide®

for intra + extra oral scanning

COL 9.S3D3.300



Lab Analog

COL 5.La3.300



Implant Pick.Up

for open-tray with integrated screw

COL TR-00013.3



Coloured Lab Screw

COL 623343



Abutment Screw

Recommended tightening torques: 20 Ncm

COL 603343



DIM-ANALOG®

COL 5.DIM.330



Sleeve (for DIM-ANALOG®)

for plaster model


DIM 5.XSLA10

COL-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**










WWW.NT.DENTAL



3,3 mm
CONELOG®

ALL DETAILS
ONLINE




	nt-Preform® incl. final Abutment Screw	COL 10.Mi3.300.05
	Scanbody 3D-Guide® for intra + extra oral scanning	COL 9.S3D3.300
	nt-Mill Lab Analog for nt-Preform®	COL 5.La3.300
	colored milling screw for nt-Preform®	COL 6.Sc3.343
	Abutment Screw Recommended tightening torques 20 Ncm	COL 603343
	DIM-ANALOG®	COL 5.DIM.330
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSLA10
	Abutment Grip for nt-Preform® Holder	W12.100.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

COL-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY

WWW.NT.DENTAL



3,3 mm

CONELOG®

ALL DETAILS
ONLINE



2-CONNECT®-Abutment Set (straight)

incl. 2-CONNECT®-Coping Female Part and 2-CONNECT®-Coping-screw

GH 0,5 mm COL-2CON 8033001S

GH 2,5 mm COL-2CON 8033002S



2-CONNECT®-Abutment

Male Part, recommended tightening torques 20 Ncm

GH 0,5 mm COL-2CON 8033001M

GH 2,5 mm COL-2CON 8033002M



2-CONNECT® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONNECT®-Coping

Female Part

2CON F-U



2-CONNECT®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 61



2-CONNECT®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONNECT®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONNECT®-lab Analog

for straight Male Parts


2CON 5.La1.800

COL-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,3 mm

CONELOG®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	COL-2CON 80330C3S
B: over edge	GH 2,5 mm	COL-2CON 80330A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	COL-2CON 80330C3M
B: over edge	GH 2,5 mm	COL-2CON 80330A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

COL 6.Ms3.343



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

COL-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



3,8 / 4,3 mm
CONELOG®

ALL DETAILS
ONLINE



Titanium Base (engaged)
incl. Abutmentscrew

GH 0,8 mm COL 80384308
GH 2,0 mm COL 80384320



Titanium Base NX (non engaged)
incl. Abutmentscrew

GH 1,0 mm COL 82384310
GH 2,0 mm COL 82384320



Scanbody 3D-Guide®
for intra + extra oral scanning

COL 9.S3D3.843



Lab Analog

COL 5.La3.843



Implant Pick.Up
for open-tray with integrated screw

COL TR-00023.8



Coloured Lab Screw

COL 623343



Abutment Screw
Recommended tightening torques 20 Ncm

COL 603343



DIM-ANALOG®

COL 5.DIM.384



Sleeve (for DIM-ANALOG®)
for plaster model


DIM 5.XSL.A20

COL-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,8 / 4,3 mm
CONELOG®

ALL DETAILS
ONLINE



nt-Preform®
incl. final Abutment Screw

GH 0,5 mm COL 10.Mi3.843.05



Scanbody 3D-Guide®
for intra + extra oral scanning

COL 9.S3D3.843



nt-Mill Lab Analog
for nt-Preform®

COL 5.La3.843



colored milling screw
for nt-Preform®

COL 6.Sc3.343



Abutment Screw
Recommended tightening torques 20 Ncm

COL 603343



DIM-ANALOG®

COL 5.DIM.384



Sleeve (for DIM-ANALOG®)
for plaster model

DIM 5.XSLA20



Abutment Grip
for nt-Preform® Holder

W12.200.00



Customized Abutment
milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

COL-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,8 / 4,3 mm
CONELOG®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm COL-2CON 8038431S

GH 2,5 mm COL-2CON 8038432S



2-CONnect®-Abutment

Male Part, recommended tightening torques 30 Ncm

GH 0,5 mm COL-2CON 8038431M

GH 2,5 mm COL-2CON 8038432M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

COL-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,8 / 4,3 mm
CONELOG®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	COL-2CON 80384C3S
B: over edge	GH 2,5 mm	COL-2CON 80384A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	COL-2CON 80384C3M
B: over edge	GH 2,5 mm	COL-2CON 80384A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

COL 6.Ms3.343



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

COL-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY

WWW.NT.DENTAL



5,0 mm
CONELOG®

ALL DETAILS
ONLINE



Titanium Base (engaged)
incl. Abutmentscrew

GH 0,8 mm COL 80500008
GH 2,0 mm COL 80500020



Titanium Base NX (non engaged)
incl. Abutmentscrew

GH 1,0 mm COL 82500010
GH 2,0 mm COL 82500020



Scanbody 3D-Guide®
for intra + extra oral scanning

COL 9.S3D5.000



Lab Analog

COL 5.La5.000



Coloured Lab Screw

COL 625000



Abutment Screw
Recommended tightening torques 20 Ncm

COL 605000



DIM-ANALOG®

COL 5.DIM.500



Sleeve (for DIM-ANALOG®)
for plaster model


DIM 5.XSL.A30

COL-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**










WWW.NT.DENTAL



5,0 mm
CONELOG®

ALL DETAILS
ONLINE




	nt-Preform® incl. final Abutment Screw	GH 0,5 mm COL 10.Mi5.000.05
	Scanbody 3D-Guide® for intra + extra oral scanning	COL 9.S3D5.000
	nt-Mill Lab Analog for nt-Preform®	COL 5.La5.000
	colored milling screw for nt-Preform®	COL 6.Sc5.000
	Abutment Screw Recommended tightening torques 20 Ncm	COL 605000
	DIM-ANALOG®	COL 5.DIM.500
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSLA30
	Abutment Grip for nt-Preform® Holder	W12.300.000
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

COL-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



5,0 mm
CONELOG®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm COL-2CON 8050001S

GH 2,5 mm COL-2CON 8050002S



2-CONnect®-Abutment

Male Part, recommended tightening torques 30 Ncm

GH 0,5 mm COL-2CON 8050001M

GH 2,5 mm COL-2CON 8050002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

COL-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



5,0 mm
CONELOG®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	COL-2CON 80500C3S
B: over edge	GH 2,5 mm	COL-2CON 80500A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	COL-2CON 80500C3M
B: over edge	GH 2,5 mm	COL-2CON 80500A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

COL 6.Ms5.000



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

COL-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL

NOTES

[illegible]

EUROPE

NT-Trading GmbH & Co. KG

G.-Braun-Str. 18
76187 Karlsruhe
Germany

E info@nt.dental

W www.nt.dental

NORTH AMERICA

NT-US, Inc.

132 Pepes Farm Road,
Unit 3
Milford, CT 06460
USA

E info@nt.dental (general information)

E order-us@nt.dental (for orders)

W www.nt.dental

2023_NT_Serienbroschüre_COL_EN_REV001