


# F-SERIES




compatible to:  
Nobel Biocare™ / Nobel Active™



3,5 mm NP




4,3 / 5,0 mm RP



**PRECISION  
MEETS QUALITY.**

**OVERVIEW  
PLATFORMS:**




3,5 mm NP



4,3 / 5,0 mm RP

---









 **MADE IN  
GERMANY**



3,5 mm NP  
Nobel Biocare™ / Nobel Active™

ALL DETAILS  
ONLINE




	<b>Titanium Base (engaged)</b> incl. Abutmentscrew	F 800
	<b>Titanium Base NX (non engaged)</b> incl. Abutmentscrew	GH 0,8 mm F 82350008 GH 1,8 mm F 82350018
	<b>Scanbody 3D-Guide®</b> for intra + extra oral scanning	F 9.S3D3.500
	<b>Lab Analog</b>	F 50
	<b>Implant Pick.Up</b> for open-tray with integrated screw	F TR-NP013.5
	<b>Abutment Screw</b> Recommended tightening torques 25 Ncm	F 60
	<b>DIM-ANALOG®</b>	F 5.DIM.350
	<b>Sleeve (for DIM-ANALOG®)</b> for plaster model	DIM 5.XSL.A10

**F-SERIES**

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**










**WWW.NT.DENTAL**



3,5 mm NP  
Nobel Biocare™ / Nobel Active™

ALL DETAILS  
ONLINE




	<b>nt-Preform®</b> incl. final Abutment Screw	GH 0,5 mm F 10.Mi3.500.05
	<b>Scanbody 3D-Guide®</b> for intra + extra oral scanning	F 9.S3D3.500
	<b>nt-Mill Lab Analog</b> for nt-Preform®	F 5.La3.500
	<b>colored milling screw</b> for nt-Preform®	F 6.Sc3.500
	<b>Abutment Screw</b> Recommended tightening torques 25 Ncm	F 60
	<b>DIM-ANALOG®</b>	F 5.DIM.350
	<b>Sleeve (for DIM-ANALOG®)</b> for plaster model	DIM 5.XSLA10
	<b>Abutment Grip</b> for nt-Preform® Holder	W12.100.00
	<b>Customized Abutment</b> milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

**F-SERIES**

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



3,5 mm NP  
Nobel Biocare™ / Nobel Active™

ALL DETAILS  
ONLINE



**2-CONnect®-Abutment Set (straight)**

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm F 800 S

GH 2,5 mm F-2CON 8035002S



**2-CONnect®-Abutment**

Male Part, recommended tightening torques 25 Ncm

GH 0,5 mm F 800 M

GH 2,5 mm F-2CON 8035002M



**2-CONnect® Scanbody 3D-Guide®**

for intra + extra oral scanning

2CON 9.S3D1.800



**LOcON COPING**

For more information, see our LOcON brochure

LOC 852001



**2-CONnect®-Coping**

Female Part

2CON F-U



**2-CONnect®-Coping-Screw**

Recommended tightening torques 20 Ncm

KS 61



**2-CONnect®-DIM-ANALOG®**

for straight Male Part, KS 61

2CON 5.DIM.180



**Sleeve (for DIM-ANALOG®)**

for plaster model

DIM 5.XSL.A30



**2-CONnect®-Implant pick-up**

for open-tray with integrated screw

2CON TR-00021.8



**2-CONnect®-lab Analog**

for straight Male Parts


2CON 5.La1.800

**F-SERIES**

Notice: Products indicated with © and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



3,5 mm NP  
Nobel Biocare™ / Nobel Active™

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm F-2CON 80350A3S



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm F-2CON 80350A3M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



### Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20  
Ncm

F 6.Ms3.500



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

F-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

MADE IN  
GERMANY

WWW.NT.DENTAL



4,3/5,0 mm RP  
Nobel Biocare™ / Nobel Active™



**Titanium Base (engaged)**  
incl. Abutmentscrew

F 810



**Titanium Base NX (non engaged)**  
incl. Abutmentscrew

GH 0,8 mm F 82435008

GH 1,8 mm F 82435018



**Scanbody 3D-Guide®**  
for intra + extra oral scanning

F 9.S3D4.350



**Lab Analog**

F 51



**Implant Pick.Up**  
for open-tray with integrated screw

F TR-RP024.3



**Abutment Screw**  
Recommended tightening torques 35 Ncm

F 61




**DIM-ANALOG®**

F 5.DIM.435












**Sleeve (for DIM-ANALOG®)**  
for plaster model


DIM 5.XSL.A20



4,3/5,0 mm RP  
Nobel Biocare™ / Nobel Active™

	<b>nt-Preform®</b> incl. final Abutment Screw	GH 0,5 mm F 10.Mi4.350.05
	<b>Scanbody 3D-Guide®</b> for intra + extra oral scanning	F 9.S3D4.350
	<b>nt-Mill Lab Analog</b> for nt-Preform®	F 5.La4.350
	<b>colored milling screw</b> for nt-Preform®	F 6.Sc4.350
	<b>Abutment Screw</b> Recommended tightening torques: 35 Ncm	F 61
	<b>DIM-ANALOG®</b>	F 5.DIM.435
	<b>Sleeve (for DIM-ANALOG®)</b> for plaster model	DIM 5.XSLA20
	<b>Abutment Grip</b> for nt-Preform® Holder	W12.200.00
	<b>Customized Abutment</b> milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1





4,3/5,0 mm RP  
Nobel Biocare™ / Nobel Active™

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm F 810 S

GH 2,5 mm F-2CON 8043002S



### 2-CONnect®-Abutment

Male Part, recommended tightening torques 35 Ncm

GH 0,5 mm F 810 M

GH 2,5 mm F-2CON 8043002M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 30 Ncm

KS 61



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

# F-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 MADE IN  
GERMANY

[WWW.NT.DENTAL](http://WWW.NT.DENTAL)



4,3/5,0 mm RP

Nobel Biocare™ / Nobel Active™

**2-CONnect®-Abutment Set (angled)**incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm F-2CON 80430A3S

**2-CONnect®-Abutment (angled)**

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm F-2CON 80430A3M

**2-CONnect® Scanbody 3D-Guide®**

for intra + extra oral scanning

2CON 9.S3D1.600

**2-CONnect®-Coping**

Female Part

2CON F-U

**2-CONnect®-Coping-Screw**

Recommended tightening torques 15 Ncm

KS 62

**Male Part Screw**for angled 2-CONnect® Abutment, recommended tightening torques 20  
Ncm

F 6.Ms4.350

**2-CONnect®-DIM-ANALOG®**

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160

**Sleeve (for DIM-ANALOG®)**

for plaster model

DIM 5.XSL.A30

**2-CONnect®-Implant pick-up**

for open-tray with integrated screw

2CON TR-00031.6

**2-CONnect®-lab Analog**

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600



## EUROPE

### **NT-Trading GmbH & Co. KG**

G.-Braun-Str. 18  
76187 Karlsruhe  
Germany

**E** [info@nt.dental](mailto:info@nt.dental)

**W** [www.nt.dental](http://www.nt.dental)

## NORTH AMERICA

### **NT-US, Inc.**

132 Pepes Farm Road,  
Unit 3  
Milford, CT 06460  
USA

**E** [info@nt.dental](mailto:info@nt.dental) (general information)

**E** [order-us@nt.dental](mailto:order-us@nt.dental) (for orders)

**W** [www.nt.dental](http://www.nt.dental)

2023\_NT\_Serienbroschüre\_F\_EN\_REV002