


H-SERIES



compatible to:
Biomet 3i® Certain®




3,4 mm NP



5,0 mm WP



4,1 mm RP













3,4 mm NP
Biomet 3i® Certain®

ALL DETAILS
ONLINE












	Titanium Base (engaged) incl. Abutmentscrew	H 800
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 0,3 mm H 82340003 GH 1,3 mm H 82340013
	Scanbody 3D-Guide® for intra + extra oral scanning	H 9.S3D3.400
	Lab Analog	H 50
	Implant Pick.Up for open-tray with integrated screw	H TR-NP013.4
	Abutment Screw Recommended tightening torques 20 Ncm	H 60
	DIM-ANALOG®	H 5.DIM.340
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSL.A10



3,4 mm NP
Biomet 3i® Certain®

ALL DETAILS
ONLINE



	nt-Preform® incl. final Abutment Screw	H 10.Mi3.400.05
	Scanbody 3D-Guide® for intra + extra oral scanning	H 9.S3D3.400
	nt-Mill Lab Analog for nt-Preform®	H 5.La3.400
	colored milling screw for nt-Preform®	H 6.Sc3.450
	Abutment Screw Recommended tightening torques 20 Ncm	H 60
	DIM-ANALOG®	H 5.DIM.340
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSLA10
	Abutment Grip for nt-Preform® Holder	W12.100.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

H-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,4 mm NP
Biomet 3i® Certain®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm H 800 S

GH 2,5 mm H-2CON 8034002S



2-CONnect®-Abutment

Male Part, recommended tightening torques 20 Ncm

GH 0,5 mm H 800 M

GH 2,5 mm H-2CON 8034002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.400



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 60



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 60

2CON 5.DIM.140



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00011.4



2-CONnect®-lab Analog

for straight Male Parts

2CON 5.La1.400

H-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,4 mm NP
Biomet 3i® Certain®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm H-2CON 80340A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm H-2CON 80340A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

H 6.Ms3.450



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

H-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY









WWW.NT.DENTAL



4,1 mm RP
Biomet 3i® Certain®

ALL DETAILS
ONLINE












	Titanium Base (engaged) incl. Abutmentscrew	H 810
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 0,3 mm H 82410003 GH 1,3 mm H 82410013
	Scanbody 3D-Guide® for intra + extra oral scanning	H 9.S3D4.150
	Lab Analog	H 51
	Implant Pick.Up for open-tray with integrated screw	H TR-RP024.1
	Abutment Screw Recommended tightening torques 20 Ncm	H 60
	DIM-ANALOG®	H 5.DIM.410
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSL.A20



4,1 mm RP
Biomet 3i® Certain®

ALL DETAILS
ONLINE



	nt-Preform® incl. final Abutment Screw	GH 0,5 mm H 10.Mi4.100.05
	Scanbody 3D-Guide® for intra + extra oral scanning	H 9.S3D4.150
	nt-Mill Lab Analog for nt-Preform®	H 5.La4.100
	colored milling screw for nt-Preform®	H 6.Sc3.450
	Abutment Screw Recommended tightening torques 20 Ncm	H 60
	DIM-ANALOG®	H 5.DIM.410
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSLA20
	Abutment Grip for nt-Preform® Holder	W12.200.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1



4,1 mm RP
Biomet 3i® Certain®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm H 810 S

GH 2,5 mm H-2CON 8041002S



2-CONnect®-Abutment

Male Part, recommended tightening torques 20 Ncm

GH 0,5 mm H 810 M

GH 2,5 mm H-2CON 8041002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.400



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 60



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 60

2CON 5.DIM.140



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00011.4



2-CONnect®-lab Analog

for straight Male Parts

2CON 5.La1.400



4,1 mm RP
Biomet 3i® Certain®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm H-2CON 80410A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm H-2CON 80410A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

H 6.Ms3.450



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

H-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY









WWW.NT.DENTAL



5,0 mm WP
Biomet 3i® Certain®

ALL DETAILS
ONLINE












	Titanium Base (engaged) incl. Abutmentscrew	H 820
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 0,3 mm H 82500003 GH 1,3 mm H 82500013
	Scanbody 3D-Guide® for intra + extra oral scanning	H 9.S3D4.150
	Lab Analog	H 52
	Implant Pick.Up for open-tray with integrated screw	H TR-WP035.0
	Abutment Screw Recommended tightening torques 20 Ncm	H 60
	DIM-ANALOG®	H 5.DIM.500
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSL.A30



5,0 mm WP
Biomet 3i® Certain®

ALL DETAILS
ONLINE



	nt-Preform® incl. final Abutment Screw	GH 0,5 mm H 10.Mi5.000.05
	Scanbody 3D-Guide® for intra + extra oral scanning	H 9.S3D4.150
	nt-Mill Lab Analog for nt-Preform®	H 5.La5.000
	colored milling screw for nt-Preform®	H 6.Sc3.450
	Abutment Screw Recommended tightening torques: 20 Ncm	H 60
	DIM-ANALOG®	H 5.DIM.500
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSLA30
	Abutment Grip for nt-Preform® Holder	W12.300.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

EUROPE

NT-Trading GmbH & Co. KG

G.-Braun-Str. 18
76187 Karlsruhe
Germany

E info@nt.dental

W www.nt.dental

NORTH AMERICA

NT-US, Inc.

132 Pepes Farm Road,
Unit 3
Milford, CT 06460
USA

E info@nt.dental (general information)

E order-us@nt.dental (for orders)

W www.nt.dental

2023_NT_Serienbroschüre_H_EN_REV002