

HIO -SERIES





compatible to:
HIOSSEN® / OSSTEM TS®



3,5 mm




4,0 - 6,0 mm



**PRECISION
MEETS QUALITY.**

**OVERVIEW
PLATFORMS:**




3,5 mm



4,0 - 6,0 mm









 **MADE IN
GERMANY**



3,5 mm
HIOSEN® OSSTEM TS®

ALL DETAILS
ONLINE




	Titanium Base (engaged) incl. Abutmentscrew	HIO 80350000
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 1,0 mm HIO 82350010 GH 2,0 mm HIO 82350020
	Scanbody 3D-Guide® for intra + extra oral scanning	HIO 9.S3D3.500
	Lab Analog	HIO 5.La3.500
	Implant Pick.Up for open-tray with integrated screw	HIO TR-00013.5
	Abutment Screw Recommended tightening torques 20 Ncm	HIO 603500
	DIM-ANALOG®	HIO 5.DIM.350
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSL.A10

HIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**










WWW.NT.DENTAL



3,5 mm
HIOSSSEN® OSSTEM TS®

ALL DETAILS
ONLINE




	nt-Preform® incl. final Abutment Screw	GH 0,5 mm HIO 10.Mi3.500.05
	Scanbody 3D-Guide® for intra + extra oral scanning	HIO 9.S3D3.500
	nt-Mill Lab Analog for nt-Preform®	HIO 5.La3.500
	colored milling screw for nt-Preform®	HIO 6.Sc3.500
	Abutment Screw Recommended tightening torques 20 Ncm	HIO 603500
	DIM-ANALOG®	HIO 5.DIM.350
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSLA10
	Abutment Grip for nt-Preform® Holder	W 12.100.000
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

HIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,5 mm
HIOSEN® OSSTEM TS®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm HIO-2CON 8035001S

GH 2,5 mm HIO-2CON 8035002S



2-CONnect®-Abutment

Male Part, recommended tightening torques 20 Ncm

GH 0,5 mm HIO-2CON 8035001M

GH 2,5 mm HIO-2CON 8035002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

HIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY

WWW.NT.DENTAL



3,5 mm
HIOSEN® OSSTEM TS®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm HIO-2CON 80350A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm HIO-2CON 80350A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

HIO 6.Ms3.500



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM.160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

HIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY









WWW.NT.DENTAL



4,0 - 6,0 mm
HIOSEN® OSSTEM TS®

ALL DETAILS
ONLINE




	Titanium Base (engaged) incl. Abutmentscrew	HIO 80406000
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 1,0 mm HIO 82406010 GH 2,0 mm HIO 82406020
	Scanbody 3D-Guide® for intra + extra oral scanning	HIO 9.S3D4.060
	Lab Analog	HIO 5.La4.060
	Implant Pick.Up for open-tray with integrated screw	HIO TR-00024.0
	Abutment Screw Recommended tightening torques 30 Ncm	HIO 604060
	DIM-ANALOG®	HIO 5.DIM.406
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSL.A20

HIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY










WWW.NT.DENTAL



4,0 - 6,0 mm
HIOSEN® OSSTEM TS®

ALL DETAILS
ONLINE




	nt-Preform® incl. final Abutment Screw	GH 0,5 mm HIO 10.Mi4.060.05
	Scanbody 3D-Guide® for intra + extra oral scanning	HIO 9.S3D4.060
	nt-Mill Lab Analog for nt-Preform®	HIO 5.La4.060
	colored milling screw for nt-Preform®	HIO 6.Sc4.060
	Abutment Screw Recommended tightening torques 30 Ncm	HIO 604060
	DIM-ANALOG®	HIO 5.DIM.406
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSLA20
	Abutment Grip for nt-Preform® Holder	W 12.200.000
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

HIO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,0 - 6,0 mm
HIOSEN® OSSTEM TS®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm HIO-2CON 8040601S

GH 2,5 mm HIO-2CON 8040602S



2-CONnect®-Abutment

Male Part, recommended tightening torques 30 Ncm

GH 0,5 mm HIO-2CON 8040601M

GH 2,5 mm HIO-2CON 8040602M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw


2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts

2CON 5.La1.800



4,0 - 6,0 mm
HIOSEN® OSSTEM TS®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm HIO-2CON 80406A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm HIO-2CON 80406A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

HIO 6.Ms4060



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM.160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

EUROPE

NT-Trading GmbH & Co. KG

G.-Braun-Str. 18
76187 Karlsruhe
Germany

E info@nt.dental

W www.nt.dental

NORTH AMERICA

NT-US, Inc.

132 Pepes Farm Road,
Unit 3
Milford, CT 06460
USA

E info@nt.dental (general information)

E order-us@nt.dental (for orders)

W www.nt.dental

2023_NT_Serienbroschüre_HIO_EN_REV001