

L-SERIES




compatible to:
Straumann Bone Level®



3,3 mm NC




4,1 / 4,8 mm RC



**PRECISION
MEETS QUALITY.**

OVERVIEW PLATFORMS:



A: over surface

B: over edge




3,3 mm NC



4,1 / 4,8 mm RC

 **MADE IN
GERMANY**



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

L 800



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,7 mm L 82330007

GH 1,7 mm L 82330017



Scanbody 3D-Guide®

for intra + extra oral scanning

L 9.S3D3.300



Lab Analog

L 50



Implant Pick.Up

for open-tray with integrated screw

L TR-NC013.3



Abutment Screw

Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.330



Sleeve (for DIM-ANALOG®)

for plaster model


DIM 5.XSL.A10

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**










WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE




	nt-Preform® incl. final Abutment Screw	GH 0,5 mm L 10.Mi3.300.05
	Scanbody 3D-Guide® for intra + extra oral scanning	L 9.S3D3.300
	nt-Mill Lab Analog for nt-Preform®	L 5.La3.300
	colored milling screw for nt-Preform®	L 6.Sc3.341
	Abutment Screw Recommended tightening torques: 35 Ncm	L 60
	DIM-ANALOG®	L 5.DIM.330
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSL.A10
	Abutment Grip for nt-Preform® Holder	W12.100.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm L 800 S

GH 2,5 mm L-2CON 8033002S



2-CONnect®-Abutment

Male Part

GH 0,5 mm L 800 M

GH 2,5 mm L-2CON 8033002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 30Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 80330A3S
B: over edge	GH 2,5 mm	L-2CON 80330C3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 80330A3M
B: over edge	GH 2,5 mm	L-2CON 80330C3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

L 6.Ms3.341



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



Titanium Base (engaged)
incl. Abutmentscrew

L 810



Titanium Base NX (non engaged)
incl. Abutmentscrew

GH 0,7 mm L 82414807

GH 1,7 mm L 82414817



Scanbody 3D-Guide®
for intra + extra oral scanning

L 9.S3D4.148



Lab Analog

L 51



Implant Pick.Up
for open-tray with integrated screw

L TR-NC024.1



Abutment Screw
Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.414



Sleeve (for DIM-ANALOG®)
for plaster model


DIM 5.XSL.A20

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



nt-Preform®
incl. final Abutment Screw

GH 0,5 mm L 10.Mi4.148.05



Scanbody 3D-Guide®
for intra + extra oral scanning

L 9.S3D4.148



nt-Mill Lab Analog
for nt-Preform®

L 5.La4.148



colored milling screw
for nt-Preform®

L 6.Sc3.341



Abutment Screw
Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.414



Sleeve (for DIM-ANALOG®)
for plaster model

DIM 5.XSL.A20



Abutment Grip
for nt-Preform® Holder

W12.200.00



Customized Abutment
milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONNECT®-Abutment Set (straight)

incl. 2-CONNECT®-Coping Female Part and 2-CONNECT®-Coping-screw

GH 0,5 mm L 810 S

GH 2,5 mm L-2CON 8041482S



2-CONNECT®-Abutment

Male Part

GH 0,5 mm L 810 M

GH 2,5 mm L-2CON 8041482M



2-CONNECT® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONNECT®-Coping

Female Part

2CON F-U



2-CONNECT®-Coping-Screw

Recommended tightening torques 30Ncm

KS 61



2-CONNECT®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONNECT®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONNECT®-lab Analog

for straight Male Parts


2CON 5.La1.800

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 84148A3S
B: over edge	GH 2,5 mm	L-2CON 84148C3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 84148A3M
B: over edge	GH 2,5 mm	L-2CON 84148C3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

L 6.Ms3.341



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

L 800



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,7 mm L 82330007

GH 1,7 mm L 82330017



Scanbody 3D-Guide®

for intra + extra oral scanning

L 9.S3D3.300



Lab Analog

L 50



Implant Pick.Up

for open-tray with integrated screw

L TR-NC013.3



Abutment Screw

Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.330



Sleeve (for DIM-ANALOG®)

for plaster model


DIM 5.XSL.A10

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**










WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE




	nt-Preform® incl. final Abutment Screw	GH 0,5 mm L 10.Mi3.300.05
	Scanbody 3D-Guide® for intra + extra oral scanning	L 9.S3D3.300
	nt-Mill Lab Analog for nt-Preform®	L 5.La3.300
	colored milling screw for nt-Preform®	L 6.Sc3.341
	Abutment Screw Recommended tightening torques: 35 Ncm	L 60
	DIM-ANALOG®	L 5.DIM.330
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSL.A10
	Abutment Grip for nt-Preform® Holder	W12.100.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm L 800 S

GH 2,5 mm L-2CON 8033002S



2-CONnect®-Abutment

Male Part

GH 0,5 mm L 800 M

GH 2,5 mm L-2CON 8033002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 30Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 80330A3S
B: over edge	GH 2,5 mm	L-2CON 80330C3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 80330A3M
B: over edge	GH 2,5 mm	L-2CON 80330C3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

L 6.Ms3.341



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



Titanium Base (engaged)
incl. Abutmentscrew

L 810



Titanium Base NX (non engaged)
incl. Abutmentscrew

GH 0,7 mm L 82414807

GH 1,7 mm L 82414817



Scanbody 3D-Guide®
for intra + extra oral scanning

L 9.S3D4.148



Lab Analog

L 51



Implant Pick.Up
for open-tray with integrated screw

L TR-NC024.1



Abutment Screw
Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.414



Sleeve (for DIM-ANALOG®)
for plaster model


DIM 5.XSL.A20

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



nt-Preform®
incl. final Abutment Screw

GH 0,5 mm L 10.Mi4.148.05



Scanbody 3D-Guide®
for intra + extra oral scanning

L 9.S3D4.148



nt-Mill Lab Analog
for nt-Preform®

L 5.La4.148



colored milling screw
for nt-Preform®

L 6.Sc3.341



Abutment Screw
Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.414



Sleeve (for DIM-ANALOG®)
for plaster model

DIM 5.XSL.A20



Abutment Grip
for nt-Preform® Holder

W12.200.00



Customized Abutment
milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm L 810 S
GH 2,5 mm L-2CON 8041482S



2-CONnect®-Abutment

Male Part

GH 0,5 mm L 810 M
GH 2,5 mm L-2CON 8041482M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 30Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 84148A3S
B: over edge	GH 2,5 mm	L-2CON 84148C3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 84148A3M
B: over edge	GH 2,5 mm	L-2CON 84148C3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

L 6.Ms3.341



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY

WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

L 800



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,7 mm L 82330007

GH 1,7 mm L 82330017



Scanbody 3D-Guide®

for intra + extra oral scanning

L 9.S3D3.300



Lab Analog

L 50



Implant Pick.Up

for open-tray with integrated screw

L TR-NC013.3



Abutment Screw

Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.330



Sleeve (for DIM-ANALOG®)

for plaster model


DIM 5.XSL.A10

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



nt-Preform®
incl. final Abutment Screw

GH 0,5 mm L 10.Mi3.300.05



Scanbody 3D-Guide®
for intra + extra oral scanning

L 9.S3D3.300



nt-Mill Lab Analog
for nt-Preform®

L 5.La3.300



colored milling screw
for nt-Preform®

L 6.Sc3.341



Abutment Screw
Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.330



Sleeve (for DIM-ANALOG®)
for plaster model

DIM 5.XSL.A10



Abutment Grip
for nt-Preform® Holder

W12.100.00



Customized Abutment
milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm L 800 S

GH 2,5 mm L-2CON 8033002S



2-CONnect®-Abutment

Male Part

GH 0,5 mm L 800 M

GH 2,5 mm L-2CON 8033002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 30Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 80330A3S
B: over edge	GH 2,5 mm	L-2CON 80330C3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 80330A3M
B: over edge	GH 2,5 mm	L-2CON 80330C3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

L 6.Ms3.341



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



Titanium Base (engaged)
incl. Abutmentscrew

L 810



Titanium Base NX (non engaged)
incl. Abutmentscrew

GH 0,7 mm L 82414807

GH 1,7 mm L 82414817



Scanbody 3D-Guide®
for intra + extra oral scanning

L 9.S3D4.148



Lab Analog

L 51



Implant Pick.Up
for open-tray with integrated screw

L TR-NC024.1



Abutment Screw
Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.414



Sleeve (for DIM-ANALOG®)
for plaster model


DIM 5.XSL.A20

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



nt-Preform®
incl. final Abutment Screw

GH 0,5 mm L 10.Mi4.148.05



Scanbody 3D-Guide®
for intra + extra oral scanning

L 9.S3D4.148



nt-Mill Lab Analog
for nt-Preform®

L 5.La4.148



colored milling screw
for nt-Preform®

L 6.Sc3.341



Abutment Screw
Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.414



Sleeve (for DIM-ANALOG®)
for plaster model

DIM 5.XSL.A20



Abutment Grip
for nt-Preform® Holder

W12.200.00



Customized Abutment
milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm L 810 S

GH 2,5 mm L-2CON 8041482S



2-CONnect®-Abutment

Male Part

GH 0,5 mm L 810 M

GH 2,5 mm L-2CON 8041482M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 30Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 84148A3S
B: over edge	GH 2,5 mm	L-2CON 84148C3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 84148A3M
B: over edge	GH 2,5 mm	L-2CON 84148C3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

L 6.Ms3.341



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

L 800



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,7 mm L 82330007

GH 1,7 mm L 82330017



Scanbody 3D-Guide®

for intra + extra oral scanning

L 9.S3D3.300



Lab Analog

L 50



Implant Pick.Up

for open-tray with integrated screw

L TR-NC013.3



Abutment Screw

Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.330



Sleeve (for DIM-ANALOG®)

for plaster model


DIM 5.XSL.A10

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**










WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE




	nt-Preform® incl. final Abutment Screw	GH 0,5 mm L 10.Mi3.300.05
	Scanbody 3D-Guide® for intra + extra oral scanning	L 9.S3D3.300
	nt-Mill Lab Analog for nt-Preform®	L 5.La3.300
	colored milling screw for nt-Preform®	L 6.Sc3.341
	Abutment Screw Recommended tightening torques: 35 Ncm	L 60
	DIM-ANALOG®	L 5.DIM.330
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSL.A10
	Abutment Grip for nt-Preform® Holder	W12.100.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm L 800 S

GH 2,5 mm L-2CON 8033002S



2-CONnect®-Abutment

Male Part

GH 0,5 mm L 800 M

GH 2,5 mm L-2CON 8033002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 30Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY

WWW.NT.DENTAL



3,3 mm NC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 80330A3S
B: over edge	GH 2,5 mm	L-2CON 80330C3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 80330A3M
B: over edge	GH 2,5 mm	L-2CON 80330C3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

L 6.Ms3.341



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



Titanium Base (engaged)
incl. Abutmentscrew

L 810



Titanium Base NX (non engaged)
incl. Abutmentscrew

GH 0,7 mm L 82414807

GH 1,7 mm L 82414817



Scanbody 3D-Guide®
for intra + extra oral scanning

L 9.S3D4.148



Lab Analog

L 51



Implant Pick.Up
for open-tray with integrated screw

L TR-NC024.1



Abutment Screw
Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.414



Sleeve (for DIM-ANALOG®)
for plaster model


DIM 5.XSL.A20

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



nt-Preform®
incl. final Abutment Screw

GH 0,5 mm L 10.Mi4.148.05



Scanbody 3D-Guide®
for intra + extra oral scanning

L 9.S3D4.148



nt-Mill Lab Analog
for nt-Preform®

L 5.La4.148



colored milling screw
for nt-Preform®

L 6.Sc3.341



Abutment Screw
Recommended tightening torques: 35 Ncm

L 60



DIM-ANALOG®

L 5.DIM.414



Sleeve (for DIM-ANALOG®)
for plaster model

DIM 5.XSL.A20



Abutment Grip
for nt-Preform® Holder

W12.200.00



Customized Abutment
milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONNECT®-Abutment Set (straight)

incl. 2-CONNECT®-Coping Female Part and 2-CONNECT®-Coping-screw

GH 0,5 mm L 810 S

GH 2,5 mm L-2CON 8041482S



2-CONNECT®-Abutment

Male Part

GH 0,5 mm L 810 M

GH 2,5 mm L-2CON 8041482M



2-CONNECT® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONNECT®-Coping

Female Part

2CON F-U



2-CONNECT®-Coping-Screw

Recommended tightening torques 30Ncm

KS 61



2-CONNECT®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONNECT®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



2-CONNECT®-lab Analog

for straight Male Parts


2CON 5.La1.800

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,1 / 4,8 mm RC
Straumann Bone Level®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 84148A3S
B: over edge	GH 2,5 mm	L-2CON 84148C3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

A: over surface	GH 2,5 mm	L-2CON 84148A3M
B: over edge	GH 2,5 mm	L-2CON 84148C3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

L 6.Ms3.341



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

L-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL

NOTES

[illegible]

EUROPE

NT-Trading GmbH & Co. KG

G.-Braun-Str. 18
76187 Karlsruhe
Germany

E info@nt.dental

W www.nt.dental

NORTH AMERICA

NT-US, Inc.

132 Pepes Farm Road,
Unit 3
Milford, CT 06460
USA

E info@nt.dental (general information)

E order-us@nt.dental (for orders)

W www.nt.dental

2023_NT_Serienbroschüre_L_EN_REV002