


# MEX -SERIES



compatible to:  
Medentis Medical® ICX®




3,7 - 4,8 mm



**PRECISION  
MEETS QUALITY.**


**OVERVIEW  
PLATFORMS:**



3,7 - 4,8 mm

---









 **MADE IN  
GERMANY**



3,7 - 4,8 mm  
Medentis Medical® ICX®

ALL DETAILS  
ONLINE




	<b>Titanium Base (engaged)</b> incl. Abutmentscrew	MEX 80374800
	<b>Titanium Base NX (non engaged)</b> incl. Abutmentscrew	GH 1,0 mm MEX 82374810 GH 2,0 mm MEX 82374820
	<b>Scanbody 3D-Guide®</b> for intra + extra oral scanning	MEX 9.S3D3.748
	<b>Lab Analog</b>	MEX 5.La3.748
	<b>Implant Pick.Up</b> for open-tray with integrated screw	MEX TR-00013.7
	<b>Abutment Screw</b> Recommended tightening torques 30 Ncm	MEX 603748
	<b>DIM-ANALOG®</b>	MEX 5.DIM.374
	<b>Sleeve (for DIM-ANALOG®)</b> for plaster model	DIM 5.XSL.A20

**MEX-SERIES**

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**










**WWW.NT.DENTAL**



3,7 - 4,8 mm  
Medentis Medical® ICX®

ALL DETAILS  
ONLINE




	<b>nt-Preform®</b> incl. final Abutment Screw	GH 0,5 mm MEX 10.Mi3.748.05
	<b>Scanbody 3D-Guide®</b> for intra + extra oral scanning	MEX 9.S3D3.748
	<b>nt-Mill Lab Analog</b> for nt-Preform®	MEX 5.La3.748
	<b>colored milling screw</b> for nt-Preform®	MEX 6.Sc4.148
	<b>Abutment Screw</b> Recommended tightening torques 30 Ncm	MEX 603748
	<b>DIM-ANALOG®</b>	MEX 5.DIM.374
	<b>Sleeve (for DIM-ANALOG®)</b> for plaster model	DIM 5.XSLA20
	<b>Abutment Grip</b> for nt-Preform® Holder	W 12.200.000
	<b>Customized Abutment</b> milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

**MEX-SERIES**

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



3,7 - 4,8 mm  
Medentis Medical® ICX®

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm MEX-2CON 8041481S

GH 2,5 mm MEX-2CON 8041482S



### 2-CONnect®-Abutment

Male Part, recommended tightening torques 30 Ncm

GH 0,5 mm MEX-2CON 8041481M

GH 2,5 mm MEX-2CON 8041482M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

# MEX-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 MADE IN  
GERMANY

[WWW.NT.DENTAL](http://WWW.NT.DENTAL)



3,7 - 4,8 mm  
Medentis Medical® ICX®

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm MEX-2CON 80414A3S



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm MEX-2CON 80414A3M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



### Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20  
Ncm

MEX 6.Ms4.148



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM.160



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

# MEX-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

MADE IN  
GERMANY

WWW.NT.DENTAL



## EUROPE

### **NT-Trading GmbH & Co. KG**

G.-Braun-Str. 18  
76187 Karlsruhe  
Germany

**E** [info@nt.dental](mailto:info@nt.dental)

**W** [www.nt.dental](http://www.nt.dental)

## NORTH AMERICA

### **NT-US, Inc.**

132 Pepes Farm Road,  
Unit 3  
Milford, CT 06460  
USA

**E** [info@nt.dental](mailto:info@nt.dental) (general information)

**E** [order-us@nt.dental](mailto:order-us@nt.dental) (for orders)

**W** [www.nt.dental](http://www.nt.dental)

2023\_NT\_Serienbroschüre\_MEX\_EN\_REV001