


NEO -SERIES




compatible to:
NEOSS Implant System® / NEOSS™




3,25 mm



3,5 - 5,5 mm













3,25 mm
NEOSS Implant System® NEOSS®

ALL DETAILS
ONLINE




	Titanium Base (engaged) incl. Abutmentscrew	NEO 80325000
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 0,5 mm NEO 82320005 GH 1,5 mm NEO 82320015
	Scanbody 3D-Guide® for intra + extra oral scanning	NEO 9.S3D3.250
	Lab Analog	NEO 5.La3.250
	Implant Pick.Up for open-tray with integrated screw	NEO TR-000232
	Abutment Screw Recommended tightening torques 30 Ncm	NEO 603200
	DIM-ANALOG®	NEO 5.DIM.325
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSL.A10

NEO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY










WWW.NT.DENTAL



3,25 mm
NEOSS Implant System® NEOSS®

ALL DETAILS
ONLINE




	nt-Preform® incl. final Abutment Screw	GH 0,5 mm	NEO 10.Mi3.250.05
	Scanbody 3D-Guide® for intra + extra oral scanning		NEO 9.S3D3.250
	nt-Mill Lab Analog for nt-Preform®		NEO 5.La3.250
	colored milling screw for nt-Preform®		NEO 6.Sc3.200
	Abutment Screw Recommended tightening torques 30 Ncm		NEO 603200
	DIM-ANALOG®		NEO 5.DIM.325
	Sleeve (for DIM-ANALOG®) for plaster model		DIM 5.XSLA20
	Abutment Grip for nt-Preform® Holder		W 12.100.000
	Customized Abutment milled by NT-iDent / Not available in the USA		TZ15.TiO.Si1

NEO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY

WWW.NT.DENTAL



3,5 - 5,5 mm
NEOSS Implant System® NEOSS®

ALL DETAILS
ONLINE



Titanium Base (engaged)
incl. Abutmentscrew

NEO 80355500



Titanium Base NX (non engaged)
incl. Abutmentscrew

GH 0,5 mm NEO 82355505
GH 1,5 mm NEO 82355515



Scanbody 3D-Guide®
for intra + extra oral scanning

NEO 9.S3D3.555



Lab Analog

NEO 5.La3.555



Implant Pick.Up
for open-tray with integrated screw

NEO TR-00013.5



Abutment Screw
Recommended tightening torques 30 Ncm

NEO 603555



DIM-ANALOG®

NEO 5.DIM.355



Sleeve (for DIM-ANALOG®)
for plaster model


DIM 5.XSL.A20

NEO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**










WWW.NT.DENTAL



3,5 - 5,5 mm
NEOSS Implant System® NEOSS®

ALL DETAILS
ONLINE




	nt-Preform® incl. final Abutment Screw	GH 0,5 mm	NEO 10.Mi3.555.05
	Scanbody 3D-Guide® for intra + extra oral scanning		NEO 9.S3D3.555
	nt-Mill Lab Analog for nt-Preform®		NEO 5.La3.555
	colored milling screw for nt-Preform®		NEO 6.Sc3.555
	Abutment Screw Recommended tightening torques 30 Ncm		NEO 603555
	DIM-ANALOG®		NEO 5.DIM.355
	Sleeve (for DIM-ANALOG®) for plaster model		DIM 5.XSLA20
	Abutment Grip for nt-Preform® Holder		W 12.200.000
	Customized Abutment milled by NT-iDent / Not available in the USA		TZ15.TiO.Si1

NEO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



3,5 - 5,5 mm
NEOSS Implant System® NEOSS®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm NEO-2CON 8035551S

GH 2,5 mm NEO-2CON 8035552S



2-CONnect®-Abutment

Male Part, recommended tightening torques 30 Ncm

GH 0,5 mm NEO-2CON 8035551M

GH 2,5 mm NEO-2CON 8035552M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 25 Ncm

KS 61



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw


2CON TR-00021.8



2-CONnect®-lab Analog

for straight Male Parts

2CON 5.La1.800



3,5 - 5,5 mm
NEOSS Implant System® NEOSS®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm NEO-2CON 80355A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm NEO-2CON 80355A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

NEO 6.Ms3.555



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

NEO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 MADE IN
GERMANY

WWW.NT.DENTAL

EUROPE

NT-Trading GmbH & Co. KG

G.-Braun-Str. 18
76187 Karlsruhe
Germany

E info@nt.dental

W www.nt.dental

NORTH AMERICA

NT-US, Inc.

132 Pepes Farm Road,
Unit 3
Milford, CT 06460
USA

E info@nt.dental (general information)

E order-us@nt.dental (for orders)

W www.nt.dental

2023_NT_Serienbroschüre_NEO_EN_REV002