


R-SERIES



compatible to:
Zimmer Dental Tapered Screw-Vent®




3,5 mm




5,7 mm



4,5 mm














3,5 mm
Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE












	Titanium Base (engaged) incl. Abutmentscrew	R 800
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 0,5 mm R 82350005 GH 1,5 mm R 82350015
	Scanbody 3D-Guide® for intra + extra oral scanning	R 9.S3D3.500
	Lab Analog	R 50
	Implant Pick.Up for open-tray with integrated screw	R TR-00013.5
	Abutment Screw Recommended tightening torques 30 Ncm	R 60
	DIM-ANALOG®	R 5.DIM.350
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSL.A10



3,5 mm
Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE




	nt-Preform® incl. final Abutment Screw	R 10.Mi3.500.05
	Scanbody 3D-Guide® for intra + extra oral scanning	R 9.S3D3.500
	nt-Mill Lab Analog for nt-Preform®	R 5.La3.500
	colored milling screw for nt-Preform®	R 6.Sc3.557
	Abutment Screw Recommended tightening torques 30 Ncm	R 60
	DIM-ANALOG®	R 5.DIM.350
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSLA10
	Abutment Grip for nt-Preform® Holder	W12.100.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



3,5 mm
Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm R 800 S

GH 2,5 mm R-2CON 8035002S



2-CONnect®-Abutment

Male Part, recommended tightening torques 30 Ncm

GH 0,5 mm R 800 M

GH 2,5 mm R-2CON 8035002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.400



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 60



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 60

2CON 5.DIM.140



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw


2CON TR-00011.4



2-CONnect®-lab Analog

for straight Male Parts

2CON 5.La1.400



3,5 mm
Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm R-2CON 80350A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm R-2CON 80350A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

R 6.Ms3.557



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY









WWW.NT.DENTAL



4,5 mm
Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE




	Titanium Base (engaged) incl. Abutmentscrew	R 810
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 0,5 mm R 82450005 GH 1,5 mm R 82450015
	Scanbody 3D-Guide® for intra + extra oral scanning	R 9.S3D4.500
	Lab Analog	R 51
	Implant Pick.Up for open-tray with integrated screw	R TR-00024.5
	Abutment Screw Recommended tightening torques 30 Ncm	R 60
	DIM-ANALOG®	R 5.DIM.450
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSL.A20

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY










WWW.NT.DENTAL




4,5 mm
Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



	nt-Preform® incl. final Abutment Screw	GH 0,5 mm R 10.Mi4.500.05
	Scanbody 3D-Guide® for intra + extra oral scanning	R 9.S3D4.500
	nt-Mill Lab Analog for nt-Preform®	R 5.La4.500
	colored milling screw for nt-Preform®	R 6.Sc3.557
	Abutment Screw Recommended tightening torques 30 Ncm	R 60
	DIM-ANALOG®	R 5.DIM.450
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSLA20
	Abutment Grip for nt-Preform® Holder	W12.200.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1



4,5 mm
Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm R 810 S

GH 2,5 mm R-2CON 8045002S



2-CONnect®-Abutment

Male Part, recommended tightening torques 30 Ncm

GH 0,5 mm R 810 M

GH 2,5 mm R-2CON 8045002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.400



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 30 Ncm

KS 60



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 60

2CON 5.DIM.140



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00011.4



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.400

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



4,5 mm
Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm R-2CON 80450A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm R-2CON 80450A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

R 6.Ms3.557



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY








WWW.NT.DENTAL




5,7 mm
Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



	Titanium Base (engaged) incl. Abutmentscrew	R 820
	Titanium Base NX (non engaged) incl. Abutmentscrew	GH 0,5 mm R 82570005 GH 1,5 mm R 82570015
	Scanbody 3D-Guide® for intra + extra oral scanning	R 9.S3D5.700
	Lab Analog	R 52
	Abutment Screw Recommended tightening torques 30 Ncm	R 60
	DIM-ANALOG®	R 5.DIM.570
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSLA30












5,7 mm

Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



	nt-Preform® incl. final Abutment Screw	GH 0,5 mm R 10.Mi5.700.05
	Scanbody 3D-Guide® for intra + extra oral scanning	R 9.S3D5.700
	nt-Mill Lab Analog for nt-Preform®	R 5.La5.700
	colored milling screw for nt-Preform®	R 6.Sc3.557
	Abutment Screw Recommended tightening torques: 30 Ncm	R 60
	DIM-ANALOG®	R 5.DIM.570
	Sleeve (for DIM-ANALOG®) for plaster model	DIM 5.XSLA30
	Abutment Grip for nt-Preform® Holder	W12.300.00
	Customized Abutment milled by NT-iDent / Not available in the USA	TZ15.TiO.Si1

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL

EUROPE

NT-Trading GmbH & Co. KG

G.-Braun-Str. 18
76187 Karlsruhe
Germany

E info@nt.dental

W www.nt.dental

NORTH AMERICA

NT-US, Inc.

132 Pepes Farm Road,
Unit 3
Milford, CT 06460
USA

E info@nt.dental (general information)

E order-us@nt.dental (for orders)

W www.nt.dental

2023_NT_Serienbroschüre_R_EN_REV002