


R-SERIES



compatible to:
Zimmer Dental Tapered Screw-Vent®




3,5 mm




5,7 mm



4,5 mm





3,5 mm

Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

R 800



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,5 mm R 82350005

GH 1,5 mm R 82350015



Scanbody 3D-Guide®

for intra + extra oral scanning

R 9.S3D3.500



Lab Analog

R 50



Implant Pick.Up

for open-tray with integrated screw

R TR-00013.5



Abutment Screw

Recommended tightening torques 30 Ncm

R 60



DIM-ANALOG®

R 5.DIM.350



Sleeve (for DIM-ANALOG®)

for plaster model


DIM 5.XSL.A10

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL












3,5 mm

Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE




| | | |
|---|---|-----------------|
|  | nt-Preform® incl. final Abutment Screw | R 10.Mi3.500.05 |
|  | Scanbody 3D-Guide® for intra + extra oral scanning | R 9.S3D3.500 |
|  | nt-Mill Lab Analog for nt-Preform® | R 5.La3.500 |
|  | colored milling screw for nt-Preform® | R 6.Sc3.557 |
|  | Abutment Screw Recommended tightening torques 30 Ncm | R 60 |
|  | DIM-ANALOG® | R 5.DIM.350 |
|  | Sleeve (for DIM-ANALOG®) for plaster model | DIM 5.XSL.A10 |
|  | Abutment Grip for nt-Preform® Holder | W12.100.00 |
|  | Customized Abutment milled by NT-iDent / Not available in the USA | TZ15.TiO.Si1 |

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



3,5 mm

Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm R 800 S

GH 2,5 mm R-2CON 8035002S



2-CONnect®-Abutment

Male Part, recommended tightening torques 30 Ncm

GH 0,5 mm R 800 M

GH 2,5 mm R-2CON 8035002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.400



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 60



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 60

2CON 5.DIM.140



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00011.4



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.400

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



3,5 mm

Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm R-2CON 80350A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm R-2CON 80350A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

R 6.Ms3.557



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



4,5 mm

Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

R 810



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,5 mm R 82450005

GH 1,5 mm R 82450015



Scanbody 3D-Guide®

for intra + extra oral scanning

R 9.S3D4.500



Lab Analog

R 51



Implant Pick.Up

for open-tray with integrated screw

R TR-00024.5



Abutment Screw

Recommended tightening torques 30 Ncm

R 60



DIM-ANALOG®

R 5.DIM.450



Sleeve (for DIM-ANALOG®)

for plaster model


DIM 5.XSL.A20

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,5 mm

Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



nt-Preform®

incl. final Abutment Screw

GH 0,5 mm R 10.Mi4.500.05



Scanbody 3D-Guide®

for intra + extra oral scanning

R 9.S3D4.500



nt-Mill Lab Analog

for nt-Preform®

R 5.La4.500



colored milling screw

for nt-Preform®

R 6.Sc3.557



Abutment Screw

Recommended tightening torques 30 Ncm

R 60



DIM-ANALOG®

R 5.DIM.450



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A20



Abutment Grip

for nt-Preform® Holder

W12.200.00



Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,5 mm

Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm R 810 S

GH 2,5 mm R-2CON 8045002S



2-CONnect®-Abutment

Male Part, recommended tightening torques 30 Ncm

GH 0,5 mm R 810 M

GH 2,5 mm R-2CON 8045002M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.400



LOcON COPING

For more information, see our LOcON brochure

LOC 852001



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 30 Ncm

KS 60



2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 60

2CON 5.DIM.140



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00011.4



2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.400

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



4,5 mm

Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm R-2CON 80450A3S



2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm R-2CON 80450A3M



2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



2-CONnect®-Coping

Female Part

2CON F-U



2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20
Ncm

R 6.Ms3.557



2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

MADE IN
GERMANY

WWW.NT.DENTAL



5,7 mm

Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



Titanium Base (engaged)

incl. Abutmentscrew

R 820



Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,5 mm R 82570005

GH 1,5 mm R 82570015



Scanbody 3D-Guide®

for intra + extra oral scanning

R 9.S3D5.700



Lab Analog

R 52



Abutment Screw

Recommended tightening torques 30 Ncm

R 60



DIM-ANALOG®

R 5.DIM.570



Sleeve (for DIM-ANALOG®)

for plaster model


DIM 5.XSL.A30

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL



5,7 mm

Zimmer Dental Tapered Screw-Vent®

ALL DETAILS
ONLINE



nt-Preform®

incl. final Abutment Screw

GH 0,5 mm R 10.Mi5.700.05



Scanbody 3D-Guide®

for intra + extra oral scanning

R 9.S3D5.700



nt-Mill Lab Analog

for nt-Preform®

R 5.La5.700



colored milling screw

for nt-Preform®

R 6.Sc3.557



Abutment Screw

Recommended tightening torques: 30 Ncm

R 60



DIM-ANALOG®

R 5.DIM.570



Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSLA30



Abutment Grip

for nt-Preform® Holder

W12.300.00



Customized Abutment

milled by NT-iDent / Not available in the USA

TZ15.TiO.Si1

R-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.
For products available in the USA, please check our website: www.nt.dental/en/

 **MADE IN
GERMANY**

WWW.NT.DENTAL

EUROPE

NT-Trading GmbH & Co. KG

G.-Braun-Str. 18
76187 Karlsruhe
Germany

E info@nt.dental

W www.nt.dental

NORTH AMERICA

NT-US, Inc.

132 Pepes Farm Road,
Unit 3
Milford, CT 06460
USA

E info@nt.dental (general information)

E order-us@nt.dental (for orders)

W www.nt.dental

2023_NT_Serienbroschüre_R_EN_REV002