

# S-SERIES




compatible to:  
Dentsply Implants® / Astra Tech  
Implant System OsseoSpeed® TX



3,5 / 4,0 mm




4,5 / 5,0 mm



**PRECISION  
MEETS QUALITY.**

**OVERVIEW  
PLATFORMS:**




3,5 / 4,0 mm



4,5 / 5,0 mm

 **MADE IN  
GERMANY**



3,5 / 4,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS  
ONLINE



### **Titanium Base (engaged)**

incl. Abutmentscrew

S 800



### **Titanium Base NX (non engaged)**

incl. Abutmentscrew

GH 0,8 mm S 82354008

GH 1,8 mm S 82354018



### **Scanbody 3D-Guide®**

for intra + extra oral scanning

S 9.S3D3.540



### **Lab Analog**

S 50



### **Implant Pick.Up**

for open-tray with integrated screw

S TR-00013.5



### **Abutment Screw**

Recommended tightening torques 25 Ncm

S 60



### **DIM-ANALOG®**

S 5.DIM.354



### **Sleeve (for DIM-ANALOG®)**

for plaster model


DIM 5.XSL.A10

# S-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



3,5 / 4,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS  
ONLINE



### nt-Preform®

incl. final Abutment Screw

GH 0,5 mm S 10.Mi3.540.05



### Scanbody 3D-Guide®

for intra + extra oral scanning

S 9.S3D3.540



### nt-Mill Lab Analog

for nt-Preform®

S 5.La3.540



### colored milling screw

for nt-Preform®

S 6.Sc3.540



### Abutment Screw

Recommended tightening torques 25 Ncm

S 60



### DIM-ANALOG®

S 5.DIM.354



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSLA10



### Abutment Grip

for nt-Preform® Holder

W12.100.00



### Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

# S-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



3,5 / 4,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm S 800 S

GH 2,5 mm S-2CON 8035402S



### 2-CONnect®-Abutment

Male Part, recommended tightening torques 25 Ncm

GH 0,5 mm S 800 M

GH 2,5 mm S-2CON 8035402M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 61



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

# S-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

MADE IN  
GERMANY

WWW.NT.DENTAL



3,5 / 4,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm S-2CON 80354A3S



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm S-2CON 80354A3M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



### Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20  
Ncm

S 6.Ms3.540



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

# S-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

MADE IN  
GERMANY

WWW.NT.DENTAL



4,5 / 5,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS  
ONLINE



### Titanium Base (engaged)

incl. Abutmentscrew

S 820



### Titanium Base NX (non engaged)

incl. Abutmentscrew

GH 0,8 mm S 82455008

GH 1,8 mm S 82455018



### Scanbody 3D-Guide®

for intra + extra oral scanning

S 9.S3D4.550



### Lab Analog

S 52



### Implant Pick.Up

for open-tray with integrated screw

S TR-00024.5



### Abutment Screw

Recommended tightening torques 25 Ncm

S 61



### DIM-ANALOG®

S 5.DIM.455



### Sleeve (for DIM-ANALOG®)

for plaster model


DIM 5.XSL.A20

# S-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



4,5 / 5,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS  
ONLINE



### nt-Preform®

incl. final Abutment Screw

GH 0,5 mm S 10.Mi4.550.05



### Scanbody 3D-Guide®

for intra + extra oral scanning

S 9.S3D4.550



### nt-Mill Lab Analog

for nt-Preform®

S 5.La4.550



### colored milling screw

for nt-Preform®

S 6.Sc4.550



### Abutment Screw

Recommended tightening torques: 25 Ncm

S 61



### DIM-ANALOG®

S 5.DIM.455



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSLA20



### Abutment Grip

for nt-Preform® Holder

W12.200.00



### Customized Abutment

milled by NT-iDent / Not available in the USA

TZ15.TiO.Si1


# S-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**





4,5 / 5,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm S 820 S

GH 2,5 mm S-2CON 8045502S



### 2-CONnect®-Abutment

Male Part, recommended tightening torques 25 Ncm

GH 0,5 mm S 820 M

GH 2,5 mm S-2CON8045502M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 61



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

# S-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



4,5 / 5,0 mm

Dentsply Implants® / Astra Tech Implant System OsseoSpeed® TX

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm S-2CON 80455A3S



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm S-2CON 80455A3M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



### Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20  
Ncm

S 6.Ms4550



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM160



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

# S-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

MADE IN  
GERMANY

WWW.NT.DENTAL

## NOTES

[illegible]

## EUROPE

### **NT-Trading GmbH & Co. KG**

G.-Braun-Str. 18  
76187 Karlsruhe  
Germany

**E** [info@nt.dental](mailto:info@nt.dental)

**W** [www.nt.dental](http://www.nt.dental)

## NORTH AMERICA

### **NT-US, Inc.**

132 Pepes Farm Road,  
Unit 3  
Milford, CT 06460  
USA

**E** [info@nt.dental](mailto:info@nt.dental) (general information)

**E** [order-us@nt.dental](mailto:order-us@nt.dental) (for orders)

**W** [www.nt.dental](http://www.nt.dental)

2023\_NT\_Serienbroschüre\_S\_EN\_REV002