

# THO -SERIES



compatible to:  
Thommen Medical® / SPI®



3,5 mm



4,0 mm




4,5 mm



5,0 mm




6,0 mm



**PRECISION  
MEETS QUALITY.**

**OVERVIEW  
PLATFORMS:**



3,5 mm

---



4,0 mm

---



4,5 mm

---



5,0 mm


---



6,0 mm

---

 **MADE IN  
GERMANY**



3,5 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### **Titanium Base (engaged)**

incl. Abutmentscrew

THO 80350000



### **Scanbody 3D-Guide®**

for intra + extra oral scanning

THO 9.S3D3.500



### **Lab Analog**

THO 5.La3.500



### **Implant Pick.Up**

for open-tray with integrated screw

THO TR-00013.5



### **Abutment Screw**

Recommended tightening torques 15 Ncm

THO 603500



### **DIM-ANALOG®**

THO 5.DIM.350



### **Sleeve (for DIM-ANALOG®)**

for plaster model


DIM 5.XSLA10

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



3,5 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### nt-Preform®

incl. final Abutment Screw

GH 0,5 mm THO 10.Mi3.500.05



### Scanbody 3D-Guide®

for intra + extra oral scanning

THO 9.S3D3.500



### nt-Mill Lab Analog

for nt-Preform®

THO 5.La3.500



### colored milling screw

for nt-Preform®

THO 6.Sc3.500



### Abutment Screw

Recommended tightening torques 15 Ncm

THO 603500



### DIM-ANALOG®

THO 5.DIM.350



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSLA10



### Abutment Grip

for nt-Preform® Holder

W 12.100.000



### Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



3,5 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm THO-2CON 8035001S

GH 2,5 mm THO-2CON 8035002S



### 2-CONnect®-Abutment

Male Part, recommended tightening torques 15 Ncm

GH 0,5 mm THO-2CON 8035001M

GH 2,5 mm THO-2CON 8035002M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.400



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 60



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 60

2CON 5.DIM.140



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00011.4



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.400

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



3,5 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm THO-2CON 80350A3S



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm THO-2CON 80350A3M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



### Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20  
Ncm

THO 6.Ms3.500



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM.160



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

# THO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

MADE IN  
GERMANY

WWW.NT.DENTAL



4,0 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### Titanium Base (engaged)

incl. Abutmentscrew

THO 80400000



### Scanbody 3D-Guide®

for intra + extra oral scanning

THO 9.S3D4.000



### Lab Analog

THO 5.La4.000



### Implant Pick.Up

for open-tray with integrated screw

THO TR-00024.0



### Abutment Screw

Recommended tightening torques 25 Ncm

THO 604060



### DIM-ANALOG®

THO 5.DIM.400



### Sleeve (for DIM-ANALOG®)

for plaster model


DIM 5.XSLA20

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



4,0 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### nt-Preform®

incl. final Abutment Screw

GH 0,5 mm THO 10.Mi4.000.05



### Scanbody 3D-Guide®

for intra + extra oral scanning

THO 9.S3D4.000



### nt-Mill Lab Analog

for nt-Preform®

THO 5.La4.000



### colored milling screw

for nt-Preform®

THO 6.Sc4.060



### Abutment Screw

Recommended tightening torques 25 Ncm

THO 604060



### DIM-ANALOG®

THO 5.DIM.400



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSLA20



### Abutment Grip

for nt-Preform® Holder

W 12.200.000



### Customized Abutment

milled by NT-iDent / Not available in the USA

TZ15.TiO.Si1


# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**





4,0 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm THO-2CON 8040001S

GH 2,5 mm THO-2CON 8040002S



### 2-CONnect®-Abutment

Male Part, recommended tightening torques 25 Ncm

GH 0,5 mm THO-2CON 8040001M

GH 2,5 mm THO-2CON 8040002M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.400



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 60



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.140



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00011.4



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.400

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



4,0 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm THO-2CON 80400A3S



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm THO-2CON 80400A3M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



### Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20  
Ncm

THO 6.Ms4.060



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM.160



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

# THO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

MADE IN  
GERMANY

WWW.NT.DENTAL



4,5 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### **Titanium Base (engaged)**

incl. Abutmentscrew

THO 80450000



### **Scanbody 3D-Guide®**

for intra + extra oral scanning

THO 9.S3D4.500



### **Lab Analog**

THO 5.La4.500



### **Implant Pick.Up**

for open-tray with integrated screw

THO TR-00034.5



### **Abutment Screw**

Recommended tightening torques 25 Ncm

THO 604060



### **DIM-ANALOG®**

THO 5.DIM.450



### **Sleeve (for DIM-ANALOG®)**

for plaster model


DIM 5.XSLA20

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



4,5 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### nt-Preform®

incl. final Abutment Screw

GH 0,5 mm THO 10.Mi4.500.05



### Scanbody 3D-Guide®

for intra + extra oral scanning

THO 9.S3D4.500



### nt-Mill Lab Analog

for nt-Preform®

THO 5.La4.500



### colored milling screw

for nt-Preform®

THO 6.Sc4.060



### Abutment Screw

Recommended tightening torques 25 Ncm

THO 604060



### DIM-ANALOG®

THO 5.DIM.450



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSLA20



### Abutment Grip

for nt-Preform® Holder

W 12.200.000



### Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



4,5 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm THO-2CON 8045001S

GH 2,5 mm THO-2CON 8045002S



### 2-CONnect®-Abutment

Male Part, recommended tightening torques 25 Ncm

GH 0,5 mm THO-2CON 8045001M

GH 2,5 mm THO-2CON 8045002M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 61



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



4,5 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm THO-2CON 80450A3S



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm THO-2CON 80450A3M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



### Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20  
Ncm

THO 6.Ms4.060



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM.160



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

# THO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 MADE IN  
GERMANY

[WWW.NT.DENTAL](http://WWW.NT.DENTAL)



5,0 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### **Titanium Base (engaged)**

incl. Abutmentscrew

THO 80500000



### **Scanbody 3D-Guide®**

for intra + extra oral scanning

THO 9.S3D5.000



### **Lab Analog**

THO 5.La5.000



### **Implant Pick.Up**

for open-tray with integrated screw

THO TR-00045.0



### **Abutment Screw**

Recommended tightening torques 25 Ncm

THO 604060



### **DIM-ANALOG®**

THO 5.DIM.500



### **Sleeve (for DIM-ANALOG®)**

for plaster model


DIM 5.XSLA30

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



5,0 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### nt-Preform®

incl. final Abutment Screw

GH 0,5 mm THO 10.Mi5.000.05



### Scanbody 3D-Guide®

for intra + extra oral scanning

THO 9.S3D5.000



### nt-Mill Lab Analog

for nt-Preform®

THO 5.La5.000



### colored milling screw

for nt-Preform®

THO 6.Sc4.060



### Abutment Screw

Recommended tightening torques: 25 Ncm

THO 604060



### DIM-ANALOG®

THO 5.DIM.500



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSLA30



### Abutment Grip

for nt-Preform® Holder

W 12.300.000



### Customized Abutment

milled by NT-iDent / Not available in the USA

TZ15.TiO.Si1


# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**





5,0 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm THO-2CON 8050001S

GH 2,5 mm THO-2CON 8050002S



### 2-CONnect®-Abutment

Male Part, recommended tightening torques 25 Ncm

GH 0,5 mm THO-2CON 8050001M

GH 2,5 mm THO-2CON 8050002M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 61



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



5,0 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm THO-2CON 80500A3S



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm THO-2CON 80500A3M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



### Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20  
Ncm

THO 6.Ms4.060



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM.160



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts


2CON 5.La1.600

# THO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

MADE IN  
GERMANY

WWW.NT.DENTAL



6,0 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### **Titanium Base (engaged)**

incl. Abutmentscrew

THO 80600000



### **Scanbody 3D-Guide®**

for intra + extra oral scanning

THO 9.S3D6.000



### **Lab Analog**

THO 5.La6.000



### **Implant Pick.Up**

for open-tray with integrated screw

THO TR-00056.0



### **Abutment Screw**

Recommended tightening torques 25 Ncm

THO 604060



### **DIM-ANALOG®**

THO 5.DIM.600



### **Sleeve (for DIM-ANALOG®)**

for plaster model


DIM 5.XSLA30

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



6,0 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### nt-Preform®

incl. final Abutment Screw

GH 0,5 mm THO 10.Mi6.000.05



### Scanbody 3D-Guide®

for intra + extra oral scanning

THO 9.S3D6.000



### nt-Mill Lab Analog

for nt-Preform®

THO 5.La6.000



### colored milling screw

for nt-Preform®

THO 6.Sc4.060



### Abutment Screw

Recommended tightening torques 25 Ncm

THO 604060



### DIM-ANALOG®

THO 5.DIM.600



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSLA30



### Abutment Grip

for nt-Preform® Holder

W 12.300.000



### Customized Abutment

milled by NT-iDent / Not available in the USA


TZ15.TiO.Si1

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



6,0 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (straight)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw

GH 0,5 mm THO-2CON 8060001S

GH 2,5 mm THO-2CON 8060002S



### 2-CONnect®-Abutment

Male Part, recommended tightening torques 25 Ncm

GH 0,5 mm THO-2CON 8060001M

GH 2,5 mm THO-2CON 8060002M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.800



### LOcON COPING

For more information, see our LOcON brochure

LOC 852001



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 20 Ncm

KS 61



### 2-CONnect®-DIM-ANALOG®

for straight Male Part, KS 61

2CON 5.DIM.180



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00021.8



### 2-CONnect®-lab Analog

for straight Male Parts


2CON 5.La1.800

# THO-SERIES

**Notice:** Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

 **MADE IN  
GERMANY**

**WWW.NT.DENTAL**



6,0 mm

Thommen Medical® SPI®

ALL DETAILS  
ONLINE



### 2-CONnect®-Abutment Set (angled)

incl. 2-CONnect®-Coping Female Part and 2-CONnect®-Coping-screw,  
incl. Male Part screw, mounting aid, angled 20°

GH 2,5 mm THO-2CON 80600A3S



### 2-CONnect®-Abutment (angled)

incl. 2-CONnect® Male Part screw, mounting aid, angled 20°

GH 2,5 mm THO-2CON 80600A3M



### 2-CONnect® Scanbody 3D-Guide®

for intra + extra oral scanning

2CON 9.S3D1.600



### 2-CONnect®-Coping

Female Part

2CON F-U



### 2-CONnect®-Coping-Screw

Recommended tightening torques 15 Ncm

KS 62



### Male Part Screw

for angled 2-CONnect® Abutment, recommended tightening torques 20  
Ncm

THO 6.Ms4.060



### 2-CONnect®-DIM-ANALOG®

for angled 2-CONnect® Abutment Male Parts

2CON 5.DIM.160



### Sleeve (for DIM-ANALOG®)

for plaster model

DIM 5.XSL.A30



### 2-CONnect®-Implant pick-up

for open-tray with integrated screw

2CON TR-00031.6



### 2-CONnect®-lab Analog

for angled 2-CONnect® Abutment Male Parts

2CON 5.La1.600

# THO-SERIES

Notice: Products indicated with ® and ™ are registered brand names of the respective manufacturer.  
For products available in the USA, please check our website: [www.nt.dental/en/](http://www.nt.dental/en/)

MADE IN  
GERMANY

WWW.NT.DENTAL

## NOTES

[illegible]

## EUROPE

### **NT-Trading GmbH & Co. KG**

G.-Braun-Str. 18  
76187 Karlsruhe  
Germany

**E** [info@nt.dental](mailto:info@nt.dental)

**W** [www.nt.dental](http://www.nt.dental)

## NORTH AMERICA

### **NT-US, Inc.**

132 Pepes Farm Road,  
Unit 3  
Milford, CT 06460  
USA

**E** [info@nt.dental](mailto:info@nt.dental) (general information)

**E** [order-us@nt.dental](mailto:order-us@nt.dental) (for orders)

**W** [www.nt.dental](http://www.nt.dental)

2023\_NT\_Serienbroschüre\_THO\_EN\_REV001