



2-CONnect®

Multiunit Abutment, available in various gingival heights and angulations for a variety of common implant systems.



Dateiplat: 2-CONnect_EN_V002 Seite 1/1



BENEFITS

- self-centering cone connection
- proven flexibility for divergent implant positions, can be additionally combined with NT LOcON (Locator® type attachment)
- tension-free fit through optional adhesive connection (female part)
- applicable for conventional processing or CAD/CAM technology
- optimized abutment design enables even distribution of occurring forces
- female part (cap) covers all common implant systems and platforms



The NT LOcON attachment is an ingenious concept that redefines flexibility and efficiency.

APPLICATION

- bars
- bridges
- prosthetics in connectivity with the NT LOcON system

TECHNICAL DETAILS

- Ti6Al4V ELI Grade 23

