



# DIM-ANALOG®

Digital Implant Model Analog - available for a variety of common implant systems within the NT-IQ implant libraries.



Dateplatt\_DIM-ANALOG\_230818\_EN\_V002 - Seite 1/1

## BENEFITS

- implant analog for digital model production within the modelbuilder CAD applications
- two-part repositionable implant analog
- exact positioning through guidance elements
- position correction possible
- design supports exact repositioning of removable gingival parts
- sleeve for optional use in conventional plaster models

## APPLICATION

- For the production of implant models in the digital workflow.

## TECHNICAL DETAILS

- Ti6Al4V ELI Grade 23
- stainless steel

