



3D-Guide®

Scanbody for high-precision transfer of the implant position. Available for a wide range of common implant systems.



Datenblatt: 3DGuide_230818_EN_V002 Seite 1/1

BENEFITS

- exact positioning via metal interface
- suitable for intraoral scanners
- the different geometries of the lateral surfaces guarantee an exact repositioning in the CAD software is even with incomplete capture
- integrated aspiration-proof screw
- x-ray detectable
- rotatable in case of tooth crowding
- series and platform can be easily identified by laser engraving

APPLICATION

- Scanbody for model scan
- Scanbody for intraoral scan

TECHNICAL DETAILS

- 2-3 µm manufacturing tolerance

WWW.NT.DENTAL

NT-Trading GmbH & Co. KG, G.-Braun-Straße 18 · 76187 Karlsruhe · Germany

NT DIGITAL
DENTAL
TECHNOLOGY