



TITANIUM BASE

For the production of customized
implant-supported hybrid abutments.
Available for a wide range of common
implant systems.



Datenblatt_Titanbasis_230818_EN_V002 Seite 1/1



BENEFITS

- engaging
- individual abutment design within the CAD software via the NT-IQ implant libraries
- large-area bonding, therefore permanent and stable connection to the base
- applicable for conventional processing or CAD/CAM technology

APPLICATION

- to bond individually manufactured superstructures from ceramic materials (zirconium dioxide, lithium disilicate) or plastics.

TECHNICAL DETAILS

- Ti6Al4V ELI Grade 23

 THE SECURE CONNECTION
PRECISE. PROVEN. RELIABLE.

