



TITANIUM BASE NX

Titanium Base (NoHex) for the production of individual bridge/bar constructions on multiple implants.



BENEFITS

- non-engaging
- can be shortened from 5 mm to 3 mm, considering different occlusal space conditions
- available in two gingival heights
- optimized emergence profile, very good adaptation to existing gingiva conditions
- specially developed for CAD/CAM technology

APPLICATION

- For individual bridge/bar constructions on several implants made of ceramic materials (zirconium dioxide, lithium disilicate) or plastics (PMMA, PEEK, composite).

TECHNICAL DETAILS

- Ti6Al4V ELI Grade 23