




NT X-tra Holder Set

The multi-part X-tra Holder Set protects the implant interface from damage during the finishing/processing phase.



Date: 2018-01-11, NTX-100-Holder-Set_230818_EN_V002 Seite 1/1

BENEFITS

- solid stainless steel grip for optimal heat dissipation during finishing/processing
- includes three inserts to cover a large number of common implant systems on three mill lab analog sizes
- easy and quick exchange of mill lab analogs and inserts

SET INCLUDES

- handle
- screwdriver for handle
- insert for handle Ø 3.5 mm
- insert for handle Ø 4.3 mm
- insert for handle Ø 5.0 mm

