


# K-SERIES



compatible to:  
Nobel Biocare™ / Brånemark®




3,5 mm NP



5,1 mm WP




4,1 mm RP



**PRECISION  
MEETS QUALITY.**

**OVERVIEW  
PLATFORMS:**



3,5 mm NP




4,1 mm RP



5,1 mm WP

 **MADE IN  
GERMANY**



3,5 mm NP  
Nobel Biocare™ Brånemark®



**Titanium Base (engaged)**  
incl. Abutmentscrew

K 800



**Scanbody 3D-Guide®**  
for intra + extra oral scanning

K 9.S3D3.500



**Lab Analog**

K 50



**Implant Pick Up**  
for open-tray with integrated screw

K TR-00013.5



**Abutment Screw**  
recommended tightening torques 35 Ncm

K 60




**DIM-ANALOG®**

K 5.DIM.350











**Sleeve (for DIM-ANALOG®)**  
for plaster model

DIM 5.XSLA10



4,1 mm RP  
Nobel Biocare™ Brånemark®

	<b>Titanium Base (engaged)</b> incl. Abutmentscrew	K 810
	<b>Scanbody 3D-Guide®</b> for intra + extra oral scanning	K 9.S3D4.100
	<b>Lab Analog</b>	K 51
	<b>Implant Pick Up</b> for open-tray with integrated screw	K TR-00024.1
	<b>Abutment Screw</b> recommended tightening torques 35 Ncm	K 61
	<b>DIM-ANALOG®</b>	K 5.DIM.410
	<b>Sleeve (for DIM-ANALOG®)</b> for plaster model	DIM 5.XSLA10



5,1 mm WP  
Nobel Biocare™ Brånemark®



**Titanium Base (engaged)**  
incl. Abutmentscrew

K 820



**Scanbody 3D-Guide®**  
for intra + extra oral scanning

K 9.S3D5.100



**Lab Analog**

K 52



**Abutment Screw**  
recommended tightening torques 35 Ncm

K 62



**DIM-ANALOG®**

K 5.DIM.510



**Sleeve (for DIM-ANALOG®)**  
for plaster model

DIM 5.XSL.A30



## EUROPE

### **NT-Trading GmbH & Co. KG**

G.-Braun-Str. 18  
76187 Karlsruhe  
Germany

**E** [info@nt.dental](mailto:info@nt.dental)

**W** [www.nt.dental](http://www.nt.dental)

## NORTH AMERICA

### **NT-US, Inc.**

132 Pepes Farm Road,  
Unit 3  
Milford, CT 06460  
USA

**E** [info@nt.dental](mailto:info@nt.dental) (general information)

**E** [order-us@nt.dental](mailto:order-us@nt.dental) (for orders)

**W** [usa.nt.dental](http://usa.nt.dental)

2024\_NT\_US\_Serienbroschüre\_K\_EN\_REV01