


DIM ANALOG[®]

*Digital Implant Model
Analog – available for
24 implant systems in
the nt-IQ library*



FEATURES:

- implant analog for digital model manufacturing in the model builder library
- two-piece repositionable implant analog
- precise positioning
- position adjustment possible
- analog design supports precise repositioning of removable gingiva mask
- Sleeve for optional usage in conventional plaster models




APPLICATION:

- for model manufacturing in digital workflow



TECHNICAL DETAILS:

- Stainless Steel
- Titanium Grade 5 ELI



DIM PARALLEL DRILL
for plaster model

W 16.DIM.000

